

SPECIAL COMPACT HT

Public price starting from £ 4,650,00

M.D.R.: **2544485**

Basid UDI-DI: **805378246CHT045A3**

Order Date

Quote

Dealer

Order reference

MODEL FEATURES

- Triride Software
- Customisable electronics
- Display with 5 levels of power and speed
- Electronic brake
- Reverse
- Cruise Control
- Assisted load
- Weight: 11 kg
- Range: up to 50 km
- Torque: 75Nm
- Lithium-ion battery 36V or 48V
- Tyre City or Cross
- ECODRIVE-Energy saving mode
- Two driving modes: Normal/Sport
- Charging of the battery during downhill and braking
- Light activated from display

Maximum compactness with a new exclusive, powerful high performance motor, combined with the sophisticated Made in Triride electronics, it guarantees maximum performance to satisfy even the most demanding users. Perfect for tackling the more demanding inclines, and more uneven routes with the Allroad tyre option.

HIGH TORQUE ALWAYS AT YOUR DISPOSAL



£4,650.00




12" City

£4,650.00



12" Cross

£4,750.00



Allroad

Carbon mudguard default

£4,750.00



Supermotard

Carbon mudguard default

WHEELCHAIR FEATURES

Brand

Model

Front tube's diameter

Seat width

Seat depth (E)

Toe to front tube distance (A)

Ground/knee height (D)

Length of stretched arms (C)

Front tubes angle (degree)

Frame length (B)

Front tubes width(F)*

User's weight

Footrest

Fixed ☐

Removable ☐

Internal ☐

External ☐

USER'S PERSONAL DATA

Name

Tel.:

E-mail:

Age Pathology:

DELIVERY ADDRESS

Delivery to:

Postcode


Address

Town

County

* The measurements for the width of the front tubes must be taken from the outside of the tubes at a height between 30 and 35cm

All fields are obligatory



N.B.: For wheelchairs with non-standard measurements or for users with special needs and dimensions the frame and driving controls can be customized with a quote (eg: dwarfism, bariatrics, etc)

Assembly, installation and instruction for use carried out by a Triride specialist technician

(on request) ☐

MTW: TREKKING POWER

Adaptation for MTW Trekking Power Kit purchase (complete and attach dedicated form)	no cost
T0013 Adaptation for MTW Trekking Power (future purchase)	£185.00

COMMANDS FOR TETRAPLEGIC USERS

Configuration for Tetraplegic Commands N. _____
(For the desired configuration fill and attach the dedicated form)

Personalised configuration (attach picture or drawing)

"For more information, see the attached page with various configurations"

TRIRIDE CONFIGURATION

Battery

BS-R-369	Battery 36V 7 USB (252 Wh)	default
BS-R-394	Battery 36V 9 USB (324Wh)	£0.00
BS-R-393-S	Battery 36V 14 USB (504 Wh)	£198.00
BS-R-368-S	Battery 48V 6.2 USB (297 Wh)	£153.00
BS-R-395-S	Battery 48V 10.5 USB (504 Wh)	£225.00

Electronics

BS-O-096	Performance Pack : IBS+ICC+IBS button	£171.00
BS-O-101	IBS - Intelligent Braking System	£68.00
BS-O-097	ICC - Intelligent Cruise Control	£68.00
BS-O-086	ECODRIVE - Increased battery range	£68.00
BS-O-113	MDC - Motion Direct Control	£132.00
BS-R-366	Electronic brake button / IBS	£44.00
BS-O-076	Analogue brake with thumb control	£95.00
BS-O-077	Automatic brake on accelerator release	£41.00
BS-O-079	Emergency brake with remote control	£95.00

Display LCD

BS-C-077	Display LCD Special with USB	default
BS-C-078	Display LCD Colour with USB	£53.00
BS-C-079	Display Colour with USB with customisable Triride programming	£126.00

customisation: _____ (21 character)

Notes



Throttle

BS-R-299	Thumb throttle ambidextrous	☐ LHS ☐ RHS	£0.00	☐
L: BS-R-349				
R: BS-R-348	Palm throttle with long lever	☐ LHS ☐ RHS	£20.00	☐
L: BS-R-333				
R: BS-R-335	Half twist throttle	☐ LHS ☐ RHS	£38.00	☐
L: T0062				
R: BS-R-324	Full twist throttle	☐ LHS ☐ RHS	£38.00	☐

Brakes

L: BS-R-319				
R: BS-R-318	Single mechanical brake	☐ LHS ☐ RHS	£0.00	☐
L: BS-O-203				
R: BS-O-273	Single mechanical parking brake	☐ LHS ☐ RHS	£14.00	☐
L: BS-O-274	(not available with hydraulic brake)			
R: BS-O-275	Single hydraulic brake	☐ LHS ☐ RHS	£122.00	☐
BS-O-218	Mechanical brake - double caliper		£149.00	☐
BS-O-194	Hydraulic brake - double caliper		£224.00	☐

N.B: The electronic brake (including IBS when required) can be actioned also by the brake levers of the device

Handlebar

BS-O-246	Straight (indicate width between 46cm and 50cm)		cm	£0.00	☐
BS-O-208	4 Bends riser			£0.00	☐
BS-O-107	4 Bends high riser			£0.00	☐

Handlebar stem

BS-O-222	Adjustable 110 mm		£0.00	☐
BS-O-220	Adjustable 130 mm		£0.00	☐
BS-O-185	Adjustable 160 mm		£81.00	☐

Custom frame

BS-O-189	Extended frame (to bring the handlebar closer)		£81.00	☐
----------	--	--	--------	--------------------------

Frame colour

BS-O-078	Titanium		default	☐
BS-O-078-P	Polished		£252.00	☐
BS-O-217-N	Black Ral 9005		£135.00	☐
BS-O-217-R	Red Ral 3020		£135.00	☐
BS-O-217-W	White Ral 9003		£135.00	☐
BS-O-217-G	Yellow Ral 1003		£135.00	☐
BS-O-217-A	Orange Ral 2004		£135.00	☐
BS-O-217-B	Blue Ral 5017		£135.00	☐
BS-O-217	RAL on request		£207.00	☐



Personalized frame wraps

BS-O-087-M	Camouflage	£504.00	☐
BS-O-087-C	Carbon Effect	£504.00	☐
BS-O-087-F	Graffiti wrap Black/White	£504.00	☐
	Personalisation: Two-tone, gradients, etc.	on request	☐

Mudguard

BS-O-245	Mudguard Black	(not for Allroad and Supermotard)	default	☐
BS-O-167	Mudguard Blue	(not for Allroad and Supermotard)	£0.00	☐
BS-O-169	Mudguard Orange	(not for Allroad and Supermotard)	£0.00	☐
BS-O-166	Mudguard Green	(not for Allroad and Supermotard)	£0.00	☐
BS-O-164	Mudguard Red	(not for Allroad and Supermotard)	£0.00	☐
BS-O-215	Mudguard Carbon	(default on Allroad and Supermotard)	£96.00	☐

Handlescrews

BS-O-301	Handlescrews - Black	£0.00	☐
BS-O-128	Handlescrews - Blue	£0.00	☐
BS-O-127	Handlescrews - Orange	£0.00	☐
BS-O-130	Handlescrews - Green	£0.00	☐
BS-O-129	Handlescrews - Red	£0.00	☐
BS-O-139	Lever handlescrews for Tetraplegic user	£5.00	☐

SPEED SETTING

We herewith advise, that using the Triride with max. speed settings higher than 4mph is only permitted on public roads in certain countries and otherwise may only be used on private property. In the absence of approval, in accordance with the road traffic laws, the Triride with speed-settings >4mph may not be driven on public roads, bicycle-lanes and pavements. You need to inform yourself about the country specific legislation regarding usage and your legal obligation to insure the Triride with the speed-settings greater than 4mph.

When placing an order for a max. speed-setting higher than 4mph please confirm with your signature below, that you have understood the above information, you will follow strictly the country specific legislation and you will intend the usage of your Triride on private property.

Up to 4mph (6 Kph) ☐

Up to 8mph (13 Kph) ☐

>8mph (13kph) ☐

Date

Customer Name

Customer signature

OPTIONAL

Attachment Kit for SECOND WHEELCHAIR - please provide details below (prices vary as per specifics)

Tilting frame ☐

Connecting arms ☐

Clamps ☐

Brand

Model

Seat width cm

Seat depth cm

Front tubes diameter mm

ACCESSORIES / OPTIONAL

BS-R-369	Battery 36V 7A USB (252Wh)	£549.00	
BS-R-394	Battery 36V 9A USB (324Wh)	£567.00	
BS-R-393	Battery 36V 14 USB (504 Wh)	£702.00	
BS-R-368	Battery 48V 6.2 USB (297 Wh)	£651.00	
BS-R-395	Battery 48V 10.5 USB (504 Wh)	£753.00	
BS-O-012	Neoprene 20/19 battery cover	£24.00	
BS-R-294	Front spare battery holder with switch	£224.00	
BS-O-083	Central arch stand	£149.00	
BS-O-278	Lateral stand with wheels	£113.00	
BS-O-229	Rearview mirror	£29.00	LHS RHS
BS-O-251	Carry Trailer for fixed frame	£279.00	Rear axle diameter mm
BS-O-019	Handlebar bag with B-Twin adaptor 10 L	£144.00	
BS-O-020	Handlebar bag with B-Twin adaptor 6 L	£144.00	
BS-O-011	Handlebar basket with fabric cover 24 L	£144.00	
BS-O-214	Leg cover	£19.00	
BS-O-292	Calf support with adaptation to tilting frame	£142.00	
BS-O-213	Triride bag	£153.00	
BS-O-114	Arms stopper rings, pair	£16.00	
BS-O-174	Laser customised engraving on frame (on inox and black frame only)	£68.00	
BS-O-156	Digital engraving	£68.00	
BS-O-157	Personalised sticker	£59.00	
BS-O-081	Integrated front light, activated from display	default	
BS-O-022	Rear light with handlebar switch (combined with the front light)	£95.00	
BS-O-232	Light - Front 1200 Lumen	£100.00	
BS-O-198	Light Front USB	£29.00	
BS-O-180	Light - Rear USB	£33.00	
BS-O-015	Rear axle extension for wheelchair (please complete dedicated form)	on request	

Prices are intended for public and before tax

For details and photos of accessories and spare parts, see the relevant "spare parts and accessories" price list or visit our website [web www.triride.uk.com](http://www.triride.uk.com)

N.B.: The Company reserves the exclusive right to bring changes, without notice, if it deems it useful for the improvement of the Product