



Model "LIGHT"

£3,495.00

Medical Devices Registration No.: 2168023

ORDER


QUOTE

ORDER REFERENCE

DATE

Please complete and send this form by e-mail to: orders@triride.uk

CODE	DESCRIPTION	
TR-LIGHT	Wheel 10", 12" or 14", programmable and customisable Triride electronics, 5 speed, reverse gear, cruise control, assisted lifting, electronic brake, display LCD, Cylindrical battery 36V (252 Wh)	
USER'S PERSONAL DATA*		
Name	<input type="text"/>	
Surname	<input type="text"/>	
Phone	<input type="text"/>	
E-mail	<input type="text"/>	
Age	<input type="text"/>	
Pathology	<input type="text"/>	
DELIVERY ADDRESS		
Delivery to	<input type="text"/>	
Street	<input type="text"/>	
Town	<input type="text"/>	
Postcode	<input type="text"/>	
WHEELCHAIR FEATURES		
Brand	<input type="text"/>	
Model	<input type="text"/>	
Front tubes diameter	<input type="text"/>	mm
Seat width	<input type="text"/>	cm
Seat depth (b)	<input type="text"/>	cm
Toe to front tube distance (a)*	<input type="text"/>	cm
Ground/knee height (d)	<input type="text"/>	cm
Length of arms from the shoulder back (c)	<input type="text"/>	cm
Front tubes angle	<input type="text"/>	degree
User's weight	<input type="text"/>	kg
Footrest	Fixed Internal	
	Removable External	
*For wheelchairs with removable footrests, measure as per "a" from toe to front tube.		



TRIRIDE CONFIGURATION

Special commands for tetraplegic users	Configuration for Tetraplegic commands _____ Please see the dedicated attachment with different configurations and prices Personalized configuration (please attach picture or drawing)				on request on request
Battery	BS-R-344	Battery 36V 7A (252 Wh) cylindrical			standard
	BS-R-369	Battery 36V 7A (252 Wh)			£135.00
	BS-R-394	Battery 36V 9 A (324Wh)			£216.00
	BS-R-393-S	Battery 36V 14 A (504 Wh)			£342.00
	BS-R-368-S	Battery 48V 6.2 A (297 Wh)			£396.00
	BS-R-395-S	Battery 48V 10.5 A (504 Wh)			£527.00
Electronics	BS-O-096	Performance Pack: IBS + IBS button + ICC + ECODRIVE			£205.00
	BS-O-101	IBS - Intelligent Braking System (included with tetra controls)			£68.00
	BS-O-097	ICC-Intelligent Cruise Control (included with tetra controls)			£68.00
	BS-O-086	ECODRIVE - Increased battery range			£68.00
	BS-O-113	MDC- Motion Direct Control			£132.00
	BS-R-366	Electronic brake button/IBS			£44.00
	BS-O-076	Analog brake with thumb control (not compatible with MDC)			£95.00
	BS-O-077	Automatic brake on accelerator release (included with push&pull)			£41.00
	BS-O-079	Emergency brake with remote control (not compatible with MDC)			£95.00
Display LCD	BS-C-076	Display LCD Basic			standard
	BS-C-010	Display LCD Mini			standard
	BS-C-077	Display LCD Special with USB			£53.00
	BS-C-011	Display LCD Special colour without USB			£53.00
	BS-C-078	Display LCD Colour with USB			£99.00
Throttle	BS-R-299	Thumb throttle	LHS	RHS	standard
	LS:BS-R-349 RS:BS-R-348	Palm throttle with long lever	LHS	RHS	£20.00
	LS: BS-R-333 RS:BS-R-335	Half twist throttle	LHS	RHS	£38.00
Brakes	LS:BS-R-319 RS:BS-R-318	Single mechanical brake	LHS	RHS	standard
	LS:BS-O-203 RS:BS-O-273	Single mechanical Parking brake	LHS	RHS	£14.00
	LS:BS-O-274 RS:BS-O-275	Single hydraulic brake	LHS	RHS	£122.00
	BS-O-218	Mechanical brake - double caliper			£149.00
	BS-O-194	Hydraulic brake - double caliper			£224.00

N.B: The electronic brake (including IBS when required) can be actioned also by the brake levers of the device					
Wheel	BS-R-280	10" wheel city tyre			£162.00
	BS-R-353	12" wheel city long life tyre (extra spokes for 12"rim)			£32.00
	BS-R-297	12" wheel cross tyre (extra spokes for 12" rim)			£38.00
	BS-R-301	14" wheel city long life tyre			standard
Handlebar	BS-O-246	Straight			standard
	BS-O-208	4 bends riser			standard
	BS-O-107	4 bends high riser			standard
Handlebar stem	BS-O-222	Adjustable 110 mm			standard
	BS-O-220	Adjustable 130 mm			standard
	BS-O-185	Adjustable 160 mm			£81.00
Custom frame	BS-O-189	Dedicated long version			£81.00
Frame colour	BS-C-078	Neutral (titanium steel colour)			standard
	BS-C-078-P	Polished			£189.00
	BS-O-217-N	Black Ral 9005			£135.00
	BS-O-217-R	Red Ral 3020			£135.00
	BS-O-217-W	White Ral 9003			£135.00
	BS-O-217-G	Yellow Ral 1023			£135.00
	BS-O-217-A	Orange Ral 2004			£135.00
	BS-O-217-B	Blue Ral 5017			£135.00
	BS-O-217	RAL on request	RAL N. _____		£185.00
Frame special wraps	BS-O-087-M	Camouflage			£504.00
	BS-O-087-C	Carbon effect			£504.00
	BS-O-087-F	Graffiti wrap			£504.00
Frame customisation	BS-O-174	Laser customised engraving on frame (on inox and black frame only)			£68.00
	BS-O-156	Digital engraving			£68.00
	BS-O-157	Custom sticker			£59.00
Mudguard Black	BS-O-245	standard	Triangular Screws - Black	BS-O-301	standard
Mudguard Blue	BS-O-167	standard	Triangular Screws - Blue	BS-O-128	standard
Mudguard Orange	BS-O-169	standard	Triangular Screws - Orange	BS-O-127	standard
Mudguard Green	BS-O-166	standard	Triangular Screws - Green	BS-O-130	standard
Mudguard Red	BS-O-164	standard	Triangular Screws - Red	BS-O-129	standard
Mudguard Carbon	BS-O-215	£96.00	Lever Handlescrews	BS-O-330	standard
			Lever Handlescrews for Tetraplegic user	BS-O-139	£5.00

OPTIONAL

Attachment Kit for SECOND WHEELCHAIR - please provide details below for quotation (prices vary as per specifics)

Tilting frame _____ Connecting arms _____ Clamps _____

Brand _____ Model _____

Seat width _____ cm Seat depth _____ cm Front tubes diameter _____ cm

Footrest: Internal External Fixed Removable

SPEED SETTING

Up to 6 Kph Up to 15 Kph >15Kph (48V battery only)

We herewith advise, that using the Triride with max. speed settings higher than 6kph is only permitted on public roads in certain countries and otherwise may only be used on private property. In the absence of approval, in accordance with the road traffic laws, the Triride with speed-settings >6 Kph may not be driven on public roads, bicycle-lanes and pavements. You need to inform yourself about the country specific legislation regarding usage and your legal obligation to insure the Triride with the speed-settings greater than 6 Kph.

When placing an order for a max. speed-setting higher than 6 Kph please confirm with your signature below, that you have understood the above information, that you will follow strictly the country specific legislation and you will intend the usage of your Triride on private property.

Date	Customer Name	Customer signature
------	---------------	--------------------

NOTES

ACCESSORIES

36V 7A Battery (252 Wh) cylindrical		BS-R-344	£468.00
36V 7A Battery (252 Wh) USB		BS-R-369	£549.00
36V 9A Battery (324 Wh) USB		BS-R-394	£567.00
36V 14 A Battery (504 Wh) USB		BS-R-393	£702.00
48V 6.2 A Battery (297 Wh) USB		BS-R-368	£651.00
48V 10.5 A Battery (504 Wh) USB		BS-R-395	£753.00
Neoprene 20/19 battery cover		BS-O-012	£24.00
Front spare battery holder with switch		BS-R-294	£225.00
Central stand		BS-O-201	£126.00
Central arch stand		BS-O-083	£149.00
Lateral stand with wheels		BS-O-278	£113.00
Rearview mirror	LHS RHS	BS-O-229	£29.00
**Carry Trailer	Rear axle diameter _____ mm	BS-O-251	£279.00
Handlebar basket with fabric cover (not compatible with double battery holder)		BS-O-011	£144.00
Leg cover		BS-O-214	£19.00
Triride bag		BS-O-213	£153.00
Arms stopper rings, pair		BS-O-114	£16.00
Integrated front light, activated from display (first mount)		BS-O-081	£95.00
Light - Front 1200 Lumen		BS-O-232	£100.00
Light - Front + USB		BS-O-198	£29.00
Light - Rear + USB		BS-O-180	£33.00
Rear axle extension for wheelchair (please complete dedicated form)		BS-O-015	on request

"For a complete list of accessories and spare parts visit our website"

* Mandatory fields for warranty purposes

** Available for fixed frame and some foldable wheelchairs

*** Prices are intended for public and before tax

N. B.:The Company reserves the exclusive right to change, without notice, if it deems is useful for the improvement of the Product