



## TRIBIKE HYBRID PLUS

**£4,500.00** \*\*\*

**ORDER**

**QUOTE**

**ORDER REFERENCE**

**DATE**

Please complete and send this form by e-mail to: [orders@triride.uk](mailto:orders@triride.uk)

### TECHNICAL FEATURES

- |                                                                                                                                                                                                                                                                                                                                                                                                                            |                                                                                                                                                                                                                                                                                                                                 |
|----------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|---------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| <ul style="list-style-type: none"> <li>1 - Electronics by Triride</li> <li>2 - Electronics pack: accelerator, brake ,motor brake, reverse, battery recharge, IBS, ICC</li> <li>3- Battery 20-19 with USB: 36V or 48V</li> <li>4 - Pedal Assist PAS</li> <li>5 - Stainless steel frame</li> <li>6 - Display LCD Special USB with 5 assistance levels</li> <li>7 - 16" or 20" aluminium rim and double inner tube</li> </ul> | <ul style="list-style-type: none"> <li>8 - Shimano 7-speed gearbox</li> <li>9 - Disc brake + V-BRAKE</li> <li>10 - Speed with engine support max 25 Km/h</li> <li>11 - Weight: 11 Kg approx.(without battery)</li> <li>12 - Adjustable pedaling position - N.B.: <b>after each adjustment check the chain length</b></li> </ul> |
|----------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|---------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|


For any special need or customisation please contact TRIRIDE UK

USER'S PERSONAL DATA*	
Name	<input style="width: 80%;" type="text"/>
Surname	<input style="width: 80%;" type="text"/>
Phone	<input style="width: 80%;" type="text"/>
E-mail	<input style="width: 80%;" type="text"/>
Weight	<input style="width: 80%;" type="text"/>
Age	<input style="width: 80%;" type="text"/>
Pathology	<input style="width: 80%;" type="text"/>

DELIVERY ADDRESS	
Delivery to	<input style="width: 80%;" type="text"/>
Street	<input style="width: 80%;" type="text"/>
Town	<input style="width: 80%;" type="text"/>
Postcode	<input style="width: 80%;" type="text"/>

### WHEELCHAIR FEATURES

Brand	<input style="width: 80%;" type="text"/>	
Model	<input style="width: 80%;" type="text"/>	
Front tubes diameter	<input style="width: 40%;" type="text"/>	mm
Seat width	<input style="width: 40%;" type="text"/>	cm
Seat depth (b)	<input style="width: 40%;" type="text"/>	cm
Toe to front tube distance (a)*	<input style="width: 40%;" type="text"/>	cm
Ground/knee height (d)	<input style="width: 40%;" type="text"/>	cm
Length of stretched arms from the back(c)	<input style="width: 40%;" type="text"/>	cm
Footrest	Fixed	Internal
	Removable	External
<p>* For wheelchairs with removable footrests, measure as per "a" from toe to front tube.</p>		



CONFIGURATION				
Wheel	16" wheel	City tyre	Cross tyre	standard
	20" wheel	City tyre	Cross tyre	standard
Battery	BS-R-394	Battery 36V 9Ah USB (324Wh)		standard
	BS-R-393-SB	Battery 36V 14 USB (504 Wh)		£252.00
	BS-R-395-SB	Battery 48V 10.5 USB (504 Wh)		£306.00
Cranck arms	TR-B-032	120mm		standard
	TR-B-021	140mm		standard
	TR-B-004	155mm		standard
Handlebar Stem	BS-O-221	Fixed 60mm		standard
	TR-B-090	Fixed 80mm 40 degrees		£37.00
	TR-B-091	Fixed 100mm 40 degrees		£37.00
	BS-C-441	Fixed 120mm 40 degrees		£37.00
Frame Colour	BS-O-217-N	Black Ral 9005		standard
	BS-O-217-R	Red Ral 3020		standard
	BS-O-217-W	White Ral 9003		£135.00
	BS-O-217-G	Yellow Ral 1023		£135.00
	BS-O-217-A	Orange Ral 2004		£135.00
	BS-O-217-B	Blue Ral 5017		£135.00
	BS-O-217	Customized RAL	RAL no. _____	£185.00
Frame special wraps	BS-O-087-M	Camouflage		£504.00
	BS-O-087-C	Carbon effect		£504.00
	BS-O-087-F	Graffiti wrap		£504.00
Frame customisation	BS-O-174	Laser customised engraving on frame (on inox and black frame only)		£68.00
	BS-O-156	Digital engraving		£68.00
	BS-O-157	Custom sticker		£59.00
Triangular handlescrews	BS-O-301	Triangular handlescrews - black		standard
	BS-O-130	Triangular handlescrews - green		standard
	BS-O-128	Triangular handlescrews - blue		standard
	BS-O-127	Triangular handlescrews - orange		standard
	BS-O-129	Triangular handlescrews - red		standard
Lever handlescrews	BS-O-330	Lever handlescrews		standard
	BS-O-139	Lever handlescrews for Tetraplegic user		£5.00

**NOTES**

---



---

**ACCESSORIES**

36V 9A Battery (324 Wh) USB		BS-R-394	£567.00
36V 14 A Battery (504 Wh) USB		BS-R-393	£702.00
48V 6.2 A Battery (297 Wh) USB		BS-R-368	£651.00
48V 10.5 A Battery (504 Wh) USB		BS-R-395	£753.00
Neoprene 20/19 battery cover		BS-O-012	£24.00
Central stand		BS-O-201	£126.00
Central arch stand		BS-O-083	£149.00
Lateral stand with wheels		BS-O-278	£113.00
Rearview mirror	LHS                      RHS	BS-O-229	£29.00
**Carry Trailer	Rear axle diameter _____mm	BS-O-251	£279.00
Leg cover		BS-O-214	£19.00
Arms stopper rings, pair		BS-O-114	£16.00
Integrated front light, activated from display ( first mount )		BS-O-081	£95.00
Light - Front 1200 Lumen		BS-O-232	£100.00
Light - Front+USB		BS-O-198	£29.00
Light - Rear+USB		BS-O-180	£33.00
Rear axle extension for wheelchair ( please complete dedicated form )		BS-O-015	on request

"For a complete list of accessories and spare parts visit our website"

\* Mandatory fields for warranty purposes

\*\* Available for fixed frame and some foldable wheelchairs

\*\*\* Prices are intended for public and before tax

**N. B.:The Company reserves the exclusive right to change, without notice, if it deems is useful for the improvement of the Product**