

mod.  
**LIGHT**

### MODEL SPECIFIC FEATURES

• Max. power with 36 V / 48 V	670 W/890 W
• Max. torque	40 N/m
• Transport weight <sup>1</sup>	from 8,5 Kg
• Wheel size	12" or 14"
• Lithium-ion battery	36 V o 48 V
• Range <sup>2</sup>	up to 50 Km
• Tyre	City or Cross
• Energy saving mode	Eco Drive

### GENERAL FEATURES

- Stainless steel frame
- Customisable electronics by Triride
- Reverse + Assisted load
- Cruise control
- Aluminum rim
- Electronic regenerative brake
- Disc brake
- LCD Display
- 5 selectable power levels
- Speed compliant with legal standards
- TÜV and CE certifications
- Customisable frame colour or pattern

### Available with customised control configuration for tetraplegic users

Medical Devices registered of the Health Ministry  
TUV Certification No. 15-PEL-0010022-TIC

<sup>1</sup>Weight refers to device excluding battery and arms for wheelchair connection.

<sup>2</sup>Range refers to device with 36 V 14 Ah or 48 V 10,5 Ah battery on a flat route at medium speed of 6 km/h (level 1) with user weight of 70 kg. Range depends on: user weight, road condition, general state of the wheelchair and tires, condition of the battery, speed and type of driving.

<sup>3</sup>USB port, to recharge devices such as smartphones or tablets, is available depending on the chosen battery or display.

Triride srl reserves the right to modify the models with updates and innovations and therefore the images may differ from the final product sold.

TRIRIDE srl - via Massimo D'Antona, 8 – 63812 Montegranaro (FM) - Italy  
Tel: +39 0733 801405 - Mob: +39 327 6266267 - Mail: info@trirideitalia.it

WWW.TRIRIDEITALIA.COM



12" / 14"



36 V / 48 V



670 W / 890 W



8,5 Kg<sup>(1)</sup>



up to 50 Km<sup>(2)</sup>



40 N/m

All models are customisable and configurable with a wide range of accessories and options.

MORE INFORMATION IS AVAILABLE ON OUR WEBSITE [WWW.TRIRIDEITALIA.COM](http://WWW.TRIRIDEITALIA.COM)

