

mod. SPECIAL HP16

MODEL SPECIFIC FEATURES

Max. power with 48 V	1500 W
• Torque	50 N/m
• Transport weight ¹	from 12,8 Kg
• Wheel size	16"
• Lithium-ion battery	48 V
• Range ²	up to 50 Km
• Tyre	City or Cross
• Energy saving mode	Eco Drive
• Dual drive profile as standard	Normal / Sport

GENERAL FEATURES

- Stainless steel frame
- Customisable electronics by Triride
- Reverse + Assisted load
- Cruise control
- Aluminum rim
- Electronic regenerative brake
- Disc brake
- LCD Display
- 5 selectable power levels
- Speed compliant with legal standards
- TÜV and CE certification
- Customisable frame colour or pattern

Available with customised control configuration for tetraplegic users

Medical Devices registered of the Health Ministry TUV Certification No. 15-PEL-0010022-TIC

- $^{\rm I}$ Weight refers to device excluding battery and arms for wheelchair connection. $^{\rm 2}$ Range refers to device with 48 V 10,5 Ah battery on a flat route at medium speed of 6 km/h (level 1) with user weight of 70 kg. Range depends on: user weight, road condition, general state of the wheelchair and tires, condition of the battery, speed and type of
- ³USB port to recharge devices such as smartphones or tablets.

Triride srl reserves the right to modify the models with updates and innovations and therefore the images may differ from the final product sold.

TRIRIDE srl - via Massimo D'Antona, 8 – 63812 Montegranaro (FM) - Italy Tel: +39 0733 801405 - Mob: +39 327 6266267 - Mail: info@trirideitalia.it













48 V

1500 W



12,8 Kg⁽¹⁾

up to 50 Km⁽²⁾

50 N/m

a wide range of accessories and options.

All models are customisable and configurable with

MORE INFORMATION IS AVAILABLE ON OUR WEBSITE WWW.TRIRIDEITALIA.COM





