

SPECIAL LIGHT SPECIAL LIGHT ALLROAD

MODEL SPECIFIC FEATURES

Max. power with 36 V/48 V

• Max. torque

• Transport weight¹

• Wheel size

· Lithium-ion battery

• Range²

Tyre

Energy saving mode

• Dual drive profile as standard

900 W/1200 W

50 N/m

from 8,5 Kg

12" or 14"

36 V or 48 V

up to 50 Km

City or Cross

Eco Drive

Normal / Sport

ALLROAD VERSION

- Tyre 12"x3.00 8
- · Carbon fender as standard

GENERAL FEATURES

- Stainless steel frame
- Customisable electronics by Triride
- Reverse + Assisted load
- Cruise control
- Aluminum rim
- Disc brake + Electronic regenerative brake
- LCD Display with 5 selectable power levels
- TÜV and CE certification
- Customisable frame colour or pattern

Available with customised control configuration for tetraplegic users

Medical Devices registered of the Health Ministry TUV Certification No. 15-PEL-0010022-TIC

¹Weight refers to device excluding battery and arms for wheelchair connection.

 2 Range refers to device with 36 V 14 Ah or 48 V 10,5 Ah battery on a flat route at medium speed of 6 km/h (level 1) with user weight of 70 kg. Range depends on: user weight, road condition, general state of the wheelchair and tires, condition of the battery, speed and type of driving.

³ USB port to recharge devices such as smartphones or tablets.

Triride srl reserves the right to modify the models with updates and innovations and therefore the images may differ from the final product sold.

TRIRIDE srl - via Massimo D'Antona, 8 – 63812 Montegranaro (FM) - Italy Tel: +39 0733 801405 - Mob: +39 327 6266267 - Mail: info@trirideitalia.it















12" / 14"

36 V / 48 V

900 W /1200 W







8,5 Kg⁽¹⁾

up to 50 Km⁽²⁾

50 N/m

All models are customisable and configurable with a wide range of accessories and options.

MORE INFORMATION IS AVAILABLE ON OUR WEBSITE WWW.TRIRIDEITALIA.COM





