

TRIWHEEL

Public price starting from **£950.00**

Order Dealer
Quote Date

Dealer's stamp

MODEL FEATURES

- Weight: 4 kg
- Wheel size 12", 14" or 16"

Non-motorised front wheel with high smoothness. Ideal for having fun tackling downhill descents

FOR THOSE WHO WANT TO PUSH



Notes

WHEELCHAIR DETAILS

Brand

Model

Front tube diameter

Seat width

Seat depth (E)

Toe to front tube distance (A)

Ground/knee height (D)

Length of stretched arms (C)

Front tubes angle (degree)

Frame length (B)

Front tubes width(F)*

User's weight

Legrest

Fixed

Internal

Removable

External

USER'S PERSONAL DATA

Name & Surname


Tel.:

E-mail:

Age Pathology:

Notes

* Measurements for the width of the front tubes must be taken from the outside of the tubes at a height between 30 and 35cm
All fields are obligatory



N.B.: For wheelchairs with non-standard measurements or for users with special needs and dimensions the frame and driving controls can be customized with a quote (eg: dwarfism, bariatrics, etc)

Assembly, installation and instruction for use carried out by a Triride specialist technician on request

TRIWHEEL CONFIGURATION

Handlebar

BS-O-246	Straight (width of your choice between 46cm and 50cm)	£ 0,00	<input type="checkbox"/>
BS-O-208	4 Bends riser	£ 0,00	<input type="checkbox"/>
BS-O-107	4 Bends high riser	£ 0,00	<input type="checkbox"/>

Handlebar stem

BS-O-222	Adjustable 110 mm	£ 0,00	<input type="checkbox"/>
BS-O-220	Adjustable 130 mm	£ 0,00	<input type="checkbox"/>
BS-O-185	Adjustable 160 mm	£ 90,00	<input type="checkbox"/>

Tyre

BS-R-353	Tyre 12" City Long Life	£ 32,00	<input type="checkbox"/>
BS-R-297	Tyre 12" Cross Long Life	£ 42,00	<input type="checkbox"/>
BS-R-301	Tyre 14" City Long Life	default	<input type="checkbox"/>
BS-R-334	Tyre 16" City Long Life	£ 32,00	<input type="checkbox"/>
BS-R-298	Tyre 16" Cross Long Life	£ 42,00	<input type="checkbox"/>

Brakes

SX: BS-R-319 DX: BS-R-318	Single mechanical brake	<input type="checkbox"/> LHS	<input type="checkbox"/> RHS	£ 0,00	<input type="checkbox"/>
SX: BS-O-203 DX: BS-O-273	Single mechanical Parking brake	<input type="checkbox"/> LHS	<input type="checkbox"/> RHS	£ 15,00	<input type="checkbox"/>
SX: BS-O-274 DX: BS-O-275	Single hydraulic brake	<input type="checkbox"/> LHS	<input type="checkbox"/> RHS	£ 135,00	<input type="checkbox"/>

Mudguard

BS-O-245	Mudguard Black	default	<input type="checkbox"/>
BS-O-167	Mudguard Blue	£ 0,00	<input type="checkbox"/>
BS-O-169	Mudguard Orange	£ 0,00	<input type="checkbox"/>
BS-O-166	Mudguard Green	£ 0,00	<input type="checkbox"/>
BS-O-164	Mudguard Red	£ 0,00	<input type="checkbox"/>
BS-O-215	Mudguard Carbon	£ 106,00	<input type="checkbox"/>

Frame colour

BS-O-078	Titanium	default	<input type="checkbox"/>
BS-O-078-P	Polished	£ 280,00	<input type="checkbox"/>
BS-O-217-N	Black Ral 9005	£ 150,00	<input type="checkbox"/>
BS-O-217-R	Red Ral 3020	£ 150,00	<input type="checkbox"/>
BS-O-217-W	White Ral 9003	£ 150,00	<input type="checkbox"/>
BS-O-217-A	Orange Ral 2004	£ 150,00	<input type="checkbox"/>
BS-O-217-B	Blue Ral 5017	£ 150,00	<input type="checkbox"/>
BS-O-217	RAL on request	£ 150,00	<input type="checkbox"/>

Personalized frame wraps

BS-O-087-M	Camouflage	£ 560,00	
BS-O-087-C	Carbon Effect	£ 560,00	
BS-O-087-F	Graffiti wrap Black/White	£ 560,00	
	Personalisation: Two-tone, gradients, etc.	on request	

Handlescrews

BS-O-301	Handlescrews - Black	default	
BS-O-128	Handlescrews - Blue	£ 0,00	
BS-O-127	Handlescrews - Orange	£ 0,00	
BS-O-130	Handlescrews - Green	£ 0,00	
BS-O-129	Handlescrews - Red	£ 0,00	
BS-O-139	Lever handlescrews for Tetraplegic user	£ 0,00	

ACCESSORIES/OPTIONAL

BS-O-083	Central arch stand	£ 165,00	
BS-O-278	Lateral stand with wheels	£ 125,00	
BS-O-229	Rearview mirror	£ 32,00	<input type="checkbox"/> LHS <input type="checkbox"/> RHS
BS-O-213	Triride bag	£ 170,00	
BS-O-114	Arms stopper rings, pair	£ 24,00	
BS-O-232	Light - Front 1200 Lumen	£ 110,00	
BS-O-198	Light Front USB	£ 32,00	
BS-O-180	Light - Rear USB	£ 36,00	

Prices are intended for public and before tax

For details and photos of accessories and spare parts, see the relevant "spare parts and accessories" price list or visit our website [web www.trirideitalia.com](http://www.trirideitalia.com)

N. B.: The Company reserves the exclusive right to bring changes, without notice, if it deems it useful for the improvement of the Product