

ENVOY RANGE

Scoters

ENVOY

“The Envoy brings a new stylish design and offers a huge 30 mile range”



CE

The Envoy range of Scooters are excellent additions to the Drive DeVilbiss range, with the performance usually associated with much larger models. These scooters enjoy a whole host of excellent features including all-round suspension, LED headlights and longer range batteries.

www.drivedevilbiss.co.uk

drive
DeVilbiss®
HEALTHCARE



ENVOY 8+

- › Max range of 30 miles*
- › Top speed of 8mph*
- › Fully functioning indicators and mirror for road use
- › Powerful 500W motor
- › Secondary cable brake
- › Sport style seat
- › Modern dashboard design
- › Black pneumatic tyres with large 5" rims
- › Hydraulic tiller adjustment
- › Delta bars as standard
- › 50ah batteries

CODE: MS058

DESCRIPTION: Envoy 8+



Length:
121cm (48")



Height:
122cm (48")



Width:
60cm (24")



Max Capacity:
160kg (25st)



Weight:
100kg (220lbs)



Max Gradient:
8°



ENVOY 6

- › Max range of 30 miles*
- › Top speed of 6mph*
- › Fully functioning indicators and mirror for road use
- › All round suspension
- › Powerful 470W motor
- › Stylish silver wheel hubs
- › Delta bars as standard
- › Black pneumatic tyres with large 5" rims
- › LED front headlight
- › 50ah batteries

CODE: MS051

DESCRIPTION: Envoy 6



Length:
121cm (48")



Height:
102cm (40")



Width:
60cm (24")



Max Capacity:
160kg (25st)



Weight:
94kg (207lbs)



Max Gradient:
8°



ENVOY 4

- › Max range of 30 miles*
- › Top speed of 4mph*
- › Powerful 350W motor
- › All round suspension
- › Stylish silver wheel hubs
- › Delta bars as standard
- › Black pneumatic tyres with large 5" rims
- › LED front headlight
- › 50ah batteries

CODE: MS050

DESCRIPTION: Envoy 4



Length:
121cm (48")



Height:
102cm (40")



Width:
60cm (24")



Max Capacity:
160kg (25st)



Weight:
94kg (207lbs)



Max Gradient:
8°


Colour Options:



Red



Blue



White (4 & 8 Only)



Silver (4 & 8 Only)

Leaflet Code: LL157 If you would like to view this product please contact your local stockist, details across.

We have a continuous product improvement policy and consequently reserve the right to amend design and specification without prior notice. All sizes and weights stated are nominal. All errors and omissions excepted. Due to the limitations of photography and the printing process colours shown may not be 100% accurate.* Depends on user mass, terrain and battery condition.