



Carbon*Lite*



CarbonLite

The CarbonLite emerges as a revolutionary solution that combines cutting-edge carbon fiber technology with exceptional portability and convenience.

New and Exclusive:

Innovative design and state-of-the-art materials, this scooter sets itself apart from the crowd, catering to the needs of individuals seeking a lightweight, durable, and feature-rich option.

Electromagnetic Brake

Safer, Accurate Breaking

20 Layers of Carbon
Film Wrapping
1/3 Lighter



Only 19kg!



Bold Red
Classic Black
Deep Blue





The Carbon*Lite* folding mobility scooter represents a game-changing innovation in the world of personal mobility. With its lightweight carbon fiber construction, impressive weight capacity, extended range and a host of other convenient features.

Rediscover the joy of independence and freedom with Carbon*Lite*.



Removeable Lithium Batteries



Powered by lithium batteries the CarbonLite delivers optimal performance and reliability. Lithium batteries are known for their high energy density and prolonged lifespan, offering a peace of mind. With a range of 15 miles on a single charge, the CarbonLite empowers users with extended mobility and independence.

Full Specification

Total Weight - 19kg

Max Range - 12 Miles

Maximum Speed - 4mph (6kmh)

Weight Capacity - 22st (140kg)

Total Weight - 19kg

Folded Dimensions - 400 x 480 x 680mm

Dimensions - 890 x 480 x 880mm

Material - Carbon Fibre

Max Gradient - 8°

Scooter Ground Clearance - 38 mm (1.5")

Batteries Required - 24V 6AH

Removable Battery - Yes

Airline Friendly Batteries - Yes

Tyre Type - Solid

Colour - Black, Red, Blue

Airline Approved Battery

020 8159 7070



sales@lifestyleandmobility.co.uk

www.lifestyleandmobility.co.uk