

Airium

25 mi.
Range

13 kg
Heaviest Part

21.4 st.
Weight Limit

The Airium offers the innovation of Lithium battery technology to a lightweight aluminium chassis, making it the easiest to lift and largest travel scooter in our range.



Lithium-ion
Battery Technology



Rear monoshock suspension



Large wheels with super soft tyres



Front LED light



Comfy rotating seat



Wide floor pan



Soft grip delta handlebars



Easily disassembles

Available in:



Plum



Teal

+44 (0) 844 257 8150
motionhealthcare.co.uk

Motion

Specifications

Range	25 miles
Heaviest Part	13 kg
Top Speed	4 mph
Seat Width	46 cm
Wheels	22.5 cm
Battery Watt Hours	793.8 Wh
Total Weight	47.31 kg
Battery Weight	7 kg
Weight Limit	21.4 stone
Battery	31.5ah
Motor	250w
Ground Clearance	5.2 cm
Max. Gradient	8°
Turning Radius	145 cm

