

LENGTH	112cm / 44"
WIDTH	62cm / 24.4"
HEIGHT	90.5cm / 35.6"
FOLDED LENGTH	N/A
FOLDED WIDTH	N/A
FOLDED HEIGHT	N/A
SEAT WIDTH	46cm / 18"
FRONT WHEEL	Solid 22.5cm / 9"
REAR WHEEL	Solid 22.5cm / 9"
TOTAL WEIGHT	93.91lbs / 42.6kg
HEAVIEST PART	29.9lbs / 13.58kg
BATTERY WEIGHT	15.4lbs / 7kg
MAX USER WEIGHT	21.4stone / 136kg
RANGE	25miles / 39km
SPEED	4mph / 6.4kph
BATTERY	25.2v 30.15ah Lithium
MOTOR	24v 250w
GROUND CLEARANCE	5.2cm / 2"
MAX GRADIENT	8°
TURNING RADIUS	145cm / 57.1"
BOX LENGTH	130cm / 51"
BOX WIDTH	55cm / 21"
BOX HEIGHT	53cm / 20"
PACKAGED WEIGHT	108.3lbs / 49.12kg