TRAVELUX TIEMPO Specification Sheet

Removable basket

360 degree swivel chair with adjustable seat height

Angle adjustable Tiller

Black puncture proof PU tyres





level and diagnostic



armrests for easier

41 cm

40 cm

50 - 55 cm*



Convenient flip-back Clever hidden storage

Available Colours



Full Specifications

Seat Width Seat Depth Seat Height (From Ground) **Total Length**

102 cm **Total Width** 51 cm **Total Height** 90 cm

Maximum Speed Range **Turning Radius Ground Clearance** Battery Charger **Total Weight**

(Including Battery) Maximum User Weight 4 mph

19 km (12 miles)** 115 cm 8 cm 2 x 12V / 12Ah 24V / 2Ah

48 kg

113 kg (17.7 Stone)

