# VelumiLi

Lithium powered scooters









### Rascal **VelumiLi** Profile:



The **Rascal VelumiLi** is a lightweight, Aluminium boot scooter with high quality State-Of-The-Art **Lithium Ferro-Phosphate Battery**. At just 28.8kg\*, our lightest Rascal ever also includes a 2 Year Warranty for peace of mind. Brimming with features the slimline frame is just 18.6 "wide yet feels sturdy with 8" dual compound tyres providing a smooth ride, comfortable & spacious seating position for travelling & stability for everyday adventures. Discover your style with the choice of 5 vibrant colours.

Featuring two easy access charging points, on the tiller or on the removable battery box, which is ideal for indoor charging. **The compact & clever styling with super light front section weighing just 10.2kg** makes it exceptional to manoeuvre, with 1150mm turning radius perfect for travelling near or far, at home or overseas.

Enjoy the convenience & impressive durability...



...Today is the day to start your journey!

With the **Rascal VelumiLi** every detail makes a statement, from the crafted white stitching on the seat, LED front light to be seen with, reflective side strip & a good size front basket which is quick to attach on or off for bitesize essentials. This 4mph lightweight scooter has a range of upto 11.6 miles, ideal for both indoor & outdoor use, it is slick & easy to dismantle, with a T Bar tiller locking feature for efficient transportation. It may be small but it is mighty & great value for money.

Relax & enjoy the ride in the knowledge that the **Lithium Ferro-Phosphate** (LiFePO4) Battery pack is warrantied for a design life of over 1000 cycles providing an increased efficiency & performance to stand up to everyday life.

\*Weight excluding Battery

Your perfect travel companion

#### Main Features:

- Lightweight, aluminium, high quality rugged frame
- British Made / Medical / Military grade Lithium Ferro-Phosphate Battery
- Airline Compliant 12.4Ah battery weighs just 4.1kgs / 9 lbs
- 1150mm turning radius for exceptional manoeuvrability
- Tiller locking for easy transportation
- Fold down backrest & simple to remove adjustable armrests
- Front basket has good size & depth
- Two easy access charging points on tiller & battery box
- Dual compound puncture proof tyres for impressive stability
- 15.7" seat height that offers real comfort & support for long journeys
- Folded Dimensions (LxWxH) 106cm x 47cm x 71cm / 41.7" x 18.6" x 27.9"



Clever, practical front handle increases accessibility



#### **Boot Transportable**






#### Air travel & CAA Lithium Battery Legislation

Airline Compliant 12.4Ah battery meets I.A.T.A. UN38.3 Passenger Aircraft Transportation requirements for overseas travel & can be taken on a plane.

(Subject to Airline Regulation & Restrictions)



Park to pavement, car to countryside, the VelumiLi is super lightweight & compact, simple to store or transport. Perfect for travels near or far!

## VelumiLi



	contraction of the contraction o	
Specification:	VelumiLi 12.4Ah	
Type Class	B/2	
Overall length	106cm / 41.7"	
Overall width	47cm / 18.6"	
Overall height	88cm/34.7"	
Floor to seat height	39 - 43cm (min-max)	
Seat width	41cm / 16.1"	
Seat depth	37cm / 14.6"	
Folded Dimensions LxWxH	106cm / 47cm / 71cm (41.7" / 18.6" / 27.9")	
Maximum carrying capacity	115kgs / 18 stone	
Battery voltage and capacity	25.6V / 6.2Ah x 2	
Maximum speed	6.4km/h / 4mph	
Ground clearance	4.5cm / 1.8"	
Maximum kerb height	5cm / 2"	
Range*	up to 18.7km / 11.6 miles	**up to 14.2km / 8.8 miles
The maximum safe slope	6°	
Turn-around width	1.25m / 4'1"	
Colour	Blue Moon / Green Lightning / Orange Sunset / Pink Thunder / Yellow Star	
Single Heaviest Component (Front Section)	10.2kgs / 22.5 lbs	

 $<sup>{}^*</sup>Range\ on\ full\ charge\ and\ flat\ ground\ -\ with\ 68kg\ User\ Weight \qquad {}^{**}Range\ with\ max\ 115kg\ User\ w$ 

### Today is the day to begin your adventure

28.8kgs / 63.5 lbs (Excluding battery)

Local Rascal Retailer:

Product weight



32.8kgs / 72.3 lbs (Including battery)





