



APEX
ALUMA-LITE
PLUS



Colours Available:



Lightweight Aluminium Travel Scooter

One of the lightest aluminum travel scooter on the market! Experience the convenience of a lightweight, compact travel scooter with airline safe batteries. Easy to put into the boot, with its heaviest piece weighing only 11.6 kg. The Apex Aluma-lite plus, with its aluminum frame, weighs only 33.2kg making it one of our lightest four wheeled scooter yet!

Specifications	
Scooter	
Weight Capacity	21 st 6 lbs (136 kg)
Overall Length	1010 mm (39.76")
Overall Width	495 mm (19.48")
Overall Height	840 mm (33")
Drive Tyres	Front & Rear - 203 mm (8") Solid
Ground Clearance	32 mm (1.25") at Motor
Heaviest Piece When Disassembled	11.6 kg - Rear Main Frame
Colour(s)	Red, Blue & Black
Mobility	
Turning Radius	1385 mm (54.5")
Max Speed	Up to 4 mph
Max. Obstacle Climb	57.5 mm (2.25")
Drivetrain	Rear-Wheel Drive Single Motor 24v DC Motor
Dual Breaking System	Regenerative & Electromechanical
Maximum Obstacle	57.5 mm (2.25")
Battery	
Battery Requirements	12 Ah Lithium Iron Phosphate 12 Ah (Standard) 30 Ah (Optional)
Battery Charger	Off-Board
Total Weight Without Batteries	33.6 kg
Battery Weight	12 Ah Battery Pack 4.25 kg Each 30 Ah Battery Pack 8.63 kg Each

Features



Suspension



Easy to Use Controls



Airline Safe Battery