



Colours Available:



Sporty Style with Advanced Comfort

The all new Colt Sport is fully loaded with sporty style and advanced features, including a DX Tiller, LED headlights and an advanced suspension seat post as well as an all new memory foam adjustable seat.

Drive at speeds up to 8 mph and see what life has in store for you!

Specifications	
Scooter	
Weight Capacity	25 st (159 kg)
Overall Length	1194 mm (47")
Overall Width	560 mm (22")
Drive Tyres	254 mm Front & Rear - (10") Solid or Pneumatic
Ground Clearance	57 mm (2.25")
Heaviest Piece When Disassembled	29.75 kg (65.5lbs)
Colour(s)	Matte Black and Stone Grey
Suspension	Front & Rear Suspension
Mobility	
Turning Radius	1370 mm (54")
Max Speed	Up to 8 mph (12.9 km/h)
Max. Obstacle Climb	50 mm (2")
Drivetrain	24v DC Motor
Dual Breaking System	Regenerative & Electromechanical
Battery	
Battery Requirements	2 x 32 Ah Deep Cycle. (optional 40Ah)
Battery Weight	10 kg (40Ah 14.75kg)
Battery Charger	Off-Board
Total Weight Without Batteries	64.5 kg (142.5lbs)

Features



Suspension Seat Post



10" Alloy Style Pneumatic or solid tyres (solid tyres only available in black)



Memory foam adjustable seat



DX Tiller and LED Lighting