

# EXCEL GALAXY II Specification Sheet

Height adjustable headrest

360 degree rotating chair with depth and height adjustable Captain seat

Height adjustable and foldable armrests

Angle adjustable Tiller

Removable basket with carry handle

Spacious footwell for optimum seating posture



New multifunction LCD display dashboard



Delta Tiller handle for easier steering



Rear LED lighting for improved visibility and extra safety



Sleek modern shroud design with alloy wheels

Available Colours



## Full Specifications

Seat Width	51 cm	Maximum Speed	8 mph
Seat Depth	46 cm	Range	56 km (35 miles)*
Seat Height (From Ground)	65 cm	Turning Radius	161 cm
Total Length	146 cm	Ground Clearance	13 cm
Total Width	65 cm	Battery	2 x 12V / 90Ah
Total Height	130 cm	Charger	24V / 8Ah
		Total Weight (Including Battery)	127 kg
		Maximum User Weight	160 kg (25 Stone)

KEY: \* = Varies depending on weather and terrain conditions



[www.vanosmedical.co.uk](http://www.vanosmedical.co.uk)