

AIRSCAPE

10 mi.
Range

14.5 kg
Heaviest Part

21 st.
Weight Limit

The Airscape offers the innovation of Lithium battery technology to a durable steel chassis, making it the most affordable Lithium scooter in our range.



Lithium-ion
Battery Technology




Quick release T-bar handle



Super soft tyres




Air safe Lithium battery



Rotating padded seat



Long floor pan



Soft grip & dual controls



Easily disassembles

Available in:



Black

+44 (0) 844 257 8150
motionhealthcare.co.uk

Motion

Specifications

Range	10 miles
Heaviest Part	14.5 kg
Top Speed	4 mph
Seat Width	43 cm
Wheels	19 cm
Battery Watt Hours	240 Wh
Total Weight	35 kg
Battery Weight	1.7 kg
Weight Limit	21 stone
Battery	10ah
Motor	270w
Ground Clearance	5 cm
Max. Gradient	6°
Turning Radius	125 cm

Accessories



Delta handlebars



In-car charger

