



Colours Available:



High Performance With Exclusive Features


The sleek, sporty Colt Twin scooter has unique dual front wheels offering the stability of a four-wheel scooter. The Colt Twin has advanced standard features such as low-profile tyres and a wraparound delta tiller.

Specifications	
Scooter	
Weight Capacity	25 st (159 kg)
Overall Length	1160 mm (46")
Overall Width	560 mm (22")
Drive Tyres	Front: 230 mm (9") Solid Rear: 250 mm (10") Solid
Ground Clearance	57 mm (2.25") at Motor
Heaviest Piece When Disassembled	27.25 kg (front section)
Colour(s)	Red, Blue
Standard Seating	Type CRS Junior Weight: 19 kg (42 lbs) Material: Black Vinyl Dimensions: 445 mm (17.5") Overall Width: 432 mm (17") Effective Depth
Mobility	
Turning Radius	1067 mm (42")
Max Speed	4 mph (6.4km/h)
Max. Obstacle Climb	50mm
Drivetrain	Rear Wheel Drive, Sealed Transaxle, 24v DC Motor
Dual Breaking System	Regenerative & Electromechanical
Battery	
Battery Requirements	2 x 30 Ah (standard)
Battery Weight	10 kg
Battery Charger	Off-Board
Total Weight Without Batteries	64.5 kg

Features



Large Amount of Leg Room



Comfortable High Back Seat With Armrests



Exclusive & Stylish Low Profile Tyres



REV 02 - 31/01/2023

