

Permobil

# Panthera X



permobil



# Strong and light-weight for maximum maneuverability

Designed for experienced manual wheelchair users who need a dynamic mobility device to fit into their busy, active lifestyles, the Panthera X is one of the world's lightest wheelchairs at 2.1 kg.

## Key benefits:

- Ultra-light weight
- Extremely strong chassis
- Extremely maneuverable
- Easy to transport
- Adjusts for comfort



## Specifications

Seat width	33 - 36 - 39 - 42 - 45 cm
Seat depth	35 - 46 cm
Seat height	Front 47 cm (w. 24" rear wheels) Back 43 cm (w. 24" rear wheels)
Seat angle	7°
Total length	82.5 cm
Total width	Seat width +20.5 cm
Front frame angle	100°
Backrest height	22 - 35 cm
Backrest angle	11.5°- (-7.5°)
Footplate distance	35 - 43 cm
Camber	2.2°
Weight without rear wheels	2.1 - 2.3 kg
Wheelchair weight	4.4 - 4.6 kg
Maximum load	100 kg
Material (chassis)	Carbon fibre
Material (upholstery)	Polyurethane covered polyester

Permobil reserves the right to change specifications and design without prior notice.

Personalise your ride

CARBON  
FIBRE

Actual product color may vary.