

Adjusted by:

Tel:

Save the form to the desktop and open it with Adobe Reader.

E-mail:

Date:

The price changes depending on the choices made in the form £:

Standard configuration includes; Chassis, backrest, castors, rearwheel Spinergy 24", alu pushrim, Schwalbe One tire, one hand X brake, Carbon sideguards.

Choice required: Max. user weights : Seat width 33-45 cm 100 kg,

Optional  Upcharge /Add-price  Marked yellow Click on "blue underlined text" for more info or image



**Original design**

**Alternative options**

**Other options**

**1) Panthera X**

Panthera X



**2) Chassis - Seat width** [Art.no.](#)

33 cm  36 cm  39 cm  42 cm  45 cm

Clear



**3) Backrest height** 22-35 cm (HH) [Art.no.](#)

22 cm  23 cm  24 cm  25 cm  26 cm  27 cm  28 cm

29 cm  30 cm  31 cm  32 cm  33 cm  34 cm  35 cm

Clear



**4) Balance** [Info](#)

**Position 1**  
Most easy balanced

**Position 2**  
Less easy balanced

**Position 3**  
Medium balanced

**Position 4**  
Least easy balanced

The positioning of the backrest on the chassis affects the balance / COG of the wheelchair. "easy balanced" = easy to practice rear wheel balance. This can later easily be changed.



**5) Castor** [Art.no.](#)

**3" X** (87 mm)

[90 mm S3](#)

Clear



**6) Footrest**

[Footrest Alu X](#)  
36004BB

[Footrest Titanium](#)  
36915BB

[Footrest Titanium](#), extended. reinforced  
36006BB [needed for Freewheel](#)

[Carbon Footplate](#)  
65060BB

[Plastic plate for footbow](#)  
36010BB

Clear



**7) Brake** [Art.no.](#)

[One hand brake right](#)

[One hand brake left](#)

Clear



### 8) Rear wheel - without pushrim

[Rear wheel sizes](#)

- [24" Spinergy X](#) black spokes
- [24" Spinergy X](#) white spokes
- Clear
- [25" Spinergy X](#) black spokes
- [25" Spinergy X](#) white spokes
- [26" Spinergy X](#) black spokes
- [26" Spinergy X](#) white spokes



### 9) Tyre - choose tyre

The price is for pair

- [Schwalbe One](#) 24", 25", 26"
- Clear
- [Marathon Plus 24", 25", 26"](#) £68
- [Right Run 24", 25", 26"](#) £64
- [Rugged 24"](#) £68
- [PU 24" 8 massive](#) £110
- reTyre does not fit PU tyres
- [reTyre skins 24"](#) £245
- [reTyre skins 25"](#) £245
- [reTyre skins 26"](#) £245
- Clear



### 10) Pushrim

The price is for pair

- [Aluminum X](#) 1904404 x2 (04=24", 05=25", 06=26")
  - Clear
  - [Titanium](#) 190DD11 (x 2)
  - [Paragrip](#) 19051DDBT x 2 £150
  - [Tetragrip](#) 19052DDBT x 2 £150
  - [Maxgrip 22", 24"](#) 1903060 x 2 (60=24", 80=22" (contains latex)) £330
- DD= Wheeldiameter in inches

### Accessories



### Side guard Carbon [art.no.](#)

- [Carbon without fender](#)
- [Carbon small fender](#)
- [Carbon large fender](#)
- Clear



£230


### Side guard S3 [video](#)

- [Side guard S3](#)



### Cushion Wedge cushion

- [5 cm](#) 572BB40
- [Wedge cushion](#) 57710BB



£105

£49



### Arm hook

- [Arm hook pair](#) 4200010

£85



### Transport protection

- [Housing](#) 9024009 £92
- [Rear end](#) 9024012 £92
- [Rear axle](#) 9024008 £92



### Other

- [Tube protection, front](#) 9900252 (2 pcs) Tube £15
- [protection, rear](#) 9900252 (1 pcs) £15
- [Freewheel](#) £715

Height 14-19 cm, art: 2600106

Chose art.no. for Freewheel according to the footrest height of the wheelchair  
Read more about [Freewheel here](#)

### Set-up adjustments [Backrest angle](#)



Note measure (A) for desired back rest angle 10-20 mm. Standard = 15 mm

A=  mm

### [Footrest height](#)



Note measure (A) for desired footrest height

A=  mm

### Extra