



Feather **FOLD**

The lightweight folding transportable mobility scooter



Feather Fold **Technical Specifications**



Class 2 Mobility Scooter
Class 2 scooters can be used on pavements.



Batteries
10Ah Lithium



Overall Size
Height: 914mm / 36 inches
Length: 1020mm / 40 inches
Width: 559mm / 22 inches



Speed Mode
4mph



Braking System
Electromagnetic



Folded Flat Size
Height: 457mm / 18 inches
Length: 1953mm / 37.5 inches
Width: 559mm / 22 inches



Range
Up to 7.5 miles*
Range may vary with battery charge/condition/option, terrain type, tyre condition & vehicle loading.



Weight Capacity
120kg / 19 stone



Tyres - Solid
Front - 152mm / 6"
Rear - 254mm / 10"



Scooter Weight
16.8kg
Excluding batteries.



Ground Clearance
60mm / 2.4 inches (unloaded)



Maximum Slope
8 degrees



Charge Time
6 - 7 hours



Seat Width & Seat Depth
380mm / 15 inches
356mm / 14 inches