




# Feather PROPEL

The lightest mobility range in the UK



## Feather Propel Technical Specifications



**Wheelchair Weight**

7.2kg  
Excluding footrests.



**Armrest Height**

229mm / 9 inches



**Overall Size**

Height: 889mm / 35 inches  
Length: 762mm / 30 inches  
Width: 635mm / 25 inches



**Frame Weight**

6.1kg




**Brakes**

Attendant handles  
& push to lock



**Folded Flat Size**

Height: 381mm / 15 inches  
Length: 737mm / 29 inches  
Width: 711mm / 28 inches



**Backrest Height**

406mm / 16 inches



**Weight Capacity**

113kg / 18 stone



**Tyres - Solid**

Front - 152mm / 6 inches  
Rear - 559mm / 22 inches



**Seat-to-Floor Height**

483mm / 19 inches



**Seat Width & Seat Depth**

457mm / 18 inches  
406mm / 16 inches



**Footrest**

Swing-away