



LuXe Mobility X500 Midi Sport: Compact, Powerful, and Designed to Elevate Your Freedom



Lifestyle & Mobility

020 8159 7070

Specifications

Overall Width **525 mm**
Overall Length **1110 mm**
Turning Radius **1500 mm**
Ground Clearance **45 mm**
Total Weight **50KG**
Max Speed **4mph**
Max Range **14 Mile**
Weight Capacity **21.5st**
Colours **Blue, Red, Purple, Silver, Green**



21.5 St User Weight

With a maximum user weight of 21.5 stone, this boot scooter offers exceptional strength without compromising on portability or convenience.

Rear Suspension

Equipped with rear suspension and pneumatic tyres, the LuXe Mobility X500 Midi Sport ensures a smoother, more comfortable ride across pavements and uneven surfaces alike.



Miles of Freedom

Impressive Range Travel up to 14 miles on a single charge with reliable, long-lasting battery performance.