

Extreme X8



Rev4.:10/2023


ORDER FORM

Errors and imissions excepted..

Billing Address Or Customer-ID			Shipping Address		Identical To Billing Address
C.-ID:			C.-ID:		
Name:			Name:		
Street:			Street:		
Town/PC:			Town/PC:		
E-Mail:			E-Mail:		
Phone:			Phone:		
Order Information			User Weight		
Order	Calculation	Built-4-Me®	Art. No.:		
Purchase Order No.:			User weight: kg		
Reimbursement Code:					
Order date: Day Month Year			Base Price £15,950		

Client Measurements

a) Top Of Head To Seat	mm
b) Top Of Shoulder To Seat	mm
c) Underarm To Seat	mm
d) Elbow To Seat	mm
e) Top Of Pelvis To Seat	mm
f) Shoulder Width	mm
g) Trunk Width	mm
i) Hip Width	mm
j) Actual Sitting Depth	mm
m) Foot Length	mm
n) Lower Leg Length	mm
o) Actual Seat Width	mm



Technical Data

User Weight Capacity	155 kg with seat lift and seat inclination	Ground Clearance	100 mm (depends on the configuration)
Base Width	710 mm	Drive Wheels & Castors	4 x 355 mm (14") Drive Wheels
Chassis length	1030 mm	Charger	220 V/10 A
Drive Control	R-Net (120A Power Module)	Battery Size	80 Ah (108 Ah optional)
Drive Motors	10 km/h, 4 x 700 W, 4-Pole Motors	Turning radius	1415 - 1470 mm
Standard configuration	Including MPS back system, MPS seat system, maintenance-free 80 Ah batteries and 10 A charger		
Tyres	Maximum tyre pressure depending on user weight from 0.25 to 0.55 bar (see table in instruction manual)		

Drive Unit

Power Wheelchair

		£
MAG010006	EXTREME X8	15,950

Seating Systems

		£
MAGREHABSEATPAN	Rehab Seat Pan	0

Seat Cushion For Rehab Seat

		£
MAGWOSEATCUSH	No seat cushion	0
MAGCMRC	MPS Removable Seat Cushion (1) (2)	0
MAGPADDEDSC	Flat padded seat cushion	235

UPH For Seat Cushion

		£
MAGSUFABRIC	Fabric	0
MAGSUVINYL	Vinyl	0

Seat Sizes

Seat width ⁽³⁾

MAGSW305MM	12" (305mm)	MAGSW405MM	16" (405mm)	MAGSW510MM	20" (510mm)	MAGSW610MM	24" (610 mm)
MAGSW330MM	13" (330mm)	MAGSW430MM	17" (430mm)	MAGSW535MM	21" (535mm)		
MAGSW355MM	14" (355mm)	MAGSW455MM	18" (455mm)	MAGSW560MM	22" (560mm)		
MAGSW380MM	15" (380mm)	MAGSW480MM	19" (480mm)	MAGSW585MM	23" (585mm)		

Seat depth ⁽³⁾

MAGSD405MM	16" (405mm)	MAGSD480MM	19" (480mm)	MAGSD560MM	22" (560mm)
MAGSD430MM	17" (430mm)	MAGSD510MM	20" (510mm)	MAGSD585MM	23" (585mm)
MAGSD455MM	18" (455mm)	MAGSD535MM	21" (535mm)	MAGSD610MM	24" (610 mm)

Seat Tilt Modules

		£
MAG2852-KIT	Fixed Frame, 3° tilt (no power functions)	0
MAG3502-KIT	Fixed Frame, 3° tilt (only lighting)	0
MAG5335	Fixed Frame, 3° tilt (with Power Options)	0
MAG5127	Powered Tilt (0 - 50°)	1950
MAG5126	Lift and Tilt (300 mm lift, 0 - 45° tilt)	2850

Seat Height

		£
MAG6552	19" (470mm)	0
MAG3323-PI-2	20" (495mm)	0
MAG3323-PI-4	21" (520mm)	0

(1) Minimum Width And Depth 405 mm

(2) Seat cushion sizes only available according to Section A

(3) Depth minus Width Must Not Exceed 50 mm

BACKREST**Backrest for Rehab Seat Pan**

		£
MAGREHABCANES	Rehab Canes	0
MAGREHABFOLD	Rehab Canes - Fold Forward	375

Recline Option

		£
--	--	---

Recline Option for Rehab seat and Rehab canes

MAG3920-S-3921-F	Fixed (90° - 105°), (90°, 97,5°, 105°)	0
MAG3924-S-3925-F	Manual Recline (90° - 170°)	395
MAG3928-S-3929-F	Power Recline (90° - 160°)	1175
MAG3932-S-3933-F	Anti-shear power recline (95° -155°)	1875

Recline option for Rehab seat and Rehab canes, Foldable to the front

MAG3922-S-3923-F	Fixed (90° - 15°), (90°, 97,5°, 105°)	0
MAG3926-S-3927-F	Manual Recline (90° - 170°)	395
MAG3930-S-3931-F	Power Recline (90° - 160°)	1175
MAG3934-S-3935-F	Anti-shear power recline (95° -155°)	1875

Insert Back Clamps

		£
MAG3492-PI	Quick Release	0

MPS Insert Back

		£
MAGMPSINSERT	MPS Insert back	0

MPS Insert Back UPH

		£
MAGFABRICBR	Fabric	0
MAGVINYLBR	Vinyl	0

MPS Insert Back – Sizes

		£
--	--	---

Backrest width (1)(2)

MAGBW405MM	16" (405mm)	0
MAGBW455MM	18" (455mm)	0
MAGBW510MM	20" (510mm)	0

Backrest height (1)

MAGBH535MM	21" (535mm)	0
MAGBH585MM	22" (585mm)	0
MAGBH635MM	23" (635mm)	0

Rehab Insert Back

		£
MAGBACKRESTINSERT	Magic Rehab insert back (1)	535

Rehab Back Insert UPH

		£
MAGFABRICBR	Fabric	0
MAGVINYLBR	Vinyl	0

Magic Rehab Insert Backrest Cushion

		£
MAGBRMRJ	Small (350 mm height)	0
MAGBRMRJXS	Small (410 mm height)	0
MAGBRMRMS	Medium (430 mm height)	0
MAGBRMRM	Medium (480 mm height)	0
MAGBRMRMXT	Medium (530 mm height)	0
MAGBRMRLS	Large (430 mm height)	0
MAGBRMRL	Large (480 mm height)	0
MAGBRMRLXT	Large (530 mm height)	0

(1) Backrest Height Minus Backrest Width Must Exceed 50 mm

(2) Back width has to match seat width selected

Price List - QPA01
REV2 (01/01/2024)

Rehab Back Insert Bolster Size			£
For Magic Rehab back cushion large/medium			
MAGBOLSTERS	Shallow (110 mm depth)		55
MAGBOLSTERSTD	Standard (160 mm depth)		80
MAGBOLSTERD	Deep (220 mm depth)		105
MAGNOBOLSTER	No bolster		0

Rehab Back Insert Lumbar			£
MAGLUMBARS	Shallow (12 mm insert)		0
MAGLUMBARSTD	Standard (25 mm insert)		0
MAGLUMBARD	Deep (35 mm insert)		0
MAGNOLUMBAR	No lumbar		0

Push Handles			£
MAG46354636	Rehab peg handles		135

Backrest Bag Hooks			£
MAG5723	Rehab Bag hooks (J-shaped) on spreader bar		35
MAG5725	Rehab Bag hooks on Push Rail		55

HEADREST

Headrest - All Types			£
MPS Headrest			
MAG2770VI	MPS Headrest, MPS Mount		0
Rehab Headrest			
MAGREHABSTD	Magic Rehab Standard		160
MAGREHABWIDE	Magic Rehab Wide		220
Headrest Mount (for REHAB headrests only)			
MAG4113	2 Axis		165
MAG4114	3 Axis		190
MAG4115	Multi Axis		215

Rehab Headrest			£
MAGFABRICHR	Fabric		0
MAGVINYLHR	Vinyl		0

ARMREST

Armrest – Types			£
MAG3572	Armrest, Single Post		0
MAG3994	Armrest, Double Post		160
MAGFLEXIARM	Armrest, Flip Back		575
MAG6503-6502	Armrest, Flip Back, with Locked Position		550
MAG3331-6500	Armrest, Flip Back, with Locking Post		600
MAG3727-6501	Armrest, Flip Back, with Locking Post and Drop Forward		700

Armrest – Heights											
MAGAH220MM	220 mm	MAGAH260MM	260 mm	MAGAH300MM	300 mm	MAGAH340MM	340 mm	MAGAH380MM	380 mm		
MAGAH230MM	230 mm	MAGAH270MM	270 mm	MAGAH310MM	310 mm	MAGAH350MM	350 mm	MAGAH390MM	390 mm		
MAGAH240MM	240 mm	MAGAH280MM	280 mm	MAGAH320MM	320 mm	MAGAH360MM	360 mm				
MAGAH250MM	250 mm	MAGAH290MM	290 mm	MAGAH330MM	330 mm	MAGAH370MM	370 mm				

Armpad – Style			£
MAGAP75-WOKOR-2	75 mm armpads		0
MAG3971-3972	MPS armpads		80
MAG3712	MPS flip-out armpads		320
MAGAP50-2	50 mm armpads		55
MAG6711	Waterfall armpads		160
MAG3775	Gutter armpads		160
MAG6750	Gutter armpads with elbow support pad		160
MAG6751	Gutter armpads with hand pad		160
MAG6752	Gutter armpads with hand pad and elbow support pad		160

Armpad – UPH			£
MAGAPVINYL	Vinyl		0
MAGAPFABRIC	Fabric		0

Armrest Bag Hooks			£
MAGHOOKA2RH	Right Side Armrest bag hooks		35
MAGHOOKA2LH	Left Side Armrest bag hooks		35

POSITIONING

Lap Belt							£
MAG4973-MPS-4974-R	Standard belt	500 mm	600 mm	700 mm	800 mm		0
MAG4973-4974-PBPAD	Standard belt with 2 posture pads	500 mm	600 mm	700 mm	800 mm		100
MAG6540	4 Point Belt	500 mm	600 mm	700 mm	800 mm		180
MAG6540-2XPBPAD-4	4 Point Belt with 2 posture pads	500 mm	600 mm	700 mm	800 mm		185
MAG4730	Retractable Belt						110
	MAGLBLH			Left sided buckle			0
	MAGLBRH			Right sided buckle			0

Harness			£
MAGHARP	Harness standard		160
MAGHARQ	Harness pediatric		160

Thoracic Support (pad and mount)			£
MAG4943	MPS Adjustable Laterals		530

Hip Pads - In Vinyl			£
MAGHPKIT	Adjustable Hip Pads		270

LEGREST

Legrest – Choice		£
MAGCPLEGREST	Fixed Centre Mount Legrest	0
MAG3528	Powered Centre Mount Legrest	1600
MAG3527-ME	Standard Swing away Legrest	0
MAGCPOWER	Manual Elevating Swing away Legrest	595
MAG3526-ME	Power Elevating Swing away Legrest (1)	1750

Footplate		£
MAG4055FIX	Centre Post: 1 piece fixed footplate	0
MAG4055FIX-2	Centre Post: 2 piece fixed footplate	0
MAG4055	Centre Post: 1 piece flip up footplate	0
MAG4055SFP	Centre Post: 1 piece flip up, SHORT	0
MAG4055-2	Centre Post: 2 piece flip up footplate	0
MAG4055FP-2	Centre Post: 2 piece flip up, SHORT	0
MAG4755	Powered Centre Post: 1 piece flip up	0
MAG5746	Powered Centre Post: 1 piece flip up WIDE	0
MAG4755-2	Powered Centre Post: 2 piece flip up	0
MAGFLIPUP	Swing-Away: 2 piece flip up (1)	0

Footplate Options

MAG4312	Swing Away: Footplate Lock	55
MAG5323	Centre Post: Metal edges	55
MAG4260	Abduction Block	55

Footplate Height

MAGLL300MM	300 mm	MAGLL340MM	340 mm	MAGLL380MM(2)	380 mm	MAGLL420MM (2)	420 mm
MAGLL310MM	310 mm	MAGLL350MM	350 mm	MAGLL390MM(2)	390 mm	MAGLL430MM (2)	430 mm
MAGLL320MM	320 mm	MAGLL360MM	360 mm	MAGLL400MM(2)	400 mm		
MAGLL330MM	330 mm	MAGLL370MM	370 mm	MAGLL410MM(2)	410 mm		

Calf Support

		£
MAGCALFSUPPORT	Calf Support	135

Knee Pads

		£
MAG4247	Knee Pads for Fixed Centre Mount	125
MAGKPCMPE	Knee Pads for Powered Centre Mount	125
MAG4660	Knee Pads for Standard Swing-Away Legrests	125
MAG4104	Knee Pads for Elevating Swing-Away Legrests	125

CONTROLLER**Controller – R-Net**

		£
MAGPG-JSMLEDL	Remote Joystick, R-Net, LED	0
MAGPG-JSM-PI	Remote Joystick, R-Net, Colour LCD	160
MAGPG-JSM2IRL	Remote Joystick, R-Net, Colour LCD Advanced, Infrared	315
MAGPG-JSMIRBT	Remote Joystick, R-Net, Colour LCD Advanced, Infrared & Bluetooth	530

Joystick Position

		£
MAGJSRIGHT	Joystick Right	0
MAGJSLEFT	Joystick Left	0
MAGPGCYPGC2.03725	Interchangeable to Left and Right	215
MAG3329JSM-RH	Joystick Mounted as Attendant Control - Right	315
MAG3329JSM-LH	Joystick Mounted as Attendant Control - Left	315

Swingback Options

		£
MAG3083-KIT-OS	Swingback Outside (Magnets)	0
MAG3083-KIT-IS	Swingback Inside (Magnets)	0
MAG3082-KIT-OS	Swingback Outside (Locking Latch)	125
MAG3082-KIT-IS	Swingback Inside (Locking Latch)	125

Joystick Mounting Position

		£
MAG3713	Height Adjustable Joystick Stem	105
MAG3729SB	Quick Release Stem	0

Joystick Knob Option

		£
MAGKNOBSTD	Standard joystick knob	0
MAGPG-KS80	Stick S80	80
MAGPG-KCHIN	Chin C38	80
MAGPG-KB32	Ball B32	80
MAGPG-KB38	Ball B38	80
MAGPG-KM48	Mushroom M48	80
MAGPG-KTBAR	T-Bar T80	80
MAGPG-KSB70	Soft Ball SB70	80

Dual Control

		£
MAG3329AM-KIT	Dual Control	290

Dual Control Mounting Location

MAGACRH	Mounting - RHS	
MAGACLH	Mounting - LHS	

5 Switch Box

		£
MAG158709	5 Switch Box (Button)	320
MAG158708	5 Switch Box (Toggle)	320

5 Switch Box Location

		£
MAG5430	Same side as joystick	0
MAG5431	Opposite side to joystick	0
MAG5433	Opposite side to joystick + socket	0
MAG5432	As Attendant	300

Omni Control		£
MAG5436	Omni 2 Input Control Module	1170
MAGPG-OMNISCANNER	Omni 2 Input Control Module (Scanning Mode)	1170
Omni Control Mounting Location		£
MAGOMNILH	Left	0
MAGOMNIRH	Right	0
BASE OPTIONS		
Tyre		£
MAG7050	Pneumatic	0
MAG7058	Neverflat for Extreme X8	225
MAG7056	Extreme X8 Puncture Res Kit	600
X8 Drive wheel rim		
MAG6671	Silver rims	0
MAG6670	Black rims	0
LED Lights		£
MAG5227	Omit lights	0
MAG2857	Lighting mounted on seat	0
MAG6530	Lighting mounted on base	0
Colour Selection		£
MAGRED	Flame Red	0
MAGCOWLBL	Ocean Blue	0
MAGBLACK	Stealth Black	0
MAGORANGE	Desert Tangerine	0
MAGYELLOW	Sunrise Yellow	0
MAGMGREEN	Forest Green	0
MAGCOWLFC	Forest Camo	580
MAGCOWLAC	Army Camo	580
MAGCOWLFL	Flames Fury	580
Batteries		£
MAG6668	Battery, 70 AH	0
MAG6669	Battery, 90 AH	275
Chargers		£
MAG160104	Charger, 10ah	0
Extreme X8 Base Options		£
MAG2899	Crash tested tie downs	0
MAG5420	X8 lifting rings	100
MAG2944	Steering lock	165
MAG4384	Extreme X8 stability roller	165

Accessories		£
MAG3145	Luggage Rack	320
MAG4471	RAM® Cup holder	55
MAG4627	Magic Mobility Steel cup holder (1)	105
MAG4628	Magic Mobility Steel cup holder (Right sided) (1)	105
MAG4705	Accessory charger (1x12V + 2xUSB) - side rail mounted, same side as joystick	160
MAG6629	Accessory charger (2xUSB) - side rail mounted, opposite to joystick	160
MAG4801	Accessory charger (2xUSB) - armrest mounted, same side as joystick	160
MAG6594	Accessory charger (2xUSB) - armrest mounted, opposite to joystick	160
MAG5360	Loud horn	55
MAG4834	A4 Accessory bag (mounted on the backrest)	105
MAG4834RH	A4 Armrest bag (mounted on the right armrest)	105
MAG4834LH	A4 Armrest bag (mounted on the left armrest)	105
MAG4834-2	A4 Armrest bags (mounted on the left and right armrest)	215
MAG5320RHI	A5 Armrest bag (mounted internally, right armrest)	80
MAG5320LHI	A5 Armrest bag (mounted internally, left armrest)	80
MAG5320RHE	A5 Armrest bag (mounted externally, right armrest)	80
MAG5320LHE	A5 Armrest bag (mounted externally, left armrest)	80
MAG48345320I	A4 & A5 Armrest bags (A4 joystick side, A5 internal facing)	160
MAG48345320	A4 & A5 Armrest bags (A4 joystick side, A5 external facing)	160
MAGCOVER	Wheelchair cover	0

Clear Tables		£
MAG2675	Tray Table, Slide In	335
MAG2676	Tray Table, Swing Away	335
MAG2677	Tray Table, Swing Away, Small	335
MAGTDS	Tray Table, Swing Away, Double Stem	435

Comments For Adaption



Extreme X8

£

Base Price

Total Of Additional Equipment

Packaging And Transport

Exclusive VAT

VAT

Total Price

Date

Signature

Day

Month

Year



Sunrise Medical GmbH

Kahlbachring 2-4

69254 Malsch/HD

Tel: +49 (0) 7253-9800

Fax: +49 (0) 7253-980111

E-Mail: kundenservice@sunrisemedical.de

www.SunriseMedical.de


Dimension Matrix MPS-SEAT

◀ BACK



Dimension Matrix REHAB-SEAT

◀ BACK



Magic Rehab Insert Size

◀ BACK

MAGIC REHAB INSERT SIZE



Chest Width 'C'	Height 'Y' (mm)
330-470mm	
Large (Xtra tall)	530
Large	480
Large (Short)	430
280-420mm	
Medium (Xtra tall)	530
Medium	480
Medium (Short)	430
260-400mm	
Junior (Xtra tall)	410
Junior	350


- Bolsters can be mounted upside down for increased armpit support
- Bolster heights can be offset up and down, up to 25mm between each side (asymmetric)

© Magic Mobility 2018

Rehab Seat and Back Sizing

◀ BACK

REHAB SEAT AND BACK SIZING



© Magic Mobility 2018

Price List - QPA01
REV2 (01/01/2024)

Colors

◀ BACK



Camo Forest

MAG5527



Desert Tangerine

MAG5523



Flame Red

MAG5522



Flames Fury

MAG5526



Forest Green

MAG5521



Ocean Blue

MAG5528



Stealth Black

MAG5525



Sunrise Yellow

MAG5524

Drive Wheels Rims & Special Tyre Options

◀ BACK

MAG5924 Silver rims
MAG5927 Black rims

MAG5927




MAG5927 Extreme X8 Puncture Res Kit




MAG4384 Extreme X8 stability roller

Seat Tilt Modules

[◀ BACK](#)

Headrest

[◀ BACK](#)

MPS HEADREST

Armrests

◀ BACK

ARMREST TYPE (1)



STANDARD

Note: 25mm of extra adjustment is added to bottom

DOUBLE STEM
For bariatric users


Note: Depending on requested armrest height the flexi arms will be assembled using slightly different parts. Therefore the range for each armrest will not be the full 220-440mm

Requested Height	Actual Range (Flexi on MPS back)	Requested Height	Actual Range (Flexi on Rehab back)
220-330mm	220-370mm	220-330mm	180-360mm
340-370mm	250-400mm	340-440mm	275-440mm
380-440mm	330-440mm		


ARMPAD ADJUSTMENTS (ALL ARMREST TYPES)






FLEXI ARMREST ADJUSTMENTS



ARMPADS



Allows the arm pads to be in tight to the user while still allowing flip up (Flexi arms only)




If MPS, 50mm, 75mm, gutter or waterfalls armpads are selected, the length of armpad fitted will generally be based on the seat depth. See below table for details. If specific armpad length is required please state this on order form

SEAT DEPTH ▶	12"	13"	14"	15"	16"	17"	18"	19"	20"	21"	22"	23"	24"	
50mm Arm Pad	255	295	305	305	305	355	355	355	380	380	380	380	380	← Armpad length
75mm Arm Pad	255	295	305	305	305	355	355	355	380	380	380	380	380	
Gutter Arm Pad	255	295	305	305	305	355	355	355	380	380	380	380	380	
Waterfall Arm Pad	255	295	305	305	305	355	355	355	380	380	380	380	380	
MPS Arm Pad	-	-	-	-	355	355	405	405	460	460	460	460	460	← Armpad length

magic
Mobility

Joystick Mounting

◀ BACK



Swingback Joystick Release Options	
Magnet	Latch
- Requires more force but no dexterity to release - Can potentially release if wheelchair is being driven erratically or over rough terrain - Releases more easily if impacted	- Requires minimal force but some dexterity to release - Will stay locked at all times unless latch is pushed - Less likely to release if impacted, potentially causing more damage - Can be locked in fully swing-back position (optional) - Custom levers can be screwed to latch to make releasing it easier for clients with less hand movement



JOYSTICK STEM OPTIONS



HEIGHT ADJUSTABLE JOYSTICK STEM



QUICK RELEASE MOUNT

Joystick Knobs

◀ BACK



FOOTPLATE SELECTION



1 PIECE FLIP UP



2 PIECE FLIP UP




1 PIECE FIXED

Fastening bolts can be removed to allow footplate to flip up if required in the future



2 PIECE FIXED

Can adjust each footplate to a different angle



SHORT FOOTPLATES




210



210


© Magis Mobility 2016

Legrest MK2 Footplate Lock




Knee Pads


◀ BACK



**Centre mount legrest
(fixed and power elevating)**
3 axes of adjustment




Power or gas swingaway legrest
Fixed kneepad attached to legrest pivot



Standard swingaway legrest
3 axes of adjustment

Rear Boot Plate

◀ BACK



Cup Holder

◀ BACK



If swingaway legrests are fitted the cup holder mount bracket may be offset higher to avoid top of legrest (depending on armrest height)



Note: Previous cup holder can be fitted if required. See option on order form.



© Magn Mobility 2016

Armrest Bag

◀ BACK



Mounts under armrest
(unless otherwise requested)

If swingback joystick is selected it is generally more suitable to fit accessory bag to opposite side

© Magic Mot

Table Options

◀ BACK

TABLE OPTIONS

Clear polycarbonate table shown below. Black HD plastic is also available (typically without folded edges)

Unless a drawing has been supplied with specifications for cut out sizes for stomach and joystick, the slide in and swingover tables will be provided as a kit

Table stem slides into opposite side to joystick




**SLIDE IN
LARGE TABLE**

Table rests on armpad opposite to stem



**DOUBLE STEM
LARGE TABLE**

For chairs with attendant control ONLY (no joystick)



SWINGAWAY SMALL TABLE
Width 320 x Depth 240



SWINGOVER LARGE TABLE
Table rests on both armpads

© Magic Mobility 2016