



RETAIL

TRU COMFORT 2 SPP CUSHION

DEALERS DETAILS

Dealers Name	
Order Date	
Address	
	Postcode
Email	
Phone	Fax

DELIVERY DETAILS

Dealers Address		
Postcode		
PO No		
Order Type	Quote	Order

DETAILS

This cushion offers skin protection and positioning with a Coolcore, 4-way stretch, breathable cover and an urethane liner.

Multi-layered foam, deep contour, CoolCore cushioncover has pockets on sides.



Weight capacities
158.75 kg: 14"-21" width,
204.11 kg: 22" width+
Overall weight: 1.8kg
Overall height 4.5"
14x14" up to 24x24"
Odd sizes in 17" and 19"
Custom sizes available

CUSHION

TRU-Comfort 2 SPP. 14"-21"W.		ST-TC2SPP			£182
TRU-Comfort 2 SPP. 22"+		ST-TC2SPP			£204
Extra Covers And Liners.	Coolcore Cover.	TC2SPP-NBLK			£99
	Polyurethane Liner.	PUL1709080			£72

		Width										
		14	15	16	17	18	19	20	21	22	24	
Depth	14											
	15											
	16											
	17											
	18											
	19											
	20											
	21											
	22											
	24											

Select your size, then complete the Part Number by adding the Width first then Depth.

ACM CONTOURED RIGIDISER *FOR USE WITH SLING SEATS ONLY.*

15"W x 16"D. <i>CUS142467-1616.</i>	£72	16"W x 18"D. <i>CUS142467.</i>	£72	18"W x 20"D. <i>CUS142466.</i>	£72
15"W x 17"D. <i>CUS142467-1616.</i>	£72	17"W x 16"D. <i>CUS142466-1816.</i>	£72	18"W x 21"D. <i>CUS142466.</i>	£72
15"W x 18"D. <i>CUS142467.</i>	£72	17"W x 21"D. <i>CUS142466.</i>	£72	19"W x 21"D. <i>CUS142466.</i>	£72
15"W x 19"D. <i>CUS142467.</i>	£72	17"W x 22"D. <i>CUS142466.</i>	£72	20"W x 18"D. <i>CUS142465.</i>	£72
15"W x 20"D. <i>CUS142467.</i>	£72	18"W x 16"D. <i>CUS142466.</i>	£72	20"W x 19"D. <i>CUS142465.</i>	£72
16"W x 16"D. <i>CUS142467.</i>	£72	18"W x 18"D. <i>CUS142466.</i>	£72	20"W x 20"D. <i>CUS142465.</i>	£72
16"W x 17"D. <i>CUS142467.</i>	£72	18"W x 19"D. <i>CUS142466.</i>	£72	20"W x 21"D. <i>CUS142465.</i>	£72

NOTES & SPECIAL REQUESTS

Special Order Solution SPP Cushion.

£350


Please detail your special order cushion request.

All special order cushion requests must be approved with the Quantum Custom Department.

Please contact Quantum Inside Sales or Custom Department with questions. Extended lead time will apply.

Cushion part numbers involved will be the cushion part number listed in previous sections modified to meet special order request stated in details below. Drawings are available for sketch details.


REAR VIEW



RIGHT SIDE VIEW



TOP VIEW



PRICE

TOTAL

2



Stealth Products
32 Wedgwood Road, Bicester,
Oxon, OX26 4UL
Phone: 01869 324600
Email: sales@quantumrehab.co.uk
Web: www.quantumrehab.co.uk
TRU COMFORT 2 - 14/02/25 V1

