

Adjusted by:

Tel:

Save the form to the desktop and open it with Adobe Reader.

E-mail:

Date:

The price changes depending on the choices made in the form £:

Standard configuration includes; chassis, standard backrest, castor, footrest alu or titanium, rear wheel spinergy 24" with aluminum pushrims and tire Schwalbe one , and one hand brake.

Max. user weights : Seat width **33-45 cm 100 kg**,

Choice required:

Click on "**blue underlined text**" for more info or image

Optional

Upcharge /Add-price

Marked yellow

Original design

Alternative options

Other options

1) Chassis

[U3 Light](#)

[U3 Light Y-Front](#)

[U3 Light Y-Front 90°](#)

[U3 Light Y-Front Low](#)

[U3 Light Front High 4°](#)

Seat width 33-45 cm, Seat depth 35-46 cm

Seat width 33-45 cm, Seat depth 35, 37.5, 40, 42.5 cm

Seat width 33-45 cm, Seat depth 35-40 cm

Seat width 33-45 cm, Seat depth 35, 37.5, 40, 42.5 cm

Seat width 33-45 cm, Seat depth 35, 37.5 cm

Clear

2) Chassis colour

White- Pearl white

Black- Anodic black

Optional colour - A delivery time of 20 days for optional colour. Fill in art.no. in the square

Clear

3) Seat width

33 cm 36 cm 39 cm 42 cm 45 cm

33 cm L 36 cm L 39 cm L 42 cm L 45 cm L L= light balanced wheelchair, only for U3 Light and U3 Light Y-front

Clear

4) Seat depth

35 cm 37.5 cm 40 cm 42.5 cm

Clear

5) Back rest Height (HH)

[U3 Light](#)

523BBHH L

25 cm

30 cm

35 cm

Clear

Tapered 525BBHH L

Curved back 524BBHH L

[For fixed backrest](#)

Tapered 6 cm

Flared 3 cm

Curved-tapered

6) Balance - center of gravity (COG) [info](#)

Position A

Least easy balanced

Clear

Position B

Medium balance

Position C

Most easy balanced

7) Castor

[90mm X](#)

2100222 (x2/ wheelchair)

Clear

[90 mm S3](#)

2100030 x2

[120 mm S3](#)

2100020 x2

8) Footrest

[Footrest X](#)

36004BB

Clear

[Footrest Titanium](#)

36915BB

[Footrest, extended, reinforced](#)

36006BB [Needed for Freewheel](#)

[Plastic plate for footrest](#)

36010BB

[Carbon footplate](#)

65050BB

9) Brake

[One hand brake, right](#)

Brake S3
4630000

[One hand brake, left](#)

[art.no.](#)

[Clear](#)

10) Rear wheel - without pushrim

All prices for rearwheels are for pair

[24" Spinergy X white](#) spokes

[art.no.](#)

25" [Spinergy X](#) black spokes

[24" Spinergy X black](#) spokes

26" [Spinergy X](#) black spokes

25" [Spinergy X white](#) spokes

[Rear wheel sizes](#)

26" [Spinergy X white](#) spokes

[Clear](#)

Tyre - choose tyre

[Schwalbe One 24", 25"](#)

Rugged 24"

All prices for tyres are for pair

Schwalbe Marathon Plus, 24", 25", 26"

Schwalbe Right Run 22", 24", 25", 26"

PU 24" (massive)

reTyre does not fit PU tyres

reTyre skins 24"

reTyre skins 25"

reTyre skins 26"

[Clear](#)

[Clear](#)

All prices for pushrims are for pair

11) Pushrim

[Aluminum X](#)

[1904404](#) x2 (04=24", 05=25", 06=26")

Paragrip

19051DDBT x 2

Tetragrip

19052DDBT x 2

Maxgrip 24"

1903060 x 2 (contains latex)

[Titanium](#)

190DD11

[Clear](#)

Accessories

Side guard [art.no.](#)

Side guard S3

[video](#)

[Clear](#)

Arm support [art.no.](#)

Arm support S3 *

[video](#)

Makes the wheelchair 15 mm wider

Side guard Carbon [art.no.](#)

Without fender

Small fender

Large fender

Push handles

Foldable S3 [Video](#)

4245000

Max. userweight 100 kg

[Clear](#)

Cushion - Wedge cushion

2,5 cm

572BB40

5 cm

575BB40

[Clear](#)

Wedge cushion

57710BB

Other

Tube protection, front

9900251

Tube protection, rear

9900251

Backrest wedge

4710012

Freewheel

Choose art.no. for Freewheel according to the footrest height of the wheelchair

Read more about Freewheel on the [homepage](#)

Seat bag

Included in order. Choose type.

[Seat mounted](#)

9900212

[Cushion mounted](#)

9900213

[Clear](#)

Set-up adjustments

[Backrest angle](#)

Note measure (A) for desired back rest angle.
10-20 mm. Standard=18 mm

A= _____mm

[Footrest height](#)

Note measure (A) for desired footrest height.

A= _____mm

Extra