

Adjusted by:

Tel:

Save the form to the desktop and  
open it with Adobe Reader.

E-mail:

Date:

The price changes depending on  
the choices made in the form £:

Standard configuration includes; chassis, seat, calf strap, backrest lock, standard colour option, backrest, carbon rear axle, castor and fork, titanium foot rest or foldable footplates, push to lock S3 brake, Panthera crosslaced steel spoke rear wheel, titanium pushrims and Schwalbe right Run tires.

Choice required:

Max. user weights : Seat width **33-42 cm 100 kg**, seat width **45 cm 150 kg**Click on "[blue underlined text](#)" for more info or image**Original design****1) Chassis**[U3](#)Seat width **33-45** cm Seat depth 40 cm

Clear

**2) [Chassis colour](#)****Black**- Anodic Black**Gray** -Metallic Gray

Clear

**3) [Seat width](#)****33 cm****36 cm****39 cm****42 cm****45 cm**

Clear

**4) [Seat depth](#)****40 cm**

Clear

**5) [Back rest Height \(HH\)](#)** [INFO](#)[Standard](#)**20 cm****35 cm**[For fixed backrest](#)

523BBHH VER2

**25 cm****40 cm**

Clear

**30 cm****45 cm**

Clear

**6) [Rear axle](#)**[Carbon std](#)

33023BB

Clear

[Carbon +15 mm](#)

33023BC\*

Use this rear axle when choosing arm support  
at Accessories and/or increased seat height

\* seat width BB+1

[Other](#) rear axle or seat height[Increase the seat height](#)**7) [Castor and fork](#)**[art.no.](#)[90 mm](#)[120 mm](#)[150 mm](#)

Clear

**8) [Footrest](#)**[Footrest Titanium](#)

36001BB

[Footrest U-shape, small heel bar](#)

3600134 (width 33)

Clear

[Titanium, extended, reinforced \\*](#)

36111BB

\* std for S3 Large

[Foldable footplates;](#) [art.no.](#)[Heel strap for footplates](#)[Plastic plate for footrest](#)

36010BB

[Carbon footplate](#)

65050BB

**9) [Brake](#)**[Brake S3](#)

4630000

Clear

**10) Rear wheel - without pushrim**[Rear wheel sizes](#)[24" std](#)[art.no.](#)22" [standard](#)[Clear](#)

---

**Tyre**[Schwalbe Right Run](#) 22", 24", 25", 26"[Clear](#)

---

**11) Pushrim**[Titanium](#)

190DD11 (x2 / wheelchair)

[Clear](#)[Aluminum X](#)

1904404 x2 (04=24", 05=25", 06=26")

DD=Wheel diameter

---

**Accessories****Side guard** [art.no.](#)[Side guard S3](#)[video](#)[Clear](#)**Anti-Tip** [Video](#)[Anti-tip](#) std seat height  
4353000[Clear](#)[Other seat height](#)[Art.no.](#)

## ***Set-up adjustments***

### **Back rest angle**

Note measure (A) for desired back rest angle.  
14-28 mm. Standard = 20 mm

A= \_\_\_\_\_ mm

### **Rear axle adjustment**

Note measure (A) for desired balance  
Standard = 60 mm. at raised rear axle = 51 mm

A= \_\_\_\_\_ mm

### **Footrest height - Foot plates**

Note measure (A) and (B) for desired footrest

A= \_\_\_\_\_ mm    B= \_\_\_\_\_ mm

## ***Extra***