

S800 BY SCOOBY®

The Modern Mobility Scooter
for an Active Lifestyle



Drive the New Way

Modern, stylish, and built for adventure—the Sterling S800 by Scoozy takes mobility to the next level.

Large wheels, advanced suspension and a 70ah lithium battery provide an unparalleled smooth and stable ride on any terrain, with a range of up to 50 km.

The intuitive bike-style controls make the S800 by Scoozy easy to navigate; it's the ideal choice for those who want to stay active and independent—whether exploring city streets or venturing into nature.

Turn heads wherever you go. With sleek and modern aesthetics, the S800 by Scoozy isn't just a mode of transport – it's a statement.



Scan the QR code with
your phone's camera
to watch video.





This Isn't Just a Mode
of Transport...

It's a Statement.

Designed for Your Journey

The S800 by Scoozy is designed to enhance every ride with **comfort, power, and reliability**. Whether you're navigating different terrains, embarking on long adventures, or tackling steep climbs, its advanced features keep you moving with ease and confidence.

Grip That Keeps You in Control

Navigate nature trails, urban environments, and uneven roads with confidence. The robust 18" wheels are designed to provide **excellent grip and stability**, ensuring a smooth and secure ride on various surfaces.



Go Further on Fewer Charges

Enjoy the freedom to explore further with an unrivalled range of up to **50 km*** on a single charge. Perfect for long rides and extended adventures.



* Range depends on the ground surface, the weight of the rider, driving behaviour, temperature and terrain.



Superior Suspension for a Smooth Ride

From the moment you sit down, experience the most comfortable ride whether tackling city streets or rugged terrain. With both front and rear suspension, **every movement is cushioned** for a softer experience. The rear suspension is also adjustable, allowing you to fine-tune the setup for your perfect level of support and comfort.



Effortless Climbing and Acceleration

Smooth yet **responsive acceleration** delivered by a high-performance motor, that provides impressive torque and power, conquering high kerbs and steep climbs with ease.



Speed – Get There Faster

With a **top speed of 18 km/h**, you'll breeze through your journey while maintaining complete control and safety.

Designed for Your Comfort

The Sterling S800 by Scoozy ensures every journey is comfortable with its **ergonomic design**, its soft, cushioned seat, supportive backrest, and easy adjustments for your perfect position.



Scan the QR code with your phone's camera to see how easy it is to adjust the seat.



Spacious Design

The S800 by Scoozy has a large open footwell for you to comfortably place your feet, reducing strain on the body and offering a smoother, more enjoyable ride over bumps and obstacles.



Adjustable Seating for Personalized Comfort

Find your perfect position with easy seating adjustments. The Sterling S800 by Scoozy adapts to suit your ideal posture, enhancing your overall driving experience.




Easy to Use Controls


Navigate effortlessly with the familiar **bicycle-style handlebar, equipped with intuitive controls**: Thumb throttle, cruise control, and hand-operated brakes give you full control at all times.

Thanks to the ergonomic design, your hands, arms, and back stay strain-free, letting you enjoy the journey as much as the destination.

The large colour display, featuring NFC-tag locking technology, to keep your S800 safe, provides a modern user experience.



Scan the QR code with your phone's camera to see how easy the S800 is to drive.





A Scooter Like No Other

Step into tomorrow with the **S800 by Scoozy**, a masterpiece of modern design and cutting-edge engineering. Crafted for those who dare to redefine their journey, this scooter boasts a **sleek, futuristic silhouette, clean lines, and premium finishes**. Its streamlined body and subtle curves turn heads at every corner, offering a ride that's as stylish as it is advanced.



Everything You Need at a Glance

Full-colour **LCD display**, **NFC locking** function and **USB Charge Port**.



Navigate Effortlessly

Navigate effortlessly with the familiar **bike-style handlebar**, featuring a thumb throttle, hand-operated brakes, and cruise control.



Personalize Your Way

Easily mount a range of **high quality accessories** on the accessory bar.



Be Seen, Be Safe

Integrated **LED light and indicators** provide reliable lighting day and night.



Long Range – Go the Distance

The 70Ah lithium battery provides a **50 km range** for worry-free travel - even in sub-zero temperatures - and allows charging anytime without performance loss.



Ergonomic Seating

Premium seating experience with **soft, cushioned seat and supportive backrest**.



Superior Traction

Robust **18" wheels** provide unbeatable grip and stability.



Smooth Ride

Advanced suspension absorbs shocks for a comfortable ride on any terrain.



Comfortable Leg Position

Spacious legroom ensures an optimal ergonomic seating position.



Stay Connected with the Sterling App*

The Sterling App offers valuable insights and support, keeping both users and their caregivers or family members informed and empowered every step of the way.



Stay Safe with Instant Alerts & Emergency Notifications:


Get immediate notifications for impacts, incidents, or falls, with automatic alerts sent to emergency contacts for quick assistance.



Maximize Battery Life & Performance: Track charging sessions, receive low battery alerts, and monitor battery health to keep your scooter running smoothly.



Effortlessly Track Your Ride & Stay Informed: Keep an eye on location, battery life, and scooter health, while reviewing trip history to stay updated on usage and maintenance.




Accessories

Personalize your S800 by Scoozy with a range of high-quality accessories. All can be effortlessly mounted on the mounting tubes, giving you plenty of options to enhance your ride.

Mounting tubes

There is a mounting tube on the back of the frame and under the seat to mount various universal mobility scooter accessories, such as a canopy, reflective warning triangle (mandatory in Germany), holder for sticks, walker, or oxygen bottle.

Scan the QR code with your phone's camera to see how easy it is to install accessories.



Enjoy quick and easy ordering with options available on the order form.



KLIckfix Adapter for Third Party Accessories

Order the KLIckfix adapter separately to mount any suitable KLIckfix accessories yourself! The KLIckfix handlebar adapter with clamp is required for basket mounting on the accessory bar.



Rear Review Mirror (LH + RH)

The mirror for LH + RH assembly from black plastic, size 105mm.



Detachable Basket (with KLIckfix adapter)

The KLIckfix Lamello 18L basket provides easy access and a secure place for your belongings.



RAM Mount Phone

The RAM Mounts U-Bolt phone holder can be placed on the front accessory bar. You can easily click your phone in and out of the X-grip clip. Place your device horizontally or vertically in the holder.

Take advantage of compatibility with a variety of third-party accessories to easily enhance your S800 with online options.



Vario Rack (Sport)*

The KLIckfix compact, lightweight racks click onto the Handlebar Adapter at two heights. Clip-on hooks adjust straps on the aluminum frame, and an adjustable central strap with a buckle and reflector secures luggage.



Reisenthel Fabric Basket*

Add color with the KLIckfix Reisenthel fabric bike basket, with a front pouch, zippered inner compartment, phone pocket, drawstring opening, and detachable pouch. Compatible with the KLIckfix Handlebar Adapter System.



Doggy Basket*

The KLIckfix Doggy is a comfortable handlebar bag for pets up to 5 kg. Includes carrying handle, reflector, zip pocket, and mesh compartments. Attaches securely to the Handlebar Adapter.



Windshield*

Consider a Polisport windscreen from a child bike seat. This windscreen has the same dimensions as those for a (mobility) scooter and fits perfectly on the S800 by Scoozy front accessory bar**.

*Ordered through 3rd party suppliers only.

**The KLIckfix adapter with basket will not fit together with the windscreen.

Technical Data

Overall Width / Length 680 x 1370 mm
(26.8 x 53.9")

Effective Seat Depth 430 mm (16.9")

Effective Seat Width 508 mm (20")

Max. Speed 18 kph (11.2 mph)

Max. User Weight 160 kg (352.7 lb)

Range of Action
(ISO 7176-15)* 50 km

Battery Size Lithium-ion 24V 70 Ah

Tyres 18" air tyres

Turning Radius 1400 mm (55.1")

Weight Incl. Batteries 100 kg (220.5 lb)

Control Thumb throttle left
or right

Kerb Climbing Max. 100 mm (3.9")

Colour Pale Green
(RAL 6021)

*Range depends on the ground surface, the weight of the rider, driving behaviour, temperature and terrain.

Sunrise Medical GmbH

Kahlbachring 2-4

69254 Malsch/HD

Tel.: +49 (0) 72 53 / 9 80-0

Fax: +49 (0) 72 53 / 9 80-222

E-Mail: international@sunrisemedical.de



S800 By Scoozy Brochure Rev. 2 2025

Now that you have read all about the Sterling S800 by Scoozy, why not arrange a free demonstration and experience it for yourself?



Scan the QR code with your phone's camera to book a free test-drive today.

For further information on the full specification, options, and accessories, please refer to the order form. All information is subject to change without notification. Please consult Sunrise Medical with any queries you may have.



If you are visually impaired, this document can be viewed in PDF format at www.sunrisemedical.eu.



www.sunrisemedical.eu



WHERE YOU WANT TO BE.

