



Wheelchair Details:


Manufacturer:

Model:

Inside Tube Measurement: _____ cm

Lateral System

Central System



Klaxon Klick Race

Klick Race Hub only - Sections A B and C

Klick Race Accessories - Sections C - G

Klick Race with Central Mounting System - Sections A, B and H

Klick Race with Lateral Mounting Systems - Sections A, B and I

A: Race Hub	Article Number	Price
Klaxon Klick Race - 20 km/h - Thumb throttle (R010014S)	<input type="checkbox"/> KXHK0005	£5,250.00
Klaxon Klick Race - 20 km/h - Tetra handlebar* (S010014T)	<input type="checkbox"/> KXHK0006	£5,500.00

1500 Watt High Performance 50Nm of torque

14" wide wheel

LCD Display

5 Speeds, plus reverse function

Maximum user weight 120 kg

B: Advanced Electronic Features	Article Number	Price
Intellegent Cruise Control	<input checked="" type="checkbox"/>	Included
Electronic Braking System	<input checked="" type="checkbox"/>	Included
Motor Direct Control	<input checked="" type="checkbox"/>	Included
HHC - Hill Holder Control (CC-2-PROG-HHC)	<input checked="" type="checkbox"/> KXAK0042	Included

C: Batteries (Charger supplied)	Article Number	Price
Race Battery for Klick 48V 11,6 Ah, range 30-50 km (CM-1-11AH)	<input checked="" type="checkbox"/> KXBK0012	Included
Race Lite for Klick 48V 5A, range 10-25 km (CM-1-5AH)	<input type="checkbox"/> KXBK0014	£0.00
Fly Race for Klick 48V 2,9 Ah, range 5-12 km (CM-1-2,9AH)	<input type="checkbox"/> KXBK0013	£0.00



D: Spare Batteries/Charger	Article Number	Price
Race Battery for Klick 48V 11,6 Ah, range 30-50 km (CM-1-11AH)	<input type="checkbox"/> KXBK0012	£700.00
Race Lite for Klick 48V 5Ah, range 10-25 km (CM-1-5AH)	<input type="checkbox"/> KXBK0014	£600.00
Fly Race for Klick 48V 2,9 Ah, range 5-12 km (CM-1-2,9AH)	<input type="checkbox"/> KXBK0013	£500.00
Klaxon Klick Power Charger 48V 11Ah (CC-88-0002)	<input type="checkbox"/> KXSK0122	£70.00



E: Adaptions	Article Number	Price
Universal Aluminium Handlebar Extension (CC-25-0006)	<input type="checkbox"/> KXSK0128	£160.00
Additional Traction Weight Kit - 2kg (pair) (CC-93-0000 x2)	<input type="checkbox"/> KXAK0009 x2	£200.00
Bullhorn handlebar grips (pair) (T-HORN-0000)	<input type="checkbox"/> KXAK0030	£20.00
Universal clamping ring for tetra handlebar* handles, right (T-HANDLE-00)	<input type="checkbox"/> KXAK0026	£80.00
Universal clamping ring for tetra handlebar* handles, left (T-HANDLE-00)	<input type="checkbox"/> KXAK0026	£80.00
Vertical Tetra handle for tetra handlebar*, right (T-HANDLE-02)	<input type="checkbox"/> KXAK0028	£80.00
Vertical Tetra handle for tetra handlebar*, left (T-HANDLE-02)	<input type="checkbox"/> KXAK0028	£80.00
Tripod handle for Tetra handlebars*, right (T-HANDLE-03)	<input type="checkbox"/> KXAK0029	£80.00
Tripod handle for Tetra handlebars*, left (T-HANDLE-03)	<input type="checkbox"/> KXAK0029	£80.00
Horizontal Tetra Handle*, right (T-HANDLE-01)	<input type="checkbox"/> KXAK0027	£60.00
Horizontal Tetra Handle*, left (T-HANDLE-01)	<input type="checkbox"/> KXAK0027	£60.00

F: Safety Kit	Article Number	Price
White front LED Light (CC-83-KIT-2024)	<input checked="" type="checkbox"/> KXAK0007	Included
Single rear-view Mirror (SPE-0000)	<input type="checkbox"/> KXAK0025	£40.00
Basic white front LED Light (LED-WHITE-0000)	<input type="checkbox"/> KXAK0021	£10.00
Basic red rear LED Light (LED-RED-0000)	<input type="checkbox"/> KXAK0020	£10.00

G: Accessories	Article Number	Price
Fat wheels for 12,7 mm (1/2 ") axle (FATWHEELS-01-12,7)	<input type="checkbox"/> KXAK0012	£320.00
Fat wheels for 12 mm axle (FATWHEELS-01-12)	<input type="checkbox"/> KXAK0011	£320.00
Shopping bag frame (T-SHOPBAG)	<input type="checkbox"/> KXAK0031	£30.00
Shopping bag front (SHOPBAG-0001)	<input type="checkbox"/> KXAK0024	£80.00
Shopping bag side (1 piece) right (SHOPBAG-0000)	<input type="checkbox"/> KXAK0023	£80.00
Shopping bag side (1 piece) left (SHOPBAG-0000)	<input type="checkbox"/> KXAK0023	£80.00
Basket linking kit (KIT_TEL_0001)	<input type="checkbox"/> KXAK0019	£140.00
Basket 10L (BASKET-0000)	<input type="checkbox"/> KXAK0003	£60.00
Basket 16L (BASKET-0001)	<input type="checkbox"/> KXAK0004	£60.00
Flybag 14" (FLYBAG00)	<input type="checkbox"/> KXAK0015	£280.00
Handlebar phone holder (PHONE HOLDER)	<input type="checkbox"/> KXAK0022	£30.00
Battery sticker - Reflex (CM-1-STICKER-R)	<input type="checkbox"/> KXAK0036	£20.00
Battery sticker - Light blue (CM-1-STICKER-RA)	<input type="checkbox"/> KXAK0037	£20.00
Battery Sticker - Blue (CM-1-STICKER-RB)	<input type="checkbox"/> KXAK0038	£20.00
Battery Sticker - Fuxia (CM-1-STICKER-RF)	<input type="checkbox"/> KXAK0041	£20.00
Battery Sticker - Red (CM-1-STICKER-RR)	<input type="checkbox"/> KXAK0039	£20.00
Battery Sticker - Green (CM-1-STICKER-RV)	<input type="checkbox"/> KXAK0040	£20.00

Tyre Options	Article Number	Price
Allroad Tyre 14" x 2.125 (CC-3-1401)	<input type="checkbox"/> KXSK0046	£0.00
Cross Tyre 14" x 2.50 (CC-3-2510)	<input type="checkbox"/> KXSK0120	£0.00
Enduro Tyre 14" x 2.125 (CC-3-1402)	<input type="checkbox"/> KXSK0047	£0.00
Allroad 14" x 2.50 (CC-3-1425)	<input type="checkbox"/> KXSK0119	£0.00

H1: Central Hooking System - Connector	Article Number	Price
	<input type="checkbox"/> KXLK0009	£0.00
	<input type="checkbox"/> KXLK0007	£0.00
Complete Connector Small (LSTD-0002)	<input type="checkbox"/> KXLK0008	£0.00
Complete Connector Standard (LSTD-0000)	<input type="checkbox"/> KXLK0011	£0.00
Complete Connector Extra Long (LSTD-0001)	<input type="checkbox"/> KXLK0012	£0.00
Complete Connector Tetra Standard (LTET-000)	<input type="checkbox"/> KXLK0013	£0.00
Complete Connector Tetra Extra Long (LTET-0001)		
Complete Connector Tetra Tetra Support Velco Grip (TEV02)		

H2: Crossbeam for Central Hooking System	Article Number	Price
		
Measure between the front frame under the seat canvas		
	<input type="checkbox"/> KXLKT0084	£0.00
Complete Crossbeam 271-290 mm (TRAV2729)	<input type="checkbox"/> KXLKT0085	£0.00
Complete Crossbeam 291-360 mm (TRAV 2936)	<input type="checkbox"/> KXLKT0086	£0.00
Complete Crossbeam 361-480 mm (TRAV3648)	<input type="checkbox"/> KXLKT0087	£0.00
Complete Crossbeam 481-540 mm (TRAV4854)	<input type="checkbox"/> KXLKT0088	£0.00

H3: Other Components for Central Hooking System	Article Number	Price
Quick Release Locking Pins with Ring Handle (pair) (DF-34-0001 x2)	<input type="checkbox"/> KXLKT0060 (x2)	£0.00

H4: Central Hooking System - Tube Clamps for Parallel Frames	Article Number	Price
Round Tube Clamps \varnothing 17mm - 18,9mm (AGP25 x2 + RID25-18-0000)	<input type="checkbox"/> KXLKT0014 x2 + KXLKT0068	£0.00
Round Tube Clamps \varnothing 19mm - 20,9mm (AGP25 x2 + RID25-20-0000)	<input type="checkbox"/> KXLKT0014 x2 + KXLKT0069	£0.00
Round Tube Clamps \varnothing 21mm - 22,9 mm (AGP25 x2 + RID25-22-0000)	<input type="checkbox"/> KXLKT0014 x2 + KXLKT0070	£0.00
Round Tube Clamps \varnothing 23mm - 22,9mm (AGP25 x2 + RID25-23-0000)	<input type="checkbox"/> KXLKT0014 x2 + KXLKT0071	£0.00
Round Tube Clamps \varnothing 24.4mm - 26,4mm (AGP25 x2)	<input type="checkbox"/> KXLKT0014 x2	£0.00
Round Tube Clamps \varnothing 27mm - 28,9mm (AGP32 x2 + RID32-28-0000)	<input type="checkbox"/> KXLKT0016 x2 + KXLKT0072	£0.00
Round Tube Clamps \varnothing 29mm - 30,9mm (AGP32 x2 + RID32-30-0000)	<input type="checkbox"/> KXLKT0016 x2 + KXLKT0073	£0.00
Round Tube Clamps \varnothing 31mm - 33,4mm (AG32 x2)	<input type="checkbox"/> KXLKT0016 x2	£0.00
Round Tube Clamps \varnothing 34mm - 36mm (AGP35 x2)	<input type="checkbox"/> KXLKT0017 x2	£0.00
Oval width 25mm x height 35mm (Progeo) (AGPOV x2)	<input type="checkbox"/> KXLKT0021 x2	£0.00
Oval width 27mm x height 33mm (Sunrise) (AGPOVSX x2)	<input type="checkbox"/> KXLKT0024 x2	£0.00
Oval Vermiren Trigo (AGPOVT x2)	<input type="checkbox"/> KXLKT0025 x2	£0.00
Oval Progeo Carbon Noir 2.0/Duke (AGPOVPC x2)	<input type="checkbox"/> KXLKT0022 x2	£0.00


H5: Central Hooking System - Tube Clamps for Convergent Frames	Article Number	Price
Round Tube Clamps \varnothing 17mm - 18,9mm (AGCP125 x2 + RID25-18-0000)	<input type="checkbox"/> KXLKT0001 x2 + KXLKT0068	£0.00
Round Tube Clamps \varnothing 19mm - 20,9mm (AGCP125 x2 + RID25-20-0000)	<input type="checkbox"/> KXLKT0001 x2 + KXLKT0069	£0.00
Round Tube Clamps \varnothing 21mm - 22,9 mm (AGCP125 x2 + RID25-22-0000)	<input type="checkbox"/> KXLKT0001 x2 + KXLKT0070	£0.00
Round Tube Clamps \varnothing 23mm - 22,9mm (AGCP125 x2 + RID25-23-0000)	<input type="checkbox"/> KXLKT0001 x2 + KXLKT0071	£0.00
Round Tube Clamps \varnothing 24.4mm - 26,4mm (AGCP125 x2)	<input type="checkbox"/> KXLKT0001 x2	£0.00
Round Tube Clamps \varnothing 27mm - 28,9mm (AGCP132 x2 + RID32-28-0000)	<input type="checkbox"/> KXLKT0002 x2 + KXLKT0072	£0.00
Round Tube Clamps \varnothing 29mm - 30,9mm (AGCP132 x2 + RID32-30-0000)	<input type="checkbox"/> KXLKT0002 x2 + KXLKT0073	£0.00
Round Tube Clamps \varnothing 31mm - 33,4mm (AG132 x2)	<input type="checkbox"/> KXLKT0002 x2	£0.00
Round Tube Clamps \varnothing 34mm - 36mm (AGP135 x2)	<input type="checkbox"/> KXLKT0003 x2	£0.00

I1: Lateral Linking System (alternative mounting option)	Article Number	Price
	<input type="checkbox"/> KXLK0006	£0.00

I2: Lateral (Side) Linking System Tube Clamps	Article Number	Price
Round Tube Clamps \varnothing 17mm - 18,9mm (AGLSX25 + AGLDX25 + RID25-18-0000)	<input type="checkbox"/> KXLKT0009 + KXLKT0004 + KXLKT0068	£0.00
Round Tube Clamps \varnothing 19mm - 20,9mm (AGLSX25 + AGLDX25 + RID25-20-0000)	<input type="checkbox"/> KXLKT0009 + KXLKT0004 + KXLKT0069	£0.00
Round Tube Clamps \varnothing 21mm - 22,9 mm (AGLSX25 + AGLDX25 + RID25-22-0000)	<input type="checkbox"/> KXLKT0009 + KXLKT0004 + KXLKT0070	£0.00
Round Tube Clamps \varnothing 23mm - 22,9mm (AGLSX25 + AGLDX25 + RID25-23-0000)	<input type="checkbox"/> KXLKT0009 + KXLKT0004 + KXLKT0071	£0.00
Round Tube Clamps \varnothing 24.4mm - 26,4mm (AGLSX25 + AGLDX25)	<input type="checkbox"/> KXLKT0009 + KXLKT0004	£0.00
Round Tube Clamps \varnothing 27mm - 28,9mm (AGLSX32 + AGLDX32 + RID32-28-0000)	<input type="checkbox"/> KXLKT0011 + KXLKT0006 + KXLKT0072	£0.00
Round Tube Clamps \varnothing 29mm - 30,9mm (AGLSX32 + AGLDX32 + RID32-30-0000)	<input type="checkbox"/> KXLKT0011 + KXLKT0006 + KXLKT0073	£0.00
Round Tube Clamps, \varnothing 31mm - 33,4mm (AGLSX32 + AGLDX32)	<input type="checkbox"/> KXLKT0011 + KXLKT0006	£0.00
Oval width 25mm x height 35mm (Progeo) (AGLSXOVS + AGLDXOVS)	<input type="checkbox"/> KXLKT0014 + KXLKT0008	£0.00
Oval width 27mm x height 33mm (Sunrise) (AGLSXOV + AGLDXOV)	<input type="checkbox"/> KXLKT0007 + KXLKT0012	£0.00

Helpful Measures for Set Up (Optional)

A. Distance, measured vertically from the highest point of the knees to the floor	cm
B. Arm length extended measured from the centre of the shoulder to the palm	cm
C. Distance, measured horizontally from front of footplate to front of foot	cm
D. Distance, measured vertically from top of shoulder to the floor	cm
E. Trunk depth, measured horizontally at navel height	cm



Notes