




Klaxon Klick Scrambler

Klick Scrambler Hub only - Sections A B and C
Klick Scrambler Accessories - Sections C - G
Klick Scrambler with Central Mounting System - Sections A, B and H
Klick Scrambler with Lateral Mounting Systems - Sections A, B and I



Wheelchair Details:	
Manufacturer:	
Model:	
Inside Tube Measurement:	cm
Lateral System	☐
Central System	☐

A: Scrambler Hub	Article Number	Price
Klaxon Klick Scrambler - 20 km/h - Thumb throttle (S002512S)	☐ KXHK0012	£4,950.00
Klaxon Klick Scrambler - 20 km/h - Tetra handlebar* (S002512T)	☐ KXHK0013	£5,200.00

- 750 Watt 50Nm of torque
- 14" wide wheel
- LCD Display
- 5 Speeds, plus reverse function
- Maximum user weight 120 kg

B: Advanced Electronic Features	Article Number	Price
Intellegent Cruise Control	☒	Included
Electronic Braking System	☒	Included
Motor Direct Control	☒	Included
HHC - Hill Holder Control (CC-2-PROG-HHC)	☒ KXAK0042	Included

C: Batteries (Charger supplied)	Article Number	Price
Race Battery for Klick 48V 11,6 Ah, range 30-50 km (CM-1-11AH)	☒ KXBK0012	Included
Race Lite for Klick 48V 5Ah, range 10-25 km (CM-1-5AH)	☐ KXBK0014	£0.00
Fly Race for Klick 48V 2,9 Ah, range 5-12 km (CM-1-2,9AH)	☐ KXBK0013	£0.00

D: Spare Batteries/Charger	Article Number	Price
Race Battery for Klick 48V 11,6 Ah, range 30-50 km (CM-1-11AH)	☐ KXBK0012	£700.00
Race Lite for Klick 48V 5Ah, range 10-25 km (CM-1-5AH)	☐ KXBK0014	£600.00
Fly Race for Klick 48V 2,9 Ah, range 5-12 km (CM-1-2,9AH)	☐ KXBK0013	£500.00
Klaxon Klick Power Charger 48V 11Ah (CC-88-0002)	☐ KXSK0122	£70.00

E: Adaptions	Article Number	Price
Universal Aluminium Handlebar Extension (CC-25-0006)	☐ KXSK0128	£160.00
Additional Traction Weight Kit - 2kg (pair) (CC-93-0000 x2)	☐ KXAK0009 x2	£200.00
Bullhorn handlebar grips (pair) (T-HORN-0000)	☐ KXAK0030	£20.00
Universal clamping ring for tetra handlebar* handles, right (T-HANDLE-00)	☐ KXAK0026	£80.00
Universal clamping ring for tetra handlebar* handles, left (T-HANDLE-00)	☐ KXAK0026	£80.00
Vertical Tetra handle for tetra handlebar*, right (T-HANDLE-02)	☐ KXAK0028	£80.00
Vertical Tetra handle for tetra handlebar*, left (T-HANDLE-02)	☐ KXAK0028	£80.00
Tripod handle for Tetra handlebars*, right (T-HANDLE-03)	☐ KXAK0029	£80.00
Tripod handle for Tetra handlebars*, left (T-HANDLE-03)	☐ KXAK0029	£80.00
Horizontal Tetra Handle*, right (T-HANDLE-01)	☐ KXAK0027	£60.00
Horizontal Tetra Handle*, left (T-HANDLE-01)	☐ KXAK0027	£60.00

F: Safety Kit

	Article Number	Price
White front LED Light (CC-83-KIT-2024)	☒ KXAK0007	Included
Single rear-view Mirror (SPE-0000)	☐ KXAK0025	£40.00
Basic white front LED Light (LED-WHITE-0000)	☐ KXAK0021	£10.00
Basic red rear LED Light (LED-RED-0000)	☐ KXAK0020	£10.00

G: Accessories

	Article Number	Price
Fat wheels for 12,7 mm (1/2 ") axle (FATWHEELS-01-12,7)	☐ KXAK0012	£320.00
Fat wheels for 12 mm axle (FATWHEELS-01-12)	☐ KXAK0011	£320.00
Shopping bag frame (T-SHOPBAG)	☐ KXAK0031	£30.00
Shopping bag front (SHOPBAG-0001)	☐ KXAK0024	£80.00
Shopping bag side (1 piece) right (SHOPBAG-0000)	☐ KXAK0023	£80.00
Shopping bag side (1 piece) left (SHOPBAG-0000)	☐ KXAK0023	£80.00
Basket linking kit (KIT_TEL_0001)	☐ KXAK0019	£140.00
Basket 10L (BASKET-0000)	☐ KXAK0003	£60.00
Basket 16L (BASKET-0001)	☐ KXAK0004	£60.00
Flybag 14" (FLYBAG00)	☐ KXAK0015	£280.00
Handlebar phone holder (PHONE HOLDER)	☐ KXAK0022	£30.00
Battery sticker - Reflex (CM-1-STICKER-R)	☐ KXAK0036	£20.00
Battery sticker - Light blue (CM-1-STICKER-RA)	☐ KXAK0037	£20.00
Battery Sticker - Blue (CM-1-STICKER-RB)	☐ KXAK0038	£20.00
Battery Sticker - Fuxia (CM-1-STICKER-RF)	☐ KXAK0041	£20.00
Battery Sticker - Red (CM-1-STICKER-RR)	☐ KXAK0039	£20.00
Battery Sticker - Green (CM-1-STICKER-RV)	☐ KXAK0040	£20.00

Tyre Options


	Article Number	Price
Allroad HT Tyre 13,3" (300-8) (CC-3-3008)	☐ KXSK0129	£0.00
Motard HT Tyre 13,3" (300-8) (CC-3-3108)	☐ KXSK0127	£0.00

H1: Central Hooking System - Connector



	Article Number	Price
Complete Connector Small (LSTD-0002)	☐ KXLK0009	£0.00
Complete Connector Standard (LSTD-0000)	☐ KXLK0007	£0.00
Complete Connector Extra Long (LSTD-0001)	☐ KXLK0008	£0.00
Complete Connector Tetra Standard (LTET-000)	☐ KXLK0011	£0.00
Complete Connector Tetra Extra Long (LTET-0001)	☐ KXLK0012	£0.00
Complete Connector TetraTetra Support Velco Grip (TEV02)	☐ KXLK0013	£0.00

H2: Crossbeam for Central Hooking System



Measure between the front frame under the canvas



	Article Number	Price
Complete Crossbeam 230-270 mm (TRAV2327)	☐ KXLKT0084	£0.00
Complete Crossbeam 271-290 mm (TRAV2729)	☐ KXLKT0085	£0.00
Complete Crossbeam 291-360 mm (TRAV 2936)	☐ KXLKT0086	£0.00
Complete Crossbeam 361-480 mm (TRAV3648)	☐ KXLKT0087	£0.00
Complete Crossbeam 481-540 mm (TRAV4854)	☐ KXLKT0088	£0.00

H3: Other Components for Central Hooking System

	Article Number	Price
Quick Release Locking Pins with Ring Handle (pair) (DF-34-0001 x2)	☐ KXLKT0060 (x2)	£0.00

H4: Central Hooking System - Tube Clamps for Parallel Frames	Article Number	Price
Round Tube Clamps Ø 17mm - 18,9mm (AGP25 x2 + RID25-18-0000)	☐ KXLKT0014 x2 + KXLKT0068	£0.00
Round Tube Clamps Ø 19mm - 20,9mm (AGP25 x2 + RID25-20-0000)	☐ KXLKT0014 x2 + KXLKT0069	£0.00
Round Tube Clamps Ø 21mm - 22,9 mm (AGP25 x2 + RID25-22-0000)	☐ KXLKT0014 x2 + KXLKT0070	£0.00
Round Tube Clamps Ø 23mm - 22,9mm (AGP25 x2 + RID25-23-0000)	☐ KXLKT0014 x2 + KXLKT0071	£0.00
Round Tube Clamps Ø 24.4mm - 26,4mm (AGP25 x2)	☐ KXLKT0014 x2	£0.00
Round Tube Clamps Ø 27mm - 28,9mm (AGP32 x2 + RID32-28-0000)	☐ KXLKT0016 x2 + KXLKT0072	£0.00
Round Tube Clamps Ø 29mm - 30,9mm (AGP32 x2 + RID32-30-0000)	☐ KXLKT0016 x2 + KXLKT0073	£0.00
Round Tube Clamps Ø 31mm - 33,4mm (AG32 x2)	☐ KXLKT0016 x2	£0.00
Round Tube Clamps Ø 34mm - 36mm (AGP35 x2)	☐ KXLKT0017 x2	£0.00
Oval width 25mm x height 35mm (Progeo) (AGPOV x2)	☐ KXLKT0021 x2	£0.00
Oval width 27mm x height 33mm (Sunrise) (AGPOVSX x2)	☐ KXLKT0024 x2	£0.00
Oval Vermiren Trigo (AGPOVT x2)	☐ KXLKT0025 x2	£0.00
Oval Progeo Carbon Noir 2.0/Duke (AGPOVPC x2)	☐ KXLKT0022 x2	£0.00

H5: Central Hooking System - Tube Clamps for Convergent Frames	Article Number	Price
Round Tube Clamps Ø 17mm - 18,9mm (AGCP125 x2 + RID25-18-0000)	☐ KXLKT0001 x2 + KXLKT0068	£0.00
Round Tube Clamps Ø 19mm - 20,9mm (AGCP125 x2 + RID25-20-0000)	☐ KXLKT0001 x2 + KXLKT0069	£0.00
Round Tube Clamps Ø 21mm - 22,9 mm (AGCP125 x2 + RID25-22-0000)	☐ KXLKT0001 x2 + KXLKT0070	£0.00
Round Tube Clamps Ø 23mm - 22,9mm (AGCP125 x2 + RID25-23-0000)	☐ KXLKT0001 x2 + KXLKT0071	£0.00
Round Tube Clamps Ø 24.4mm - 26,4mm (AGCP125 x2)	☐ KXLKT0001 x2	£0.00
Round Tube Clamps Ø 27mm - 28,9mm (AGCP132 x2 + RID32-28-0000)	☐ KXLKT0002 x2 + KXLKT0072	£0.00
Round Tube Clamps Ø 29mm - 30,9mm (AGCP132 x2 + RID32-30-0000)	☐ KXLKT0002 x2 + KXLKT0073	£0.00
Round Tube Clamps Ø 31mm - 33,4mm (AG132 x2)	☐ KXLKT0002 x2	£0.00
Round Tube Clamps Ø 34mm - 36mm (AGP135 x2)	☐ KXLKT0003 x2	£0.00

I1: Lateral Linking System (alternative mounting option)	Article Number	Price
--	----------------	-------




Lateral Linking System (CL-KIT-0000)	☐ KXLK0006	£0.00
--------------------------------------	-----------------------------------	-------

I2: Lateral (Side) Linking System Tube Clamps	Article Number	Price
Round Tube Clamps Ø 17mm - 18,9mm (AGLSX25 + AGLDX25 + RID25-18-0000)	☐ KXLKT0009 + KXLKT0004 + KXLKT0068	£0.00
Round Tube Clamps Ø 19mm - 20,9mm (AGLSX25 + AGLDX25 + RID25-20-0000)	☐ KXLKT0009 + KXLKT0004 + KXLKT0069	£0.00
Round Tube Clamps Ø 21mm - 22,9 mm (AGLSX25 + AGLDX25 + RID25-22-0000)	☐ KXLKT0009 + KXLKT0004 + KXLKT0070	£0.00
Round Tube Clamps Ø 23mm - 22,9mm (AGLSX25 + AGLDX25 + RID25-23-0000)	☐ KXLKT0009 + KXLKT0004 + KXLKT0071	£0.00
Round Tube Clamps Ø 24.4mm - 26,4mm (AGLSX25 + AGLDX25)	☐ KXLKT0009 + KXLKT0004	£0.00
Round Tube Clamps Ø 27mm - 28,9mm (AGLSX32 + AGLDX32 + RID32-28-0000)	☐ KXLKT0011 + KXLKT0006 + KXLKT0072	£0.00
Round Tube Clamps Ø 29mm - 30,9mm (AGLSX32 + AGLDX32 + RID32-30-0000)	☐ KXLKT0011 + KXLKT0006 + KXLKT0073	£0.00
Round Tube Clamps, Ø 31mm - 33,4mm (AGLSX32 + AGLDX32)	☐ KXLKT0011 + KXLKT0006	£0.00
Oval width 25mm x height 35mm (Progeo) (AGLSXOVS + AGLDXOVS)	☐ KXLKT0014 + KXLKT0008	£0.00
Oval width 27mm x height 33mm (Sunrise) (AGLSXOV + AGLDXOV)	☐ KXLKT0007 + KXLKT0012	£0.00

Helpful Measures for Set Up (Optional)

A. Distance, measured vertically from the highest point of the knees to the floor	cm
B. Arm length extended measured from the centre of the shoulder to the palm	cm
C. Distance, measured horizontally from front of footplate to front of foot	cm
D. Distance, measured vertically from top of shoulder to the floor	cm
E. Trunk depth, measured horizontally at navel height	cm



Notes