

Deliver to:  
 Name:  
 Street:  
 Town: Postal code:  
 E-Mail:  
 Tel: Fax:



Invoice to:  
 Name:  
 Street:  
 Town: Postal code:  
 E-Mail:  
 Tel: Fax:

Order date:  
 Purchase order:  
 Number required:

Marked for:  
 Sunrise order no:

- = Standard  = Cost option
- = No cost option

**Max. user weight:** 125 kg  
**USER WEIGHT**  kg

Mandatory information to ensure correct configuration

**Part. No. : 84600001** **QPP01**

<b>Base Price</b>
<b>£2,195</b>



**CLIENT MEASUREMENTS** (all measurements without cushion)

**SEAT WIDTH SELECTION** - measured from outside frame tube to outside frame tube (cm)

28	30	32	34	36	38	40	42	44	46	48	50
<input type="checkbox"/>											

**SEAT DEPTH SELECTION** - measured in cm

32*	34*	36	38	40	42	44	46
<input type="checkbox"/>							
<b>Short Frame</b>				<b>Medium Frame</b>		<b>Long Frame</b>	

- \* not available with frame widths 34cm or greater
- \* not available with swing-away frame

**FRONT SEAT HEIGHT SELECTION:** in dependance of front castors and frame versions

**Front seat height in cm**

		Fork type (mm)		Front seat height in cm														
		40	41	42	43	44	45	46	47	48	49	50	51	52	53			
<b>Front castor</b>	<b>4"</b>	Hemi Version	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>												
		Carbotecture (111 mm)	<input type="checkbox"/>															
		Aluminium (123 mm) (1)		<input type="checkbox"/>														
		Aluminium (143 mm) (1)				<input type="checkbox"/>												
		Aluminium (168 mm) (1)							<input type="checkbox"/>									
	<b>5"</b>	Hemi Version	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>												
		Carbotecture (111 mm)			<input type="checkbox"/>													
		Aluminium (123 mm) (1)				<input type="checkbox"/>												
		Aluminium (143 mm) (1)					<input type="checkbox"/>											
	<b>6"</b>	Hemi Version			<input type="checkbox"/>													
		Carbotecture (111 mm)						<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>								
		Aluminium (123 mm) (1)							<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>						
Aluminium (143 mm) (1)									<input type="checkbox"/>									
<b>7"</b>	Hemi Version							<input type="checkbox"/>										
	Aluminium (168 mm) (1)								<input type="checkbox"/>									

Hemi version only with SA frame

**REAR SEAT HEIGHT SELECTION:** in dependance of rear wheel size

**Rear seat height in cm**

Rear wheels	40	41	42	43	44	45	46	47	48	49	50
<b>22"</b>	<input type="checkbox"/>	<input type="checkbox"/>		<input type="checkbox"/>	<input type="checkbox"/>		<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>		
<b>24"</b>				<input type="checkbox"/>	<input type="checkbox"/>		<input type="checkbox"/>				

The rear seat height can only be between 0 mm and 110 mm below the selected front seat height.

(1) compatible with power add on devices

						£	
<b>SEAT UPHOLSTERY AND CUSHIONS</b>							
<b>Seat slings</b>							
LFT020014	● Seat sling, nylon, 1 utility bag						
LFT020000	○ Without seat sling						
<b>Seat cushion, cover with zipper, black</b>							
LFT020101	□ Seat cushion, foam (medium-soft), cover with zipper, black, 3cm						£90
LFT020102	□ Seat cushion, foam (medium-soft), cover with zipper, black, 5cm						£90
LFT020103	□ Seat cushion, Latex (very soft), cover with zipper, black, 5cm						£90
<b>Jay Cushion - for complete Jay product range see Jay order forms</b>							
JAY020002	□ JAY Basic						£75
JAY020001	□ JAY Easy Visco C						£180
JAY020007	□ JAY Easy Fluid Cushion						£415
JAY020004	□ JAY LITE						£345
<b>Seat accessories</b>							
LFT090018	○ Lap belt with buckle, metal						
<b>BACKREST</b>							
<b>Backrest height range</b>	<b>25-30</b>	<b>30-35</b>	<b>35-40</b>	<b>40-45</b>	<b>45-47.5</b>		
Back height setting at (cm)	○ 25	○ 30	○ 35	○ 40	○ 45		
	○ 27.5	○ 32.5	○ 37.5	○ 42.5	○ 47.5		
	○ 30	○ 35	○ 40	○ 45			
LFT030001	● Backrest tube, aluminium						
LFT030101	□ Backrest, half-folding (to the back)					(from back height 45-47.5cm)	£90
LFT030031	□ Backrest, fold-down (to the front)					(from back height 30-47.5cm)	£90
LFT030012	○ Backrest, angle-adjustable						
	○ -12° (78 EU)	○ -8° (82 EU)		○ -4° (86 EU)			
	○ 0° (90 EU)	○ 4° (94 EU)		○ 8° (98 EU)	○ 12° (102 EU)		
<b>Push handles</b>							
LFT030201	● Push handles long						
LFT030203	□ Push handles, fold-down (1)					(only from back height 30cm; from 35cm with angle adj back)	£170
LFT030204	□ Push handles, height-adjustable (1) (2)					(only from back height 30cm)	£230
LFT030205	□ Push handles, height-adjustable for Jay backs (1) (2)					(only from back height 30cm)	£230
LFT030200	○ Without push handles (1)						
<b>Backrest upholstery (black)</b>							
LFT030316	○ EXO Backrest sling, synthetic fibre, padded, adjustable, 1 accessory bag, black (Evo Exo Version)						
	Colour, binding ● Black ○ Orange ○ Silver ○ Red ○ Blue						
LFT030317	□ EXO PRO Backrest sling, synth fibre, padded, adj, breathable, spacer fabric, 1 accessory bag, black (Evo Exo Version)						£195
	Colour, binding ● Black ○ Orange ○ Silver ○ Red ○ Blue						
LFT030300	○ Omit Back Upholstery						
<b>Backrest accessories</b>							
LFT030402	□ Stabilising bar, automatic folding (2)					(from sw 360 mm)	£190
LFT030411	□ Headrest, medium height, depth and angle-adjustable, 20 x 11cm						£270
LFT030410	□ Headrest, tall height, depth and angle-adjustable, 23 x 11cm						£270
(1)	Not in combination with backrest, half-folding (to the back)						
(2):	Not in combination backrest backrest, fold-down (to the front)						

				£
<b>JAY BACKS plus ACCESSORIES - for complete Jay product range see Jay order forms</b>				
<b>JAY 3 Back - Deep Contour</b>				
Upper Thoracic:	<input type="checkbox"/>	JAY030021	Medium 46cm	£870
JAY030004	<input type="checkbox"/>	Jay Flip Back Headrest (J3EUHDFH & J3EUHDLSP)		(only with JAY 3 backs) £345
JAY030006	<input type="checkbox"/>	Jay Std. Headrest (J3EUHDXH & J3EUHDLSP)		(only with JAY 3 backs) £235
<b>JAY J3 Carbon Back - Shallow Contour (with quick-release hardware)</b>				
Low Thoracic:	<input type="checkbox"/>	JAY060001	Short - 17cm	<input type="checkbox"/>
				JAY060002
				Tall - 24cm
				£820
Mid Thoracic:	<input type="checkbox"/>	JAY060003	Short - 30cm	<input type="checkbox"/>
				JAY060004
				Tall - 38cm
				£820
<b>JAY Zip Back</b>				
M (20 cm)	<input type="checkbox"/>	JAY070004	SW 28cm	£725
M (20 cm)	<input type="checkbox"/>	JAY070005	SW 30-34cm	£725
T (25cm)	<input type="checkbox"/>	JAY070007	SW 28cm	£750
T (25cm)	<input type="checkbox"/>	JAY070008	SW 30-34cm	£765
X (30cm)	<input type="checkbox"/>	JAY070010	SW 28cm	£780
X (30cm)	<input type="checkbox"/>	JAY070011	SW 30-34cm	£795
<b>ARMREST AND SIDEGUARDS</b>				
<b>Armrest</b>				(choose only one option, not in combination with sideguards; only with 0° and 1° camber)
LFT040050	<input type="checkbox"/>	<input type="radio"/>	Single-post sideguard, armpad (25cm), height-adjustable, removable	£195
LFT040051	<input type="checkbox"/>	<input type="radio"/>	Single-post sideguard, armpad (36cm), height-adjustable, removable	£195
LFT040052	<input type="checkbox"/>	<input type="radio"/>	Single-post sideguard, armpad (33cm), depth adjustable, height-adjustable, removable	(3)
<b>Desk armrest with long and short pad</b>				(1)
LFT040001	<input type="checkbox"/>	Desk sideguard, armpad (25cm), fixed, flip-up, removable		£215
LFT040006	<input type="checkbox"/>	Desk sideguard, armpad (36cm), fixed, flip-up, removable		£215
LFT040002	<input type="checkbox"/>	Desk sideguard, armpad (25cm), height-adjustable, flip-up, removable		£265
LFT040007	<input type="checkbox"/>	Desk sideguard, armpad (36cm), height-adjustable, flip-up, removable		£265
LFT040005	<input type="checkbox"/>	Desk sideguard, flip-up, removable		£215
LFT040020	<input type="checkbox"/>	Release assistance for tetraplegics		£40
<b>Sideguards</b>				(not in combination with armrest)
LFT040107	<input type="checkbox"/>	<input type="radio"/>	Sideguard, composite, removable, 21cm high (31cm long), black	(4)
LFT040103	<input type="checkbox"/>	<input type="radio"/>	Sideguard, aluminium, fender, cold weather protection (in frame colour)	(2) (available with safari brake)
LFT040102	<input type="checkbox"/>	<input type="checkbox"/>	Sideguard, aluminium, cold weather protection (in frame colour)	(2) (not with safari brake) £220
LFT040055	<input type="checkbox"/>	Armpad, height-adjustable, short		(not with safari brake) £220
LFT040104	<input type="checkbox"/>	<input type="checkbox"/>	Sideguard, carbon fibre, fender, cold weather protection (black)	(2) (not with safari brake) £410
<b>TRAYS</b>				
LFT090111	<input type="checkbox"/>	Tray, wood, push-fit, SW 22-36cm		only with LFT040055 £185
LFT090112	<input type="checkbox"/>	Tray, plastic, swing-away, left-mount		(only with height adj. armrests with pad) £185
LFT090113	<input type="checkbox"/>	Tray, plastic, swing-away, right-mount		(only with height adj. armrests with pad) £185
<b>FRAME</b>				
LFT010010	<input type="checkbox"/>	<input type="radio"/>	Frame FF 75° (integrated footrest)	(from LLL 26cm)
LFT010011	<input type="checkbox"/>	<input type="radio"/>	Frame FF 85° (integrated footrest)	(from LLL 26cm)
LFT010023	<input type="checkbox"/>	<input type="radio"/>	Frame FF 75° (integrated footrest), abducted (25 mm on each side)	(not for SW 48/50; from LLL 26cm, Max User Weight 100kg)
LFT010016	<input type="checkbox"/>	<input type="checkbox"/>	Frame FF 85° (integrated footrest), inset (20 mm on each side), active	(from SW 34; from LLL 26cm) £195
LFT010015	<input type="checkbox"/>	<input type="radio"/>	Frame SA (removable footrest)	(not for SD 32/34)
LFT090100	<input checked="" type="checkbox"/>	Tie down brackets for transportation		Std
(1) not in combination with backrest, half-folding (to the back), fold down back only from BH 375 mm only 1° & 2° camber				
(2) Effective seat width (measured between the sideguards) is seat width + 20 mm				
(3) with abduction only from SW 34				
(4) COG 83 not possible with abducted frame; not with safari brake				

			£
<b>CENTRE OF GRAVITY (COG)</b> measured from front of back post; choose only one option			
COG23	<input type="radio"/>	Centre of gravity 2.3cm	
COG43	<input type="radio"/>	Centre of gravity 4.3cm	
COG63	<input type="radio"/>	Centre of gravity 6.3cm	
COG83	<input type="radio"/>	Centre of gravity 8.3cm	
<b>FRAME ACCESSORIES</b>			
LFT090010	<input type="checkbox"/>	Transit wheels <i>(not with crutch holder, tip assist and anti tip)</i>	£255
<b>Tip assist (step tube) &amp; crutch holder</b>			
LFT090001	<input type="checkbox"/>	Tip assist, left	£35
LFT090002	<input type="checkbox"/>	Tip assist, right	£35
LFT090003	<input type="checkbox"/>	Tip assist, pair <i>(choose either tip assist or anti- tip tubes, not both)</i>	£70
LFT090009	<input type="checkbox"/>	Crutch holder, loop <i>(only in combination with tip assist)</i>	£90
<b>Anti-tip tubes, swing-away</b>			
LFT090036	<input type="radio"/>	Anti-tip tubes, aluminium, lightweight, swing-away (forwards), pair <i>(not with tip assist)</i>	£210
LFT090004	<input type="checkbox"/>	Anti-tip tubes, swing-away, left	£115
LFT090005	<input type="checkbox"/>	Anti-tip tubes, swing-away, right	£115
LFT090006	<input type="checkbox"/>	Anti-tip tubes, swing-away, pair <i>(choose either tip assist or anti- tip tubes, not both)</i>	£235
<b>PREPARATION FOR ATTITUDE ADD-ON BIKE</b> - only with FF-Frames			
LFT090211	<input type="checkbox"/>	without ATTITUDE order <i>(docking pins assembled)</i>	£130
LFT090212	<input type="radio"/>	with ATTITUDE order <i>(docking pins assembled; ATTITUDE order seperately)</i>	
<b>HANGERS &amp; FOOTPLATES</b>			
<b>Aluminium hanger</b>			
LFT050002	<input type="radio"/>	Hanger angle 70°	
LFT050001	<input type="radio"/>	Hanger angle 80°	
LFT050041	<input type="radio"/>	Hanger angle 90° <i>(not with 6", 7" castor)</i>	
LFT050006	<input type="radio"/>	Hanger, aluminium, swing-away, LH	LFT050007 <input type="radio"/> Hanger, aluminium, swing-away, RH
<i>pick both</i>	<input type="radio"/>	Hanger, aluminium, swing-away, pair	
<b>Composite hanger</b>			
LFT050002	<input type="radio"/>	Hanger angle 70°	
LFT050001	<input type="radio"/>	Hanger angle 80°	
LFT050003	<input type="radio"/>	Hanger, composite, swing-away, LH	LFT050004 <input type="radio"/> Hanger, composite, swing-away, RH
<i>pick both</i>	<input type="radio"/>	Hanger, composite, swing-away, pair	
<b>Elevating Legrest (ELR) (1)</b> <i>(not with footboards, platform; only SA frame)</i>			
LFT050009	<input type="checkbox"/>	Hanger, aluminium, flip-up (0° - 90°), swing-away, LH	£170
LFT050010	<input type="checkbox"/>	Hanger, aluminium, flip-up (0° - 90°), swing-away, RH	£170
<i>pick both</i>	<input type="checkbox"/>	Hanger, aluminium, flip-up (0° - 90°), swing-away, pair <b>(1)</b>	£340
<b>Amputee Support Pad</b>			
LFT050014	<input type="checkbox"/>	Amputee support, left	£170
LFT050013	<input type="checkbox"/>	Amputee support, right	£170
LFT050012	<input type="checkbox"/>	Amputee support, pair	£340
<b>Footboards (colour black)</b>			
LFT050020	<input checked="" type="radio"/>	Footboard, divided, composite, angle-adjustable, flip-up (sideways), heel loop, black	
LFT050050	<input type="radio"/>	Footboard, divided, aluminium, angle-adjustable, flip-up (sideways), heel loop, black	
LFT050102	<input type="radio"/>	Footboard, platform, aluminium, angle-adjustable, flip-up (sideways), heel loop, black	
	LFT050126	<input type="checkbox"/> Swing-away mechanism for platform footboard <i>(only from LLL 380 mm)</i>	£195
	LFT050062	<input type="checkbox"/> Locking system, platform footplate, aluminium	£70
LFT050131	<input type="radio"/>	Footboard, platform, lightweight, plastic, angle & depth-adjustable flip up (sideways), calf strap, black	
LFT050132	<input type="checkbox"/>	Footboard,platform,lightweight,carbon fibre,angle &depth-adjustable flip up (sideways),calf strap,black	£405
LFT050127	<input type="checkbox"/>	Adjustable side protection, to create a side-stop on the footboard	£25
<b>(1)</b> not in combination with footboard platform; not with SW 28-36			

**LOWER LEG LENGTH (LLL) - LLL has to be 2.5cm lower than front seat height; measurements in cm**

- 22     23     24     25     26     27     28     29     30     31  
 32     33     34     35     36     37     38     39     40     41  
 42     43     44     45     46     47     48     49     50

SA Frame: LLL 22-29 from SW34 with divided footboard / from SW32 with platform footboard

FF Frame: LLL 26-33 from SW34 (SW30 with abduction) with divided footboard/from SW32 (SW28 with abduction) with platform footboard

**FORK***(compare with your choice on the first page)*

LFT080041	<input type="radio"/>	Castor fork, Aluminium, Black	
		LFT080042 <input type="radio"/> Colour Cast. Fork Silver	
		LFT080043 <input type="radio"/> Colour Cast. Fork Orange	
LFT080013	<input type="radio"/>	Castor fork, Carbotecture © <input checked="" type="radio"/> Black <input type="radio"/> Orange <input type="radio"/> Silver	<i>(only 111 mm)</i>
LFT080101	<input type="checkbox"/>	Expansion kit, designed for hemiplegic users	£145

**CASTOR***(compare with your choice on the first page)*

LFT080301	<input type="radio"/>	Castor, solid (hard), rim: composite, black	
LFT080302	<input type="radio"/>	Castor, solid (soft), rim: composite, black	<i>(only for 5" and 6")</i>
LFT080305	<input type="radio"/>	Castor, pneumatic, rim: composite, black	<i>(only for 4" and 5")</i>
LFT080304	<input type="checkbox"/>	Castor, solid (hard), rim: aluminium	<i>(only for 4" and 5")</i> £100
LFT080306	<input type="checkbox"/>	Castor, solid (soft), rim: aluminium, silver	<i>(only for 4" and 5")</i> £100
		Colour, rim:    LFT080311 <input type="radio"/> Black    LFT080315 <input type="radio"/> Blue	
		LFT080313 <input type="radio"/> Orange    LFT080314 <input type="radio"/> Red	

**REAR WHEEL***(compare with your choice on the first page; choose only one wheel type)*

LFT070001	<input type="radio"/>	Rear wheel, Universal, 36 spokes (silver), cross-wise spokes	
LFT070002	<input type="radio"/>	Rear wheel, Design, 36 spokes (silver), straight spokes	<i>(only 24")</i>
LFT070003	<input type="radio"/>	Rear wheel, Lightweight, 24 spokes (silver), straight spokes	
LFT070004	<input type="checkbox"/>	Rear wheel, Proton, 24 spokes (black), straight spokes	<i>(Max User Weight 100kg)</i> £625
LFT070005	<input type="checkbox"/>	Rear wheel, Spinergy, 18 spokes (black), straight spokes, hub: silver	<i>(only 24")</i> £580
LFT070008	<input type="checkbox"/>	Rear wheel, drum brakes, 36 spokes (silver), cross-wise spokes <i>(22" &amp; 24"; only with stainless steel handrim, COG 43 mm; 0° camber) (3)</i>	£470
LFT070401	<input type="checkbox"/>	Rear wheel, one-hand operation, 36 spokes, cross-spoked, left actuation	<i>(not with camber)</i> £870
LFT070402	<input type="checkbox"/>	Rear wheel, one-hand operation, 36 spokes, cross-spoked, right actuation	<i>(not with camber)</i> £870
LFT070000	<input type="radio"/>	Without rear wheel	

**TYRES***(choose only one tyre type)*

LFT070107	<input checked="" type="radio"/>	Tyre, pneumatic, slick, Schwalbe One	<i>(only 24")</i>
LFT070102	<input type="radio"/>	Tyre, puncture proof, tread	<i>(not with lightweight wheels)</i>
LFT070101	<input type="radio"/>	Tyre, pneumatic, fine profile, Schwalbe Right Run	
LFT070104	<input type="checkbox"/>	Tyre, pneumatic, puncture proof, tread, Schwalbe Marathon Plus Evolution	£120

**HANDRIM**

LFT070202	<input type="radio"/>	Handrims, aluminium, anodised, silver	
LFT070201	<input type="radio"/>	Handrims, aluminium, powder-coated, clear	
LFT070204	<input type="radio"/>	Handrims, stainless steel, silver	
LFT070203	<input type="checkbox"/>	Handrims, aluminium, powder-coated, coloured	£80
LFT070207	<input type="checkbox"/>	Handrims, Supergrip ®	<i>(22" only with universal and drum brake wheels)</i> £295
LFT070208	<input type="checkbox"/>	Handrims, Maxgripp ®	<i>(only handrim mounting wide)</i> £295
LFT070212	<input type="checkbox"/>	Handrims, Ellipse 3R, black	<i>(ergonomic handrims with rubber-strip, 24" and 25" only)</i> £235
LFT070200	<input type="radio"/>	Without handrims	

**(1)** not in combination with 7" castors, Carbotecture fork; SA frame only**(2)** only with 0-4° camber**(3)** RSH for 22":40-45cm; for 24": 42-47cm; axle plate silver

										£
<b>HANDRIM MOUNTING</b>										
LFT070317	<input type="radio"/> Handrim mount, narrow (17 mm) <span style="float: right;"><i>(not with Maxgrepp)</i></span>									
LFT070327	<input type="radio"/> Handrim mount, wide (27 mm)									
LFT070307	<input type="radio"/> Handrim mount, super narrow (7 mm) <span style="float: right;"><i>(Only with 24" stainless steel)</i></span>									
<b>WHEEL CAMBER</b> <span style="float: right;"><i>(choose only one option; Armrests only with 0° &amp; 1°)</i></span>										
LFT070600	<input type="radio"/> Camber 0° <span style="float: right;"><i>(no camber brackets provided)</i></span>									
LFT070601	<input type="radio"/> Camber 1°		LFT070603	<input type="radio"/> Camber 3°		LFT070605	<input type="radio"/> Camber 5°			
LFT070602	<input checked="" type="radio"/> Camber 2°		LFT070604	<input type="radio"/> Camber 4°		LFT070606	<input type="radio"/> Camber 6°			
<b>SPOKE GUARDS</b> <span style="float: right;"><i>(not with spinergy wheels)</i></span>										
LFT070831	<input type="checkbox"/> Transparent		LFT070813	<input type="checkbox"/> Glaktica		LFT070828	<input type="checkbox"/> Space			£160
LFT070823	<input type="checkbox"/> Princess		LFT070821	<input type="checkbox"/> Pirates						£160
LFT070816	<input type="checkbox"/> Happy Lions		LFT070814	<input type="checkbox"/> Graffiti						£160
LFT070819	<input type="checkbox"/> Monkeys		LFT070802	<input type="checkbox"/> Basketball						£160
LFT070820	<input type="checkbox"/> Oriental		LFT070812	<input type="checkbox"/> Flying Water						£160
LFT070829	<input type="checkbox"/> Superhero		LFT070833	<input type="checkbox"/> Zebra						£160
LFT070817	<input type="checkbox"/> Hot Wheels		LFT070818	<input type="checkbox"/> Jungle Fever						£160
LFT070822	<input type="checkbox"/> Pool		LFT070811	<input type="checkbox"/> Feather						£160
LFT070827	<input type="checkbox"/> Silk		LFT070806	<input type="checkbox"/> Chrome						£160
LFT070830	<input type="checkbox"/> Tiger		LFT070803	<input type="checkbox"/> Basketball 2						£160
<b>WHEEL LOCKS</b> <span style="float: right;"><i>(choose only one option)</i></span>										
LFT060001	<input type="radio"/> Standard brake									
LFT060003	<input type="radio"/> Compact brake									
LFT060013	<input type="radio"/> Push to lock brake, integrated in sideguard, Safari									
LFT060010	<input type="checkbox"/> Brake lever extension, knee-lever brake, fold-down									£50
LFT060002	<input type="radio"/> Knee-lever brake									
LFT060016	<input type="checkbox"/> Push to lock brake, one-arm, right actuation									£210
LFT060015	<input type="checkbox"/> Push to lock brake, one-arm, left actuation									£210
<b>ACCESSORIES</b>										
LFT090000	<input checked="" type="radio"/> Tool kit									
LFT090025	<input type="checkbox"/> Air pump									£20
LFT090024	<input type="checkbox"/> High pressure pump (0 - 11 bar)									£95
LFT090030	<input type="checkbox"/> Back pack									£80
LFT090026	<input type="checkbox"/> Mobile phone pocket									£60
<b>COLOURS</b> <sup>(1)</sup> Please consult our colour catalogue to review your colour choice										
	Frame	Handrim	Hub/rim £230	Fork £75		Frame	Handrim	Hub/rim £230	Fork £75	
02 Brilliant Silver	<input type="radio"/>	<input type="checkbox"/>		<input type="checkbox"/>	29 Black	<input type="radio"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	
03 Silver Matt	<input type="radio"/>	<input type="checkbox"/>		<input type="checkbox"/>	31 Matt Black	<input type="radio"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	
06 Sun Yellow	<input type="radio"/>	<input type="checkbox"/>		<input type="checkbox"/>	32 Signal White	<input type="radio"/>	<input type="checkbox"/>		<input type="checkbox"/>	
14 Blue	<input type="radio"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	33 Pink	<input type="radio"/>	<input type="checkbox"/>		<input type="checkbox"/>	
15 Brilliant Marine	<input type="radio"/>	<input type="checkbox"/>		<input type="checkbox"/>	34 Dark Chocolate Metallic	<input type="radio"/>	<input type="checkbox"/>		<input type="checkbox"/>	
16 Matt Marine	<input type="radio"/>	<input type="checkbox"/>		<input type="checkbox"/>	35 Skyblue	<input type="radio"/>	<input type="checkbox"/>		<input type="checkbox"/>	
19 Brilliant Purple	<input type="radio"/>	<input type="checkbox"/>		<input type="checkbox"/>	36 Cyan-Metallic	<input type="radio"/>	<input type="checkbox"/>		<input type="checkbox"/>	
21 Red Ruby	<input type="radio"/>	<input type="checkbox"/>		<input type="checkbox"/>	37 Marsred (Orange)	<input type="radio"/>	<input type="checkbox"/>		<input type="checkbox"/>	
22 Brilliant Red	<input type="radio"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	40 Kiwi	<input type="radio"/>	<input type="checkbox"/>		<input type="checkbox"/>	
23 Matt Red	<input type="radio"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	41 Flamingo-Glitter	<input type="radio"/>	<input type="checkbox"/>		<input type="checkbox"/>	
24 Tutti Frutti Orange	<input type="radio"/>	<input type="checkbox"/>		<input type="checkbox"/>	42 Matt-Anthracite	<input type="radio"/>	<input type="checkbox"/>		<input type="checkbox"/>	
27 Graphite Metallic	<input type="radio"/>	<input type="checkbox"/>		<input type="checkbox"/>	46 Mamba Green	<input type="radio"/>	<input type="checkbox"/>		<input type="checkbox"/>	
<i>(1) Handrim: only with handrim colour powder coated; Hub/rim: only for universal wheel; Fork: not with Carbotecture</i>										

£

**STICKER COLOURS**

LFT100550	<input checked="" type="radio"/> Sticker colour - black	
LFT100551	<input type="radio"/> Sticker colour - silver	

**ANODISED COLOUR AXLE PLATE**

(not in combination with one hand operation)

LFT100580	<input checked="" type="radio"/> Axle plate anodised , black	
LFT100581	<input type="radio"/> Axle plate anodised , silver	
LFT100582	<input type="radio"/> Axle plate anodised, orange	

**24 " LIGHT WEIGHT WHEEL COLOURS**

	<input checked="" type="radio"/> Black	
LFT100336	<input type="checkbox"/> Anodised, orange	£120
LFT100337	<input type="checkbox"/> Anodised, blue	£120
LFT100338	<input type="checkbox"/> Anodised, red	£120

**Comments for special build** - please apply for a quote; lead time conditions my change.

<b>Quickie Life T - Basic Price</b>	<b>£2,195</b>
Option price	
Packaging and transport	
Exclusive VAT	
VAT.	
<b>Total Price</b>	

**Please follow this measuring method:**

The following method is recommended on the basis of extensive practical experience. Deviations are possible in individual instances.

**Seat width (SW):**

- a) For the user: Allow space between sideguard and hip on each side of the chair (approx. 10 - 20 mm altogether)
- b) On the wheelchair: From outside frame tube to outside frame tube.

**Seat depth (SD):**

- a) For the user: When sitting upright: leave space for 2 - 3 fingers between front edge of the seat upholstery and back of knee. (approx. 30 - 60 mm). Care must be taken with adj. upholstery.
- b) On the wheelchair: From back tube to front of seat upholstery.

**Front seat height:**

- a) For the user: When sitting with 90° kneeangle: From back of knee to bottom of shoe sole.
- b) On the wheelchair: From front of seat tube to the ground. 125 kg  
(min. 30 mm longer than user measurement)  
Remember to include the cushion height.

**Rear seat height:**

- b) On the wheelchair: From back of seat tube to the ground.  
Note: Max. seat height difference rear to front is 130 mm

**Back height:**

- a) For the user: From the seat, including cushion, to be no more than 10 - 20 mm below the lower shoulder blade.
- b) On the wheelchair: From back of seat tube to upper edge of backrest upholstery.  
Note: The user's ability to sit upright has a major influence on the back height and it can therefore differ considerably

**Back angle:**

Important: The back height depends on the degree of disability.  
0° is the angle 90° to the ground. Negative angles are < 90°  
and positive angles are > 90°.

**Technical Data Overview - Life T**

<b>Max. overall height:</b>	1035 mm		
<b>Max. overall length:</b>	1075 mm		
<b>Max. turning circle:</b>	860 mm		
<b>Max. safe slope:</b>	10°		
<b>Wheelchair weight:</b>	from 10 kg		
<b>Dismantled weight:</b>	from 7 kg		
<b>Weight of heaviest part:</b>	2,1 kg		
<b>Max. overall width:</b>	<b>Camber</b>	<b>22" rear wheel</b>	<b>24" rear wheel</b>
	For 0° camber:	SW + 200 mm	SW + 200 mm
	For 3° camber:	SW + 230 mm	SW + 240 mm
	For 6° camber:	SW + 250 mm	SW + 260 mm

Product is **crash tested according to ISO 7176-19** and therefore approved for transportation in a vehicle  
For transport without a user the chair must be safe clamped in the vehicle

All changes in the configuration might change the seat heights.  
Manufacturing tolerances of +/- 10 mm

## BACKREST



Backrest tube,  
aluminium



Backrest, half folding (to  
the back)



Backrest, angle  
adjustable



Backrest, fold down (to  
the front)



Push handles long



Push handles,  
fold-down



Push handles,  
height-adjustable



Without push  
handles,



Stabilizer bar,  
automatic folding



Headrest, height, depth and  
angle adjustable



EXO backrest sling



EXO PRO backrest sling

## JAY BACKS



Jay 3 Back  
Mid Contour



Jay Zip Back

## ARMREST AND SIDEGUARDS



Single post height  
adjustable



Single-post sideguard,  
armrest, height-adjustable,  
remove



Desk sideguard, armrest,  
fixed, flip-up, remove



Desk sideguard, armrest,  
height-adjustable, flip-up,  
remove



Desk sideguard, flip-up,  
remove



Composite sideguard,  
removeable (black)



Aluminium sideguard, wet  
weatherprotection (in  
frame colour)



Aluminium sideguard,  
fender, wet  
weatherprotection (in  
frame colour)



Carbon fibre sideguard,  
fender, wet  
weatherprotection (black)

## FRAME



Frame FF 105 °  
(integrated footrest)



Frame FF 95 °  
(integrated footrest)



Frame FF 105°  
(integrated footrest),  
abducted (20 mm on  
each side)



Frame FF 95 °  
(integrated footrest),  
20 mm inset, active



Frame SA  
(removable footrest)

## FRAME ACCESSORIES



Transit wheels



Tip assist



Crutch holder,  
loop



Anti-tip tubes, swing-  
away, left

## HANGERS AND FOOTPLATES



Aluminium hanger, 90°



Composite hanger  
110°/100°



Elevating legrest  
(ELR)



Amputee Support Pad



Footboard, divided,  
composite, angle-  
adjustable, flip-up  
(sideways), heel loop,  
black



Footboard, divided,  
aluminium, angle-  
adjustable, flip-up  
(sideways), heel loop,  
black



Footboard, platform,  
lightweight, plastic,  
angle & depth-  
adjustable flip up  
(sideways), calf strap,  
black



Footboard, platform,  
lightweight, plastic,  
angle & depth-adjustable  
flip up (sideways), calf  
strap



Footboard, angle and depth  
adjustable, flip up, aluminium with  
locking system

## FORKS



Fork Aluminium,  
black, orange & silver



Expansion kit, designed for  
hemiplegic users

## REAR WHEELS



Rear wheel, Universal,  
36 spokes, cross-wise  
spokes



Rear wheel, Design, 36  
spokes, straight spokes



Rear wheel, Lightweight,  
24 spokes, straight  
spokes



Rear wheel, Proton,  
24 spokes, straight  
spokes



Rear wheel,  
Spinergy, 18  
spokes (black),  
straight spokes;  
hub: silver



Rear wheel,  
drum brakes, 36  
spokes, cross-  
wise spokes

## TYRES



Tyre, puncture proof, tread



Tyre, pneumatic, fine  
profile, Schwalbe Right  
Run



Tyre, pneumatic, puncture proof,  
tread, Schwalbe Marathon Plus  
Evolution



Tyre, pneumatic, slick,  
Schwalbe One

## HANDRIMS



Supergrip



Max Grepp



Ellipse 3R Handrim

## WHEEL LOCKS



Standard brake



Knee-lever brake



Compact brake



Push to lock brake,  
integrated in sideguard,  
Safari



Push to lock brake, one-  
arm, right actuation



Brake lever extension,  
knee-lever brake, fold-  
down

## ACCESSORIES



Back pack



High pressure pump



Air pump



Mobile phone  
pocket

**ANODIZED COLOUR AXLE PLATE**



Axle plate  
anodized, orange



Axle plate  
anodized, silver



Axle plate  
anodized, black

**24 " LIGHT WEIGHT WHEEL COLOURS**



*Lightweight wheel, anodized orange*



*Lightweight wheel, anodized blue*



*Lightweight wheel, anodized red*