

Q700® M HD



ORDER FORM

Billing Address Or Customer-ID		Shipping Address		<input type="checkbox"/> Identical To Billing Address
C.-ID:		C.-ID:		
Name:		Name:		
Street:		Street:		
Town/PC:		Town/PC:		
E-Mail:		E-Mail:		
Phone:		Phone:		
Order Information		User Weight		
<input type="checkbox"/> Order <input type="checkbox"/> Calculation <input type="checkbox"/> Built-4-Me®		User Weight Capacity: 250 kg		
Purchase Order No.:		User Weight:		kg
Art. No.: 2250070M-300 Reimbursement Code: xxx				
Order Date:	Day	Month	Year	

Seating Specifications

Typical WC Weight	Seat Type	SEDEO PRO
Without Batteries 140 kg	Seat Size Range	Width 500 - 710 mm, Depth 460 - 560 mm
With Batteries 165 kg	Seat To Floor Height Range	415 mm - 490 mm

MWD Base Specifications

User Weight Capacity	250 kg	Drive Wheels	14" Drive Wheels (Black)
Base Size	Width Min. 620 mm	Tires / Castors	6" Front / 6" Rear (Black)
Drive Control	R-Net	Charger	10 A Dual Mode Charger
Driving Range	Up To 25 km	Battery Size	60 And 80 Ah Batteries
Drive Motors	4-Pole 6/10kph	Colour Accents	7 Colour Accents
Ground Clearance	70 mm	Frame Colour	Black

Q700 M HD Drive Bases (Max User Weight 250kg)		£
<input type="checkbox"/> SHP000080	Q700 M HD (MWD Base Chair)	8450

Standard Seat Size		£
Seat Cushion Width - Omit or SEDEO Cushions		
<input type="checkbox"/> SW50	19.5" (500 mm) (6)	
<input type="checkbox"/> SW54	21" (540 mm) (6)	
<input type="checkbox"/> SW56	22" (560 mm) (6)	
<input type="checkbox"/> SW61	24" (610 mm) (6)	535
<input type="checkbox"/> SW66	26" (660 mm)(6)	535
Select SW66 & Single Post Armrest	28" (710 mm)	535
Seat Depth		
<input type="checkbox"/> SD46	18" (460 mm)	
<input type="checkbox"/> SD48	19" (480 mm)	
<input type="checkbox"/> SD50	19.5" (500 mm)	
<input type="checkbox"/> SD52	20.5" (520 mm)	
<input type="checkbox"/> SD54	21" (540 mm) (1)	
<input type="checkbox"/> SD56	22" (560 mm) (1)	

Seat Cushion Type		£
<input type="checkbox"/> SHP020110	Omit Seat Cushion	0

Seat Rail		£
<input type="checkbox"/> SHP020131	Seat Rail, Long HD, For Seat Depth 460 - 560 mm	0

SEDEO Seat Cushion		£
<input type="checkbox"/> SHP020073	SEDEO, Seat Cushion, Profiled, Height 70 mm	160
<input type="checkbox"/> SHP020077	SEDEO, Seat Cushion, Form, Height 110 mm	260

SEDEO Seat Cushion Cover		£
<input type="checkbox"/> SHP020518	Seat Cushion Cover, Incontinence Protection	0

SEDEO Seat Boards		£
<input type="checkbox"/> SHP020082	SEDEO, Seat Board, Flat, Depth-Adjustable, Cushions Guide At Rear	130

JAY Seat Cushions		£
<input type="checkbox"/> JAY020011	JAY Balance, Fluid, Microclimatic Cover, Until Seat Width 500 mm	655
<input type="checkbox"/> JAY020013	JAY Balance, Fluid, Microclimatic Cover, Until Seat Width 560 mm	703

Seat Heights (2)		£
<input type="checkbox"/> SH415	415 mm (3)	0
<input type="checkbox"/> SH44	440 mm	0
<input type="checkbox"/> SH465	465 mm	0
<input type="checkbox"/> SH49	490 mm	0

Powered Seating Options		£
<input type="checkbox"/> SHP130011	Fixed Tilt Module (Adjustable 0°, 3°, 6°, 9°, Factory Setting 6°)	0
<input type="checkbox"/> SHP130004	Posterior Tilt, Powered, 0° - 30°	1340
<input type="checkbox"/> SHP130012	Posterior Tilt, Powered, -5° - 25° (4)	1500
<input type="checkbox"/> SHP130007	Posterior Tilt, Powered, 0° - 50° (5)	1630
<input type="checkbox"/> SHP130005	Posterior Tilt, Powered, -5° - 45° (4)(5)	1790

(1) Only with JAY cushion widths of 560mm upwards

(2) Measured from front of seat pan to floor. Obtain quote to verify seat to floor heights

(3) Only with 60Ah batteries

(4) Not with seat height 415 or 440mm

(5) Not with SHP030215 Backrest powered recline

(6) Can be made 2" wider via Single Post Armrest SHP040021 or SHP040019

(7) Only with 26" seat pan and Single Post Armrest SHP040021 or SHP040019

Backrest Frame Type			£
<input type="checkbox"/>	SHP030083	Backrest, Tool Angle Adjustable, 90° - 120°	0
<input type="checkbox"/>	SHP030215	Backrest, HD Powered Recline, 90° - 130°	750
Backrest Frame Option			£
<input type="checkbox"/>	SHP030098	Backrest, Flip Down, Quick Release	185
SEDEO Pro Backrest Size			£
Back Cushion Width		Back Cushion Height	
<input type="checkbox"/>	BW50	19.5" (500 mm), effective back width 460-	<input type="checkbox"/> BH50 20" (50cm for effective back height 50 - 56cm)
<input type="checkbox"/>	BW54	21" (540 mm), effective back width 500-540mm	<input type="checkbox"/> BH56 22" (56cm for effective back height 56 - 62cm)
<input type="checkbox"/>	BW56	22" (560 mm), effective back width 520-560mm	
<input type="checkbox"/>	BW61	24" (610 mm), effective back width 560-610mm	
<input type="checkbox"/>	BW66	26" (660 mm), effective back width 610-660mm	
SEDEO Pro Back Cushion			£
<input type="checkbox"/>	SHP030109	SEDEO PRO, Back Cushion, Profiled, Adjustable Straps	0
<input type="checkbox"/>	SHP030117	SEDEO Pro, Back Cushion, Flexible, Extra Deep Contour, Active, Adjustable Straps	475
SEDEO Pro Back Cushion Cover			£
<input type="checkbox"/>	SHP030108	Back Cushion Cover, Comfort	0
Headrest Receiver			£
<input type="checkbox"/>	SHP030120	Headrest Receiver, Backrest Frame Mount	0
<input type="checkbox"/>	SHP030119	Headrest Receiver, Back Cushion Mount (4)	0
Headrest Pad			£
<input type="checkbox"/>	SHP030430	SEDEO, Headrest Pad, Small, (LxH) 240 x 150 mm (1)	95
<input type="checkbox"/>	SHP030431	SEDEO, Headrest Pad, Small, With Cover, (LxH) 240 x 150 mm (1)	110
<input type="checkbox"/>	SHP030432	SEDEO, Headrest Pad, Flexible, (LxH) 400 x 110 mm (1)	170
<input type="checkbox"/>	WHT030402	Whitmyer, Headrest Pad, Plush, (LxH) 254 x 127 mm (2)	160
<input type="checkbox"/>	WHT030404	Whitmyer, Headrest Pad, Adjust-A-Plush, Small, (LxH) 355 x 76 mm (2)	260
Headrest Stem (3)			£
<input type="checkbox"/>	SHP030433	SEDEO, Headrest Mount	140
<input type="checkbox"/>	WHT030408	Whitmyer, Headrest Mount, LINX ²	175
<input type="checkbox"/>	WHT030409	Whitmyer, Headrest Mount, AXYS	145
<input type="checkbox"/>	WHT030411	Whitmyer, Headrest Mount, COBRA XTRA Flip Back	220
Push Handles			£
<input type="checkbox"/>	SHP030204	Push Handles, Height Adjustable	210

(1) Only use with SEDEO headrest mount

(2) Only use with Whitmyer headrest mount

(3) Choose a Whitmyer headrest mount and headrest receiver in case you would like to add a head control in the special controls section

(4) Only with SHP030109 SEDEO Pro, Back Cushion, Profiled

Armrest Types			£
<input type="checkbox"/>	SHP040021	Single Post Armrest, Removeable (Select for 2" wider seat width) (1) (2)	0
<input type="checkbox"/>	SHP040019	Single Post Armrest, Flip Back (Select for 2" wider seat width) (1) (2)	220
<input type="checkbox"/>	SHP040070	Flip Back Armrest	235
<input type="checkbox"/>	SHP040065	HD+ Flip Back Armrest (Reccommended for users over 205kg)	320

Side Guards			£
Side Guards for Single Post Flup Back or Removable			
<input type="checkbox"/>	SHP040033	Side Guard (3)	0
<input type="checkbox"/>	SHP040034	Side Guard, Integral Foam (4)	0
<input type="checkbox"/>	SHP040035	Side Guard, Wedge Shaped, Padded (5)	275
Side Guards for Flip Back Armrests			
<input type="checkbox"/>	SHP040036	Side Guard , Mounted on to Side Rails (6)	0

Arm Pads			£
<input type="checkbox"/>	SHP040023	Arm Pad, Integral Foam, 300 mm, Right	0
<input type="checkbox"/>	SHP040024	Arm Pad, Integral Foam, 300 mm, Left	0
<input type="checkbox"/>	SHP040025	Arm Pad, Integral Foam, 400 mm, Right	0
<input type="checkbox"/>	SHP040026	Arm Pad, Integral Foam, 400 mm, Left	0
<input type="checkbox"/>	SHP040027	Arm Pad, Integral Foam, 450 mm, Right	0
<input type="checkbox"/>	SHP040028	Arm Pad, Integral Foam, 450 mm, Left	0
<input type="checkbox"/>	SHP040204	Arm Pad, Comfort, Padded, PU Coated, 400 mm, Right (7)	75
<input type="checkbox"/>	SHP040203	Arm Pad, Comfort, Padded, PU Coated, 400 mm, Left (7)	75
<input type="checkbox"/>	SHP040031	Arm Pad, Hemiplegia, Integral Foam, Right (8)	275
<input type="checkbox"/>	SHP040032	Arm Pad, Hemiplegia, Integral Foam, Left (9)	275

Thigh Supports & Pommel			£
<input type="checkbox"/>	SHP020342	Knee Support, Lateral, 120 x 120 mm, Swing Away, Right (10)	185
<input type="checkbox"/>	SHP020341	Knee Support, Lateral, 120 x 120 mm, Swing Away, Left (10)	185
<input type="checkbox"/>	SHP020113	Thigh Support, Lateral, Small, 120 x 120 mm, Swing Away, Right (11)	185
<input type="checkbox"/>	SHP020114	Thigh Support, Lateral, Small, 120 x 120 mm, Swing Away, Left (11)	185
<input type="checkbox"/>	SHP020121	Thigh Support, Lateral, Large, 180 x 120 mm, Swing Away, Right (11)	185
<input type="checkbox"/>	SHP020122	Thigh Support, Lateral, Large, 180 x 120 mm, Swing Away, Left (11)	185
<input type="checkbox"/>	SHP030605	Abduction Pommel (12)	255

Fixed Centre Mount Legrest (Split Piece Footplates)			£
<input type="checkbox"/>	SHP050140	Fixed Centre Mount Legrest (13)	0
<input type="checkbox"/>	SHP050072	Footplate, Dual, Height & Angle-Adjustable 260 x 150mm	0
<input type="checkbox"/>	SHP050092	Calf Pad, Integral Foam, Pair	0

(1) Not with seat height 415 mm

(2) Not with Seat height 440 AND 130004 Powered Tilt 0° - 30° / 130007 Powered Tilt 0° - 50°

(3) Not with SHP040070 Armrest, Reclining, Multi adjustable, Flip back; SHP040065 Armrest Reclining, Multi-Adjustable, Flip back with extra support bracket

(4) Seat Width -20mm, Not with SHP040070 Armrest, Reclining, Multi adjustable, Flip back; SHP040065 Armrest Reclining, Multi-Adjustable, Flip back with extra support bracket

(5) Seat Width -80 mm, Only With SHP040021 Armrest, Single Post, Removable

(6) Only With SHP040070 Armrest, Reclining, Multi-Adjustable, Flip Back

(7) With SHP040019 Flip Back Armrest, Seat Width +20 mm

(8) Not With SHP110052 Joystick Mount, Right, SHP090110 Rear View Mirror, Right

(9) Not With SHP110053 Joystick Mount, Left, SHP090118 Rear View Mirror, Left

(10) Not With SHP050173-175, Legrest, Swing Away, SHP050080 / SHP050172 SEDEO PRO, Hanger Receivers Only, SHP090122 Transfer Bar

(11) Not With Seat To Floor Height 415 mm, Not With SHP040035 Side Guard, Wedge Shaped, Padded

(12) Not With Fixed or Power Centre Mount Legrests

(13) In Combination With Tilt, Choose SHP050092 Calf Pad, Integral Foam, Pair

(14) Actuator Control & Lights Module R-net ISM-6L required

Powered Centre Mount		£
<input type="checkbox"/>	SHP050076 Centre Mount Legrest, Powered, 90° - 180°, Elevate & Extend (13)	1470
<input type="checkbox"/>	SHP050071 Centre Mount Legrest, Powered, 90° - 180°, Elevate & Extend, Footplate to Floor Switch (13) (14)	1470
Footplates		
<input type="checkbox"/>	SHP050072 Footplate, Dual, Height & Angle Adjustable, 260 x 150 mm	0
<input type="checkbox"/>	SHP050073 Footplate, Dual, Height & Angle Adjustable, Large, 310 x 170 mm	0
Calf Pad		
<input type="checkbox"/>	SHP050092 Calf Pad, Padded, Pair	0

SEDEO Pro Swing Away Legrest		£
<input type="checkbox"/>	SHP050173 SEDEO PRO, Legrest HD, 70°, Swing Away	0
<input type="checkbox"/>	SHP050174 SEDEO PRO, Legrest HD, Manual Articulating, Swing Away	485
<input type="checkbox"/>	SHP050175 SEDEO PRO, Legrest HD, Power Articulating, Swing Away & Auto Connect	805
Footplates		
<input type="checkbox"/>	SHP050090 Footplate, For Seat Width 480 - 520 mm	0
<input type="checkbox"/>	SHP050091 Footplate, For Seat Width 540 - 560 mm	0
<input type="checkbox"/>	SHP050081 Footplates, Coupled (1) (2)	55
Hanger Receivers		
<input type="checkbox"/>	SHP050170 SEDEO PRO, Hanger Receivers, High, For Seat Cushion Thicknes 70mm And Higher	0
<input type="checkbox"/>	SHP050171 SEDEO PRO, Hanger Receivers, Low, For Seat Cushion Thicknes 50mm And Lower	0
Calf Pad / Strap		
<input type="checkbox"/>	SHP050082 Calf Support, Integral Foam, Depth-Adjustable, Pair	100
<input type="checkbox"/>	SHP050083 Calf Support, Padded, Depth-Adjustable, Pair	135
<input type="checkbox"/>	SHP050084 Heel Loop	0
<input type="checkbox"/>	SHP050085 Calf Strap (2)	0
<input type="checkbox"/>	SHP050200 Knee Pad, Integral Foam	75
Amputee Support		
<input type="checkbox"/>	SHP050142 Amputee Support, Integral Foam, Right	140
<input type="checkbox"/>	SHP050143 Amputee Support, Integral Foam, Left	140
<input type="checkbox"/>	SHP050013 Amputee Support, Padded, Right (3)	140
<input type="checkbox"/>	SHP050014 Amputee Support, Padded, Left (3)	140

R-Net Controls		£
<input type="checkbox"/>	SHP110015 Remote Joystick, R-Net, LED, Actuators	0
<input type="checkbox"/>	SHP110016 Remote Joystick, R-Net, LED, Actuators And Lights	0
<input type="checkbox"/>	SHP110017 Remote Joystick, R-Net, Colour LCD	140
<input type="checkbox"/>	SHP110018 Remote Joystick, R-Net, Colour LCD, Heavy Duty	1200
<input type="checkbox"/>	SHP110013 Remote Joystick, R-Net, High Resolution Colour LCD, Advanced, Infrared Control Included	280
<input type="checkbox"/>	SHP110023 Remote Joystick, R-Net, High Resolution Colour LCD, Advanced, Infrared And Bluetooth Control Included	490
R-Net Power Module & ISM		
<input type="checkbox"/>	SHP110003 Power Module, R-Net, 90 A	0
<input type="checkbox"/>	SHP110032 Actuator Control & Lights Module, R-Net, ISM-6L	0
R-Net Gyro		
<input type="checkbox"/>	SHP110009 Traction Control System For R-Net (GYRO Technology)	400

Joystick Mount		£
<input type="checkbox"/>	SHP110050 Joystick Mounting, Fixed	0
<input type="checkbox"/>	SHP110056 Joystick Mounting, Parallel, Swing Away, Capsulated	80

(1) Choose Also A Footplate

(2) Only For SHP050173 SEDEO PRO, Legrest, 70°, Swing Away

(3) Extra On A Set Of Swing Away Legrests, Attached Los

Joystick Mounting Location			£
<input type="checkbox"/>	SHP110052	Joystick Mounting Location, Right	0
<input type="checkbox"/>	SHP110053	Joystick Mounting Location, Left	0
Attendant Controls			£
<input type="checkbox"/>	SHP110021	Dual Control	425
<input type="checkbox"/>	SHP110020	Joystick Mounted As Attendant Control (1)	270
Attendant Mounting Location			
<input type="checkbox"/>	SHP110025	Attendant Mounting Location, Right	0
<input type="checkbox"/>	SHP110026	Attendant Mounting Location, Left	0
Joystick Knobs			£
<input type="checkbox"/>	SHP090505	Joystick Knob, Golf Ball, Ø 42 mm	40
<input type="checkbox"/>	SHP090506	Joystick Knob, T-Bar, 82 mm Wide	40
<input type="checkbox"/>	SHP090507	Joystick Knob, Soft Ball, Ø 70 mm	40
<input type="checkbox"/>	SHP090508	Joystick Knob, Ball, Ø 35 mm	40
<input type="checkbox"/>	SHP090045	Joystick Knob, Allround, Ball, Ø 48 mm	40
<input type="checkbox"/>	SHP090509	Joystick Knob, Chin Cup, Ø 37 mm	40
<input type="checkbox"/>	SHP090512	Joystick Knob, Stick, 80 x 15 mm	40
<input type="checkbox"/>	SHP090511	Joystick Knob, Mushroom, Ø 45 mm	40
LED Lights & Indicators			£
<input type="checkbox"/>	SHP110201	LED Light And Indicator StVZO Approved (2)	400
<input type="checkbox"/>	SHP110202	Fender Lights (2) Not for Road Use	160
Motor Packages			£
<input type="checkbox"/>	SHP180009	4-Pole Motor, 6 kph / 4 mph	0
<input type="checkbox"/>	SHP180010	4-Pole Motor, 8 kph / 5 mph	390
Drive Wheels & Castor Wheels			£
<input type="checkbox"/>	SHP070034	Drive Wheel, 14", Pneumatic, Knobbly, Black	0
<input type="checkbox"/>	SHP070033	Drive Wheel, 14", Solid, Knobbly, Black	0
<input type="checkbox"/>	SHP080310	Castor Wheel, 6", Solid, Black	0
Accent Colours			£
<input type="checkbox"/>	SHP140005	Apple Green (RAL 110 70 77)	0
<input type="checkbox"/>	SHP140004	Black (RAL 9005)	0
<input type="checkbox"/>	SHP140006	Orange (RAL 060 60 80)	0
<input type="checkbox"/>	SHP140007	Silver (RAL 9006)	0
<input type="checkbox"/>	SHP140002	Red (RAL 040 40 67)	0
<input type="checkbox"/>	SHP140003	Blue (RAL 260 40 45)	0
<input type="checkbox"/>	SHP140001	White (RAL 9010)	0
Batteries			£
<input type="checkbox"/>	SHP160001	Battery, 60 Ah, Gel, Maintenance-Free	0
<input type="checkbox"/>	SHP160005	Battery, 80 Ah, Gel, Maintenance-Free	160
Chargers			£
<input type="checkbox"/>	SHP160103	Charger, Automatic, 220 V / 10 A, UK	0

(1) Choose Also A Remote Joystick

(2) SHP110032 Actuator Control & Lights Module, R-Net, ISM-6L Required

Lap Belts			£
<input type="checkbox"/>	SHP090128	Positioning Belt, Carseat Buckle, Short, 50 x 1700 mm	0
<input type="checkbox"/>	SHP090129	Positioning Belt, Carseat Buckle, Short, 50 x 2700 mm	0

Accessories			£
<input type="checkbox"/>	SHP090110	Rear View Mirror, Right	60
<input type="checkbox"/>	SHP090118	Rear View Mirror, Left	60
<input type="checkbox"/>	SHP090009	Crutch Holder (1)	70
<input type="checkbox"/>	SHP090119	Brackets for Bag on Backrest Frame (1)	35
<input type="checkbox"/>	SHP090120	Bag for Armrest (2)	110
<input type="checkbox"/>	SHP090122	Transfer Bar (3)	155

Special Controls & Switches

5-Switch Box			£
<input type="checkbox"/>	SHP110043	#1 Ctrl+5 Switch Box With Free Assignable Functions (4)	865
#1 Ctrl+5 Switch Box Mounting Location			
<input type="checkbox"/>	SHP110180	#1 Ctrl+5 Switch Box Mounting Location, Behind Joystick	0
<input type="checkbox"/>	SHP110181	#1 Ctrl+5 Switch Box Mounting Location, Opposite Joystick	0
#1 Ctrl+5 Switch Type			
<input type="checkbox"/>	SHP110175	#1 Ctrl+5 Switch Type With Push Buttons	0
<input type="checkbox"/>	SHP110176	#1 Ctrl+5 Switch Type With Toggle Switches	0

Omni2 Input Control Module			£
<input type="checkbox"/>	SHP110116	OMNI2 Input Control Module, R-Net, Bluetooth, Infrared (5)	1170
<input type="checkbox"/>	SHP110117	OMNI2 Input Control Module, R-Net, Scanning Mode, Bluetooth, Infrared (5)	1170

Omni 2 Module Mount (11)			£
<input type="checkbox"/>	SHP110119	OMNI2 Module, Fixed Mount, Right	0
<input type="checkbox"/>	SHP110120	OMNI2 Module, Fixed Mount, Left	0
<input type="checkbox"/>	SHP110123	OMNI2 Module, Swing Away Mount, Right (6)	130
<input type="checkbox"/>	SHP110124	OMNI2 Module, Swing Away Mount, Left (7)	130

Hand Control (12)			£
<input type="checkbox"/>	SHP110301	Hand Control, Switch-It MicroGuide	3145
<input type="checkbox"/>	SHP110302	Hand Control, Switch-It MicroPilot	3495
<input type="checkbox"/>	SHP110304	Hand Control, Switch-It VersaGuide	1120
<input type="checkbox"/>	SHP110305	Hand Control, Switch-It VersaGuide EZ (8)	1330
<input type="checkbox"/>	SHP110306	Hand Control, Allround Joystick	1120
<input type="checkbox"/>	SHP110307	Hand Control, Allround Joystick Light (9)	1270
<input type="checkbox"/>	SHP110308	Hand Control, Allround Joystick Heavy Duty	2095
<input type="checkbox"/>	SHP110309	Hand Control, Allround Joystick, R-Net Bus Enabled (10)	1190
<input type="checkbox"/>	SHP110310	Hand Control, Allround Joystick Light, R-Net Bus Enabled (9) (10)	1400
<input type="checkbox"/>	SHP110311	Hand Control, Allround Joystick Heavy Duty, R-Net Bus Enabled (10)	2235

(1) Not With SHP030215 Powered Recline, 90° - 130°

(2) Not With SHP040070 Armrest, Reclining, Multi-Adjustable, Flip Back and SHP040065 Armrest Reclining, Multi-Adjustable, Flip Back with extra support bracket

(3) Not With SHP020342/SHP020341 Knee Support, Lateral, Swing Away, Right/Left & SHP050173/SHP050174/ SHP050175 SEDEO PRO, Legrest, Swing Away & SHP050080/ SHP050172 SEDEO PRO, Hanger Receivers Only

(4) If combined with attendant control, Ctrl+5 will be mounted to the armrest

(5) Remote Full VGA Display With Infrared, Separate Input Box With Bluetooth, USB Charger And Speaker

(7) Not With SHP090110 Rear View Mirror Right

(8) Not With SHP090118 Rear View Mirror Left

(9) For Switch Selection See Section E

(10) Not With SHP090512 Joystick Knob, Stick, 80 x 15 mm, SHP090511 Joystick Knob, Mushroom, Ø 45 mm, SHP090506 Joystick Knob, T-Bar, 82 mm Wide

(11) Choose OMNI2 Mount And Input Device in Sections Below And Choose Switches in Section E

(12) If Used With Standard Joystick Selection, It Must Be Mounted Opposite Side

Hand Control Mounting		£
<input type="checkbox"/>	SHP110342 Hand Control Mounting, Link-It To The Space Disk, Hand Tray (1)	700
Hand Control Mounting		
<input type="checkbox"/>	SHP110350 Hand Control Mounting, Link-It To The Armrest, Fixed	210
<input type="checkbox"/>	SHP110353 Hand Control Mounting, Link-It To The Armrest, Swing Away	220
Hand Control Mounting Location		
<input type="checkbox"/>	SHP110362 Hand Control Mounting Location, Right	0
<input type="checkbox"/>	SHP110361 Hand Control Mounting Location, Left	0
Head Control (2)		£
Proportional Input Devices		
<input type="checkbox"/>	SHP110370 Head Control, Switch-It Dual Pro, 3-Way Proportional, Facial Pads	5950
<input type="checkbox"/>	SHP110371 Head Control, Switch-It Dual Pro, 3-Way Proportional, Spot Pads	5950
Chin Control		£
<input type="checkbox"/>	SHP110470 Chin Control, Switch-It MicroGuide	3145
<input type="checkbox"/>	SHP110471 Chin Control, Switch-It MicroPilot	3495
<input type="checkbox"/>	SHP110472 Chin Control, Switch-It VersaGuide	1120
<input type="checkbox"/>	SHP110473 Chin Control, Switch-It VersaGuide EZ	1330
<input type="checkbox"/>	SHP110474 Chin Control, Allround Joystick	1120
<input type="checkbox"/>	SHP110475 Chin Control, Allround Joystick Light (3)	1330
<input type="checkbox"/>	SHP110476 Chin Control, Allround Joystick, R-Net Bus Enabled (4)	1190
<input type="checkbox"/>	SHP110477 Chin Control, Allround Joystick Light, R-Net Bus Enabled (3) (4)	1400
Chin Control Mounting		£
<input type="checkbox"/>	SHP110480 Chin Control Mounting, Link-It To Chin Harness, Small (5)	630
<input type="checkbox"/>	SHP110481 Chin Control Mounting, Link-It To Chin Harness, Large (5)	700
<input type="checkbox"/>	SHP110483 Chin Control Mounting, Link-It To Manual Chin Boom, Backrest Mount, Right	560
<input type="checkbox"/>	SHP110484 Chin Control Mounting, Link-It To Manual Chin Boom, Backrest Mount, Left	560
<input type="checkbox"/>	SHP110485 Chin Control Mounting, Link-It To Power Chin Boom, Backrest Mount, Right	2795
<input type="checkbox"/>	SHP110486 Chin Control Mounting, Link-It To Power Chin Boom, Backrest Mount, Left	2795
<input type="checkbox"/>	SHP110487 Chin Control Mounting, Link-It To Headrest, Swing Away, Left	490
<input type="checkbox"/>	SHP110488 Chin Control Mounting, Link-It To Headrest, Swing Away, Right	490
Link-It Extension Kit		
<input type="checkbox"/>	SHP110590 Link-It Extension Kit, Base Module, 100 mm Extension Tube, 5 Cable Clips	170
<input type="checkbox"/>	SHP110591 Link-It Extension Kit, Base Module, 150 mm Extension Tube, 5 Cable Clips	210
Sip 'N Puff (6)		£
<input type="checkbox"/>	SHP110510 Sip 'N Puff With Link-It Hardware, Swing Away Mount To The Headrest (7)	700

(1) Only For Micro Series Joysticks

(2) For Switch Selection See next page. Choose A Whitmyer Headrest Mount

(3) Only Use With Default Joystick Knob Or Joystick Knob, Allround, Ball Ø 48 mm

(4) OMNI2 Module Not Required, Obligated To Have Master R-Net Module

(5) If You Need Switches Mounted On The Harness, You Can Only Use Allround Joystick

(6) For Switch Selection See next page

(7) Headrest To Be Selected

#1 Switches - Section H (1)			£
<input type="checkbox"/>	SHP110611	#1 Single Switch, Ø 64 mm, Black	120
<input type="checkbox"/>	SHP110612	#1 Single Switch, Ø 64 mm, Red	120
<input type="checkbox"/>	SHP110613	#1 Single Switch, Ø 64 mm, Green	120
<input type="checkbox"/>	SHP110523	#1 Egg Switch – Black	135
<input type="checkbox"/>	SHP110524	#1 MicroLite Switch – Black	195
<input type="checkbox"/>	SHP110651	#1 Single Switch, Ø 30 mm, Black	120
<input type="checkbox"/>	SHP110654	#1 Single Switch, Ø 30 mm, Blue	120
<input type="checkbox"/>	SHP110653	#1 Single Switch, Ø 30 mm, Green	120
<input type="checkbox"/>	SHP110652	#1 Single Switch, Ø 30 mm, Red	120
<input type="checkbox"/>	SHP110655	#1 Single Switch, Ø 30 mm, Yellow	120
<input type="checkbox"/>	SHP110534	#1 Twister Switch – Black	195
<input type="checkbox"/>	SHP110535	#1 Twister Switch – Green	195
<input type="checkbox"/>	SHP110536	#1 Twister Switch – Red	195
<input type="checkbox"/>	SHP110537	#1 Satellite Switch – Black (7)	210
<input type="checkbox"/>	SHP110538	#1 Satellite Switch – Green (7)	210
<input type="checkbox"/>	SHP110539	#1 Satellite Switch – Red (7)	210

#1 Switch Mounting			£
<input type="checkbox"/>	SHP110560	#1 Switch Mounting, Link-It To The Hand Control (2)	210
<input type="checkbox"/>	SHP110566	#1 Switch Mounting, Link-It To The Armrest (Opposite Site As The Hand Control) (3)	210
<input type="checkbox"/>	SHP110561	#1 Switch Mounting, Link-It To The Headrest	280
<input type="checkbox"/>	SHP110563	#1 Switch Mounting, Link-It To The Seat Rail (4)	210
<input type="checkbox"/>	SHP110565	#1 Switch Mounting, Link-It To Chin Boom (5)	210
<input type="checkbox"/>	SHP110569	#1 Omit Switch Mounting (6)	0

#1 Switch Mount Location			£
<input type="checkbox"/>	SHP110572	#1 Switch Mounting Location, Right	0
<input type="checkbox"/>	SHP110571	#1 Switch Mounting Location, Left	0

#2 Switches - Section H (1)			£
<input type="checkbox"/>	SHP110631	#2 Single Switch, Ø 64 mm, Black	120
<input type="checkbox"/>	SHP110632	#2 Single Switch, Ø 64 mm, Red	120
<input type="checkbox"/>	SHP110633	#2 Single Switch, Ø 64 mm, Green	120
<input type="checkbox"/>	SHP110088	#2 Egg Switch – Black	120
<input type="checkbox"/>	SHP110089	#2 MicroLite Switch – Black	195
<input type="checkbox"/>	SHP110671	#2 Single Switch, Ø 30 mm, Black	120
<input type="checkbox"/>	SHP110674	#2 Single Switch, Ø 30 mm, Blue	120
<input type="checkbox"/>	SHP110673	#2 Single Switch, Ø 30 mm, Green	120
<input type="checkbox"/>	SHP110672	#2 Single Switch, Ø 30 mm, Red	120
<input type="checkbox"/>	SHP110675	#2 Single Switch, Ø 30 mm, Yellow	120
<input type="checkbox"/>	SHP110095	#2 Twister Switch – Black	195
<input type="checkbox"/>	SHP110096	#2 Twister Switch – Green	195
<input type="checkbox"/>	SHP110097	#2 Twister Switch – Red	195
<input type="checkbox"/>	SHP110111	#2 Satellite Switch – Black (7)	210
<input type="checkbox"/>	SHP110112	#2 Satellite Switch – Green (7)	210
<input type="checkbox"/>	SHP110113	#2 Satellite Switch – Red (7)	210

(1) Select 1 Switch And Mounting Position Accordingly. For A Second Switch Please Use #2 Switch Selection Part Below. Single Function For Special Input Devices
 (2) Not With HD Allround Joystick Or Space Disk Mount
 (3) Armrest Rail Mount

(4) Knee, If 2 Switches Selected Only Left And Right Mount
 (5) No Left / Right Switch Selection Required
 (6) Switches Attached Lose
 (7) Only For Allround Joystick, Maximum 2 Selectable

#2 Switch Mounting			£
<input type="checkbox"/>	SHP110098	#2 Switch Mounting, Link-It To The Hand Control (1)	210
<input type="checkbox"/>	SHP110099	#2 Switch Mounting, Link-It To The Armrest (Opposite Site As The Hand Control) (2)	210
<input type="checkbox"/>	SHP110100	#2 Switch Mounting, Link-It To The Headrest	280
<input type="checkbox"/>	SHP110104	#2 Switch Mounting, Link-It To The Seat Rail (3)	210
<input type="checkbox"/>	SHP110105	#2 Switch Mounting, Link-It To Chin Boom (4)	210
<input type="checkbox"/>	SHP110106	#2 Omit Switch Mounting (5)	0

#2 Switch Mount Location			£
<input type="checkbox"/>	SHP110108	#2 Switch Mounting Location, Right	0
<input type="checkbox"/>	SHP110107	#2 Switch Mounting Location, Left	0

Multi-Function Switches (7)			£
<input type="checkbox"/>	SHP110550	Mini 2 Button (Stereo Plug) (8)	385
<input type="checkbox"/>	SHP110551	Mini 2 Button (Dual Mono Plug) (8)	385
<input type="checkbox"/>	SHP110552	Mini 4 Button InLine (DB9 Connector) (8)	645
<input type="checkbox"/>	SHP110554	4-Way Toggle Switched Joystick	645

Accessories			£
<input type="checkbox"/>	SHP110580	Mono Extension Cable (3.5 mm Jack Socket) (2 m)	140
<input type="checkbox"/>	SHP110581	Stereo Extension Cable (3.5 mm Jack Socket) (2 m)	180
<input type="checkbox"/>	SHP110582	Extension Cable DB9 (1 m)	170
<input type="checkbox"/>	SHP110583	Y-Cable, 2 Female Mono To 1 Male Stereo	210

Safety Devices			£
<input type="checkbox"/>	SHP090147	Remote Stop Switch (30 m Activation Range, Fit Via OBC Port)	980

Comments For Personalisations



Q700 M HD with SEDEO PRO		£
Base Price		
Total Of Additional Equipment		
Packaging And Transport		
Exclusive VAT		
VAT		
Total Price		

Date			Signature
Day	Month	Year	
<input type="text"/>	<input type="text"/>	<input type="text"/>	<input type="text"/>

[Send E-Mail](#)

(1) Not With HD Allround Joystick Or Space Disk Mount
 (2) Armrest Rail Mount
 (3) Knee, If 2 Switches Selected Only Left And Right Mount
 (4) No Left / Right Switch Selection Required
 (5) Switches Attached Lose

(6) Only For Allround Joystick, Maximum 2 Selectable
 (7) Used In Combination With Assignable Button Programming
 (8) Attached Lose, When Ordered, No Link-It Fixation Possible